MIKRO 220 Robotic

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These instructions for initial operation and operation are not for users but are only meant for the manufacturer of the automated system.

For more information see chapter 1 on page 39.

Andreas Hettich GmbH & Co. KG
Föhrenstraße 12, D-78532 Tuttlingen / Germany
Phone +49 (0)7461 / 705-0
Fax +49 (0)7461 / 705-1125
info@hettichlab.com, service@hettichlab.com
www.hettichlab.com

© 2016 by Andreas Hettich GmbH & Co. KG

All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.

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Standards which apply to this device

The device is an incomplete high-tech centrifuge. It is subject to extensive testing and certification procedures according to the following standards in their respectively valid version:

Electrical and mechanical safety for design and final testing:

Standard series: IEC 61010 (conform to standards of DIN EN 61010)

- IEC 61010-1 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements" (Pollution Degree 2, Excess-voltage category II)
- IEC 61010-2-020 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-020: Particular requirements for laboratory centrifuges"

Electromagnetic Compatibility:

- EN 61326-1 "Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements"

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Contents

| | | |
|--------|---|----|
| 1 | How to use these instructions for initial operation and operation | 39 |
| 2 | Symbol meanings | 39 |
| 3 | Intended use | 40 |
| 4 | Remaining risks | 40 |
| 5 | Technical specifications | 41 |
| 6 | Notes on safety | 42 |
| 7 | Transport and storage | 43 |
| 7.1 | Transport | 43 |
| 7.2 | Storage | 43 |
| 8 | Scope of delivery | 44 |
| 9 | Unpacking the centrifuge | 44 |
| 10 | Removing the transport securing device | 44 |
| 11 | Initial operation | 45 |
| 12 | Interface | 45 |
| 13 | Opening and closing the hatch | 46 |
| 13.1 | Opening the hatch | 46 |
| 13.2 | Closing the hatch | 46 |
| 14 | Opening and closing the lid | 47 |
| 14.1 | Opening the lid | 47 |
| 14.2 | Closing the lid | 47 |
| 15 | Installation and removal of the rotor | 48 |
| 16 | Loading the rotor | 48 |
| 17 | Control and display elements | 49 |
| 17.1 | Control knob | 49 |
| 17.2 | Keys and setting options | 49 |
| 18 | Enter centrifugation parameters | 51 |
| 18.1 | Running time | 51 |
| 18.2 | Starting runtime timing | 51 |
| 18.3 | Speed (RPM) | 51 |
| 18.4 | Relative centrifugal force (RCF) and centrifuging radius (RAD) | 52 |
| 18.5 | Start-up and run-down parameters | 52 |
| 18.5.1 | Start-up stage and start-up time | 52 |
| 18.5.2 | Braking stage and run-down time | 52 |
| 18.5.3 | Brake switch-off speed | 52 |
| 19 | Programming | 53 |
| 19.1 | Inputting or changing programs | 53 |
| 19.2 | Calling up programs | 53 |
| 19.3 | Write protection for programs | 53 |
| 19.4 | Program linkage | 54 |
| 19.4.1 | Activating or deactivating program linkage | 54 |
| 19.4.2 | Linking programs or changing a program linkage | 54 |
| 19.4.3 | Calling up program linkage | 55 |

| | | |
|--------|--|----|
| 19.5 | Automatic temporary storage | 55 |
| 20 | Centrifugation | 55 |
| 20.1 | Centrifugation with preselected time | 55 |
| 20.2 | Continuous operation | 56 |
| 21 | Change the settings during the centrifugation run | 56 |
| 22 | Integral RCF | 56 |
| 22.1 | Query integral RCF | 56 |
| 22.2 | Activate or deactivate the integral RCF display | 57 |
| 23 | Emergency Stop | 57 |
| 24 | Cycle counter | 58 |
| 24.1 | After starting the first centrifugation run, enter the maximum permitted number of running cycles | 58 |
| 24.2 | Reset cycle counter to "0" and enter the maximum permitted number of cycles | 58 |
| 24.3 | Activating the cycle counter | 59 |
| 25 | Activate or deactivate the function, "Dual time mode" | 60 |
| 26 | Activating or deactivating start-up and run-down times | 60 |
| 27 | Acoustic signal | 61 |
| 28 | Displayed centrifugation data after switch-on | 61 |
| 29 | Centrifuge address | 61 |
| 30 | Query the operating hours, centrifugation runs and cycle counter | 62 |
| 31 | Query the system information | 62 |
| 32 | Immediate display of centrifugation data after switch-on | 62 |
| 33 | Key-operated switch | 63 |
| 34 | Setting rotor position 1 | 64 |
| 35 | Relative centrifugal force (RCF) | 64 |
| 36 | Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm ³ | 65 |
| 37 | Rotor recognition | 65 |
| 38 | Maintenance and servicing | 66 |
| 38.1 | Centrifuge (housing, lid and centrifuging chamber) | 66 |
| 38.1.1 | Surface cleaning and care | 66 |
| 38.1.2 | Surface disinfection | 66 |
| 38.1.3 | Removal of radioactive contaminants | 66 |
| 38.2 | Rotor and accessories | 67 |
| 38.2.1 | Cleaning and care | 67 |
| 38.2.2 | Disinfection | 67 |
| 38.2.3 | Removal of radioactive contaminants | 67 |
| 38.2.4 | Rotors and accessories with limited service lives | 67 |
| 38.3 | Autoclaving | 67 |
| 38.4 | Centrifuge containers | 67 |
| 39 | Faults | 68 |
| 40 | Anhang / Appendix | 70 |
| 40.1 | Rotoren und Zubehör / Rotors and accessories | 70 |

1 How to use these instructions for initial operation and operation

- These instructions are instructions for initial operation and operation for a incomplete centrifuge.
- These instructions for initial operation and operation are not for users but are only meant for the manufacturer of the automated system.
- Before integrating the centrifuge into the automated system, read and observe the interface description AH5680-01EN and the instructions for initial operation and operation.

2 Symbol meanings

Symbol on the device:

Attention, general hazard area.

Before using the device, make sure you read the instructions for initial operation and operation and observe the safety information.

Symbol in this document:

Attention, general hazard area.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol on the device:

Beware of squeezing the hands.

Symbol on the device and in this document:

Beware of biohazard.

Symbol on the device:

Switch positions of the key-operated switch.

Symbol on the device:

RS232 interface.

Symbol in this document:

This symbol refers to important circumstances.

Symbol on the device and in this document:

Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).

Applies in the countries of the European Union, as well as in Norway and Switzerland.

3 Intended use

This device is an incomplete centrifuge.

The centrifuge is only intended for use in an automated system with loading/unloading by a robot. The centrifuge is controlled by the automated system via the interface. A centrifugation run can only be started via the interface.

The centrifuge was designed only for the separation of materials or mixtures with a density of no more than 1.2 kg/dm³.

The centrifuge is only meant for this purpose.

Another use or one which goes beyond this, is considered to be non-intended. The company Andreas Hettich GmbH & Co. KG is not liable for damage resulting from this.

Observing the instructions for initial operation and operation and complying with the inspection and maintenance work is also a part of the intended use.

The manufacturer of the automated system is responsible for the safety of the overall system including the centrifuge integrated in the system.

4 Remaining risks

The device is built according to the state-of-the-art and the recognized safety regulations. If used and handled improperly, there could be life-threatening danger to the user or third parties, or the device could be impaired or there could be other property damage. The device is only to be used for its intended purpose and only when it is in safe working condition.

Malfunctions which could affect safety must be corrected immediately.

5 Technical specifications

| | | |
|---|--|--------------|
| Manufacturer | Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen | |
| Model | MIKRO 220 Robotic | |
| Type | 2350 | 2350-01 |
| Mains voltage ($\pm 10\%$) | 200-240 V 1~ | 100-127 V 1~ |
| Mains frequency | 50-60 Hz | |
| Connected load | 480 VA | |
| Max. capacity | 24 x 2 ml | |
| Allowed density | 1.2 kg/dm ³ | |
| Speed (RPM) | 13000 | |
| Force (RCF) | 18516 | |
| Kinetic energy | 8000 Nm | |
| Obligatory inspection (BGR 500) | no | |
| Ambient conditions (EN / IEC 61010-1) | <p>Indoors only</p> <p>Up to 2000 m above sea level</p> <p>5°C to 40°C</p> <p>Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.</p> | |
| – Set-up site | | |
| – Altitude | | |
| – Ambient temperature | | |
| – Humidity | | |
| – Excess-voltage category (IEC 60364-4-443) | II | |
| – Pollution degree | 2 | |
| Device protection class | I | |
| Not suitable for use in explosion-endangered areas. | | |
| EMC | | |
| – Emitted interference, Interference immunity | EN / IEC 61326-1, Class B; FCC Class B | |
| Noise level (dependent on rotor) | ≤ 67 dB(A) | |
| Dimensions | | |
| – Width | 410 mm | |
| – Depth | 414 mm | |
| – Height | 267 mm | |
| Weight | approx. 27 kg | |

6 Notes on safety

No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.

- The centrifuge may only be operated in an automated system with loading and unloading by a robot.
- The centrifuge should be integrated into the automated system on a good, stable base.
- The manufacturer of the automated system is responsible for the safety of the overall system including the centrifuge integrated in the system.
- Before using the centrifuge absolutely check the rotor for firm placement.
- When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
- Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.
- The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.
- With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.
- For centrifuges without temperature control, when the room temperature is increased and/or if the device is frequently used, the centrifuging chamber could be heated up. Therefore, it can't be ruled out that the sample material might be changed due to the temperature.

- **Only the manufacturer of the automated system may integrate the centrifuge into the automated system and put the centrifuge into operation. Before carrying out these activities, read and observe the interface description AH5680-01EN and the instructions for initial operation and operation.**
- In addition to the instructions for initial operation and operation and the binding regulations for accident prevention, the recognized professional rules for safe and proper work are to be observed. The instructions for initial operation and operation are to be supplemented with instructions based on the national regulations for accident prevention and environmental protection of the country where the device is used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room, the centrifuge must warm up for at least 3 hours in the warm room before it may be connected to the mains.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"). Before centrifuge vessels are used which are not listed in the chapter "Appendix, Rotors and accessories", the user must make sure they can be used by asking the manufacturer.
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.

For hazardous substances centrifuge containers with special screw caps must strictly be used. In addition to the screw cap centrifuge containers, for materials in hazard category 3 and 4 a biosafety system must be used (see the World Health Organisation's "Laboratory Biosafety Manual").

Without the use of a biosafety system the centrifuge is not microbiologically sealed in the sense of the EN / IEC 61010-2-020 standard.

No biosafety systems are available for this centrifuge.

- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:
EN / IEC 61010-1 and EN / IEC 61010-2-020 as well as their national deviations.
- The safety and reliability of the centrifuge are only guaranteed if:
 - the centrifuge is operated according to the instructions for initial operation and operation.
 - the electrical installation at the set-up site of the centrifuge meets the requirements of the EN / IEC regulations.
 - the tests for device safety required in the respective countries, e.g. in Germany in acc. with DGUV (Deutsche Gesetzliche Unfallversicherung or Engl.: German Statutory Accident Insurance) regulation 3, are carried out by an expert.
- With centrifuges for robotic use please pay attention the notes of the key operated switch.

7 Transport and storage

7.1 Transport

Before transporting the device, the transport securing device must be installed.

When the device and accessories are transported, the following ambient conditions must be complied with:

- Ambient temperature: -20°C to +60°C
- Relative humidity: 20% to 80%, non-condensing

7.2 Storage

The device and the accessories may only be stored in closed and dry rooms.

When the device and accessories are stored, the following ambient conditions must be complied with:

- Ambient temperature: -20°C to +60°C
- Relative humidity: 20% to 80%, non-condensing

8 Scope of delivery

The following accessories are delivered with the centrifuge:

- 1 Power cord
- 1 Hex. pin driver (5 mm)
- 2 Key for key-operated switch
- 1 Test programme Robot
- 2 Label for voltage and frequency
- 1 instructions for initial operation and operation
- 1 Notes on moving the equipment safely

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

9 Unpacking the centrifuge

- Lift the carton upward and remove the padding.


-  Do not lift the centrifuge at the operating unit or at the holder of the operating unit. Observe the weight of the centrifuge, refer to chapter "Technical specifications".

Lift the centrifuge on both sides with an appropriate number of helpers and place it on the laboratory table.

10 Removing the transport securing device

-  It is imperative that the transport securing device be removed. Keep the transport securing device in a safe place since it must be installed again before transporting the device. The device may only be transported with the transport securing device installed.

To protect the device from damage during transport, the motor is fixed in place. This transport securing device must be removed when the device is put into operation.

- Remove the three screws (a) and spacer sleeves (b).
- Remove the three transport securing devices (c).

 The transport securing device is installed in the opposite order.

11 Initial operation

The centrifuge may only be operated in an automated system with loading and unloading by a robot. The centrifuge is controlled by the automated system via the interface. A centrifugation run can only be started via the interface.

The manufacturer of the automated system is responsible for the safety of the overall system including the centrifuge integrated in the system.

- Remove the transport securing device; see the "Removing the transport securing device" chapter.
- **Set up and level the centrifuge in the automated system so that it is stable. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to EN / IEC 61010-2-020.**

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

- Ventilation openings may not be blocked.
A distance of 300 mm must be maintained from the ventilation slots and openings of the centrifuge.
- Connect the RS232 centrifuge interface to the automated system with an RS232 connecting cable (not included in the scope of delivery).
- Check whether the mains voltage tallies with the statement on the type plate.
- Connect the centrifuge with the power cord to a standard mains socket. For connection ratings refer to Chapter "Technical specifications".
- Switch on the mains switch (switch position "I"). The LEDs in the keys are blinking.
The following displays appear one after the other:
 1. The centrifuge model

On the MIKRO 220 Robotic, **MIKRO 220 POS** is displayed.

2. The program version
3. The rotor code (Rotor), the maximum speed of the rotor (Nmax) and a centrifuging radius (R) of the last rotor detected by the rotor detector.
4. The centrifugation data of the last used program or program 1.

12 Interface

The device is equipped with an RS232 interface.

The RS232 interface is labeled with the  symbol.

The centrifuge can be controlled and data queried via this interface.

The LED in the  key lights up during data communication.

13 Opening and closing the hatch

Only authorized skilled personnel may operate the device in the "TEACH" key position.
The key is to be kept in a safe place so that it is protected from unauthorized access.

The hatch can only be opened and closed via the operating unit (key position "TEACH") or via the interface (key position "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

The hatch opens and closes via the motor.

13.1 Opening the hatch

The hatch can only be opened if the rotor is stopped.

If the hatch is blocked during opening, the hatch drive stops.

Afterwards, **POS-ERROR 45 OVERCURRENT** is displayed (see "Malfunctions" chapter).

- Turn the key to the "TEACH" position. For example, **Teach Open =Close** might be displayed.
- Keep the  key pressed until the hatch is completely open. **Teach =Open Close** is displayed (=: position of the hatch).

13.2 Closing the hatch

To prevent injuries to the hands of the operator and damage to the robot arm, the drive stops the hatch and opens it again completely if the hatch is blocked for a short time during closing.

Afterwards, **POS-ERROR 45 OVERCURRENT** is displayed (see "Malfunctions" chapter).

- Turn the key to the "TEACH" setting. **Teach =Open Close** is displayed, for example.
- Keep the  key pressed until the hatch is completely closed. **Teach Open =Close** is displayed (=: Position of the hatch).


14 Opening and closing the lid

14.1 Opening the lid

The lid may only be opened when the rotor is at rest.
Always open the lid completely, because it will not held in position.

- Look through the window (a) in the lid to be sure that the rotor is at standstill.
- Unscrew the screw (b) at the lid in a counter clockwise direction. The left LED in the button **OPEN / STOP** extinguishes.
- Open the lid.

14.2 Closing the lid


Do not bang the lid shut.

- Close the lid and tighten the screw (b) at the lid by turning in a clockwise direction. The left LED in the button **OPEN / STOP** lights up.

15 Installation and removal of the rotor

After each installation of the rotor, the rotor-position 1 must be set again, see chapter "Set rotor-position 1".

- Look through the window in the lid to be sure that the rotor is at standstill.
- Clean the motor shaft (C) and the rotor drilling (A), and lightly grease the motor shaft afterwards. Dirt particles between the motor shaft and the rotor hinder a perfect seating of the rotor and cause an irregular operation.
- Place the rotor vertically on the motor shaft. The two carrier pins (B) at the underside of the rotor may not lay on the carrier (D) when tightening up the rotor.
- Tighten the rotor tension nut with the supplied wrench by turning in a clockwise direction.
- Check the rotor for firm seating.
- Loosening the rotor: Loosen the tension nut by turning in a counter clockwise direction, and turning until the working point for lifting. After passing the working point for lifting the rotor is loosened from the motor shaft cone. Turn the tension nut until the rotor is able to be lifted from the motor shaft.


16 Loading the rotor


Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2).

- Check the rotor for firm seating.
- The rotors must be loaded symmetrically. The centrifuge containers have to be distributed evenly on all rotor positions. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

Example:

Rotor is evenly loaded

Not permitted!
Rotor is not evenly loaded


- The centrifuge containers may only be filled outside of the centrifuge.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.
- When loading the rotor, no liquid may enter the rotor or the centrifuging chamber.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.
- The weight of the permissible filling quantity is specified on each rotor. This weight may not be exceeded.

17 Control and display elements

The following description of the operation of the centrifuge is only intended for the manufacturer of the automated system. This description may not be passed on to the user of the automated system.

The user may not manually operate the centrifuge. The only exception is performing an emergency stop, see chapter "Emergency stop".

17.1 Control knob

For setting the individual parameters.

Turning anticlockwise reduces the value. Turning clockwise increases the value.

17.2 Keys and setting options

- Running time, parameters **t/hms**.
 - h**: hours. Adjustable from 1 h to 99 h, in 1 hour increments.
 - m**: minutes. Adjustable from 1 min to 59 min, in 1 minute increments.
 - s**: seconds. Adjustable from 1 s to 59 s, in 1 second increments.
- Continual running " ∞ "
- Set runtime timing This setting can only be made if the function "Dual time mode" is activated. See chapter "Activate or deactivate function "Dual time mode". This function is activated ex works. You can decide whether the runtime is timed immediately after the start of the centrifugation run or as soon as the set speed is reached.
 - Timing begins at Start** = The runtime begins to count immediately after the centrifugation run starts.
 - Timing begins at Speed** = The runtime only begins to count after the set speed is reached. This is indicated on the display by the symbol **F**, to the left next to the time.

- Revolution, parameters **rpm**. Adjustable from 50 rpm to a maximum rotor speed of the rotor (Nmax), in increments of 10. For maximum rotor speed, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

- Relative centrifugal force, parameter **RCF**. The RCF is shown in brackets $\rangle \langle$. The LED in the button lights up. The figure that is adjustable is the figure that produces a revolution speed between 50 rpm and the maximum rotor speed (Nmax). Adjustable in increments of 1.
- Centrifuging radius, parameters **RAD**. Adjustable from 10 mm to 330 mm in 1 millimeter increments. For centrifuging radius, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
- Requesting information from integral RCF. You can only query integral RCF if the integral RCF display is activated. See chapter "Activate or deactivate integral RCF display".

Start-up and run-down parameters

- Run-in steps, parameters \swarrow . Step 9 = shortest run-in time, ... Step 1 = longest run-in time.
- Start-up time, parameter $\swarrow t$. Can be set in 1-second steps. The adjustable time range depends on the set speed. It is only possible to set start-up times if these have been activated. See the chapter "Activating or deactivating start-up and run-down times".
- Braking stages, parameter \searrow . 1-9 = Linear braking curve Step 9 = shortest run-out time, ... Step 1 = long run-out time, step 0 = non-braked run-out.
- Run-down time, parameter $\searrow t$. Can be set in 1-second steps. The adjustable time range depends on the set speed. It is only possible to set run-down times if these have been activated. See the chapter "Activating or deactivating start-up and run-down times".

- Brake switch-off revolution, parameters **N Brake**.
Adjustable from 50 rpm to the maximal rotor speed (Nmax) in increments of 10.
Once the rotor speed has been reached, non-braked run-out begins.

◦PROG

- Call up programs and program linkages; parameter **RCL**.
Programs: Program locations 1 to 99. Program linkages: Program locations A to Z.
- Save programs and program linkages; parameter **STO**.
99 programs can be saved (Program locations 1 to 99). Program location 0 is used for temporary storage for centrifuge data from the last centrifugation run. No programs can be saved at this program location. It is possible to save 25 program linkages (program locations A to Z; there is no program location J). A program linkage can consist of 20 programs.
- Link programs; parameter **EDIT**.
- Call up the "Machine menu" (keep the key press for 8 seconds).
- Scroll forwards in the menus.

◦START

- Saving entries and changes.
- In the "Machine Menu", call up the submenus.
- No centrifugation run can be started. After pressing the key, **START button blocked** is briefly displayed. A centrifugation run can only be started via the interface. The LED in the key will be illuminated during the centrifugation run as long as the rotor is turning.

STOP
OPEN

- End centrifugation run.
The rotor decelerates with the preset run-down parameters. The right-hand LED in the button lights up until the rotor is stationary. Once the rotor is stationary the left-hand LED flashes in the button. Pressing the button twice triggers the EMERGENCY STOP.
- Unlock the lid.
The left-hand LED in the button goes out.
- Exit the parameter input and the menus.

18 Enter centrifugation parameters

Parameters can be input via the operating unit (key position "0") or via the interface (key position "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

If no key is pressed for 8 seconds long after the selection or during the input of parameters, the previous values will be shown in the display. The input of parameter then has to be executed again.

If several parameters are input, the **[START]** key must be pressed after setting the last parameter in order to accept all changes.

If parameters are changed, the program place number is displayed in parentheses (). This means that the centrifugation data in the display no longer corresponds to the centrifugation data from the program place that has been saved.

You can no discontinue entering parameters at any time by pressing the key **[OPEN/STOP]**. In this case the adjustments are not saved.

18.1 Running time

In order to set continual running, the minutes, seconds and hours have to be set to zero. Continual running is represented in the display by the following symbol, "∞".

- Press the **[TIME]** key. The parameters **t/hms** is displayed. The minutes (**m**) are shown in parentheses ⟨⟩, and can be changed.
- Use the adjusting knob **○** to set the value you want.
- Press the **[TIME]** key. The seconds (**s**) are shown in parentheses ⟨⟩ and can be changed.
- Use the adjusting knob **○** to set the value you want.
- Press the **[TIME]** key. The hours (**h**) are shown in parentheses ⟨⟩ and can be changed.
- Use the adjusting knob **○** to set the value you want.
- To apply the setting to the display, either press the **[START]** key or press the **[TIME]** key as often as is required until the centrifugation data are displayed.

18.2 Starting runtime timing

The start to runtime timing can only be set if the function, "Dual time mode" has been activated. See chapter "Activate or deactivate function "Dual time mode". This function is activated ex works.

- Press the **[TIME]** key as often as is required until **Timing begins at start** or **Timing begins at speed** is displayed.
- Select the desired setting with the turning knob **○**.
Timing begins at Start = The runtime begins to count immediately after the centrifugation run starts.
Timing begins at Speed = The runtime only begins to count after the set speed is reached.
This is indicated on the display by the symbol **⏏**, to the left next to the time.
- Press the **[TIME]** or **[START]** keys to apply the setting to the display.

18.3 Speed (RPM)

- Press the **[RPM]** key. The parameters **RPM** is displayed.
- Use the adjusting knob **○** to set the value you want.
- Press the **[RPM]** or **[START]** keys to apply the setting to the display

18.4 Relative centrifugal force (RCF) and centrifuging radius (RAD)

Relative centrifugal force (RCF) depends on the centrifuging radius (RAD). Before RCF is set, the centrifuging radius has to be set.

- Press the RCF key as often as is required until the parameters **RAD** and **RCF** are displayed and the value of the parameter, **RAD** is displayed in parentheses, $\langle \rangle$ e.g. **RAD = $\langle 146 \rangle$ RCF = 3695**. The LED is lit in the key.
- Use the adjusting knob \odot to set the centrifuging radius you want.
By changing the centrifuging radius, the value adjusts automatically to the RCF.
- Press the RCF key again. The value of the parameter, **RCF** is displayed in $\langle \rangle$ parentheses, e.g. **RAD = 146 RCF = $\langle 3695 \rangle$** .
- Press the PROG key to save the set RCF value. See the "Inputting or changing programs" chapter.
- Press the RCF or START keys to apply the setting to the display.

Only by storing (STO) the set RCF value will the RPM value resulting from that be accepted.

18.5 Start-up and run-down parameters

The set start-up and run-down parameters are displayed.

x: 1-9 = start-up stage, t = start-up time

y: 1-9 = braking stage, 0 = unbraked run-down, t = run-down time

18.5.1 Start-up stage and start-up time

It is only possible to set start-up times if these have been activated. See the chapter "Activating or deactivating start-up and run-down times".

- Press the S key. The parameter S or $\text{S}t$ is displayed.
 S = start-up stage, $\text{S}t$ = start-up time
Press the TIME key to switch between the start-up stage and start-up time.
- Set the desired stage or time with the rotary knob \odot .
- If necessary, press the S key to set the next parameter.
- To apply the setting to the display, either press the START key or press the S key as often as is required until the centrifugation data are displayed.

18.5.2 Braking stage and run-down time

No B-brake stages can be set for this device. The activation of the B-brake stages in the "Settings" menu is not possible (parameter **B-Ramp = off**). B-brake stages are similar to an exponential braking curve.
It is only possible to set run-down times if these have been activated. See the chapter "Activating or deactivating start-up and run-down times".

- Press the B key until the parameter B or $\text{B}t$ is displayed.
 B = braking stage, $\text{B}t$ = run-down time
Press the TIME key to switch between the braking stage and run-down time.
- Set the desired stage or time with the rotary knob \odot .
- If necessary, press the B key to set the next parameter.
- To apply the setting to the display, either press the START key or press the B key as often as is required until the centrifugation data are displayed.

18.5.3 Brake switch-off speed

- Press the B key as often as necessary until the parameter **N Brake** is shown.
- Use the adjusting knob \odot to set the value you want.
- Press the B or START keys to apply the setting to the display.

19 Programming

Programming can be done via the operating unit (key position "0") or via the interface (key position "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

If parameters are changed, the program place number is displayed in parentheses (). This means that the centrifugation data in the display no longer corresponds to the centrifugation data from the program place that has been saved.

19.1 Inputting or changing programs

- Set the parameters you want (see chapter "Enter centrifugation parameters").
- Press the **PROG** key as often as necessary until the **STO** parameter is shown.
- Use the adjusting knob \odot to set the program place you want.

If a "+" is shown after the program location, then this data is write-protected. In this case, the write protection has to be removed first before saving (see the chapter on "Write protection for programs").

- Press the **START** key to have the settings saved to the program place you want. Briefly, **Program store ...** is displayed as confirmation.

The previous program location data will be overwritten when the new data is saved.

If **"Protected !!"** is displayed, then the data at the program location is write-protected and it will not be saved.

19.2 Calling up programs

- Press the **PROG** key. The parameters **RCL** is displayed.
- Use the adjusting knob \odot to set the program place you want.

If a "+" is shown after the program location, then this data is write-protected.

- Press the **START** key. Briefly, **Program recall ...** is displayed as confirmation. The centrifugation data of the selected program place is shown.

19.3 Write protection for programs

The programs can be protected against unintentional modification.

The write protection can be activated/deactivated as follows when the rotor is at a standstill:

- Call up the desired program (see the chapter on "Calling up programs").
- Press the **PROG** key. The parameters **RCL** is displayed.
- Press and hold the **PROG** button for eight seconds. The parameter **STO** is displayed. After eight seconds, **Set Protection = 1-** (for example) appears in the display.
- Set \odot "+" or "-" with the adjusting knob.
 - + = program is write-protected,
 - = program is not write-protected.
- Press the **START** button to save the setting.

19.4 Program linkage

"Program linkage" can be used to link several centrifuge operations together.

A program linkage is only possible if this has been activated (parameter **Multi programs = on**; see the chapter on "Activating or deactivating program linkage").

19.4.1 Activating or deactivating program linkage

The program link can be activated/deactivated as follows when the rotor is at a standstill:

This setting can only be made via the operating unit (key position "0" or "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

You can scroll backwards through the menu by pressing the $\overline{T/^\circ C}$ key.

The operation can be cancelled at any time by pressing the $\overline{\text{OPEN/STOP}}$ key. In this case, the settings are not stored.

- Press and hold the $\overline{\text{PROG}}$ button for eight seconds.
After eight seconds, ***** Machine Menu ***** appears in the display.
- Press the $\overline{\text{PROG}}$ key as often as necessary until the **->Settings** is displayed.
- Press the $\overline{\text{START}}$ key. The **SOUND / BELL = off/on** is displayed.
- Press the $\overline{\text{PROG}}$ key as often as necessary until the **Multi programs = off/on** is displayed.
- Set **off** or **on** with the adjusting knob \odot .
off = program linkage deactivated,
on = program linkage activated.
- Press the $\overline{\text{START}}$ key to save the setting.
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the $\overline{\text{OPEN/STOP}}$ key once to exit the "Settings" menu or press the $\overline{\text{OPEN/STOP}}$ key twice to exit the "Machine Menu".

19.4.2 Linking programs or changing a program linkage

It is possible to save 25 program linkages (program locations A to Z; there is no program location J). A program linkage can consist of no more than 20 programs.

The speed adjustment from one program to the next one is always done in a program linkage with the start-up parameter of the next program.

No continuous run programs or programs with start-up and run-down times (parameters $\swarrow t$ and $\searrow t$) may be linked.

No centrifuge parameters may be modified in a program linkage. A parameter modification is only possible in the individual programs.

The $\overline{\text{TIME}}$ button can be used during the centrifuge operation to call up the total run time of the program linkage

(e. g. $\Sigma=00:05:30$) and the run time of the currently running program (e.g. **t B.02=00:01:00**).

1. Press the $\overline{\text{PROG}}$ button repeatedly until the parameter **EDIT A...Z** is displayed.
2. Use the \odot adjusting knob to set the desired program location where the program linkage is to be saved.
3. Press the $\overline{\text{START}}$ button. The program location of the program linkage and the first program of the program linkage will be shown (e.g. **EDIT B.01 = 01**).
4. Use the \odot adjusting knob to set the first program of the program linkage.
5. Press the $\overline{\text{PROG}}$ button. The next program of the program linkage will be shown (e.g. **EDIT B.02 = END**).
6. Use the \odot adjusting knob to set the next program of the program linkage.
7. Press the $\overline{\text{PROG}}$ button. The next program of the program linkage will be shown (e.g. **EDIT B.03 = END**).
8. Repeat steps 6 and 7 until all programs have been set.
9. Set **END** with the \odot adjusting knob (turn adjusting knob anti-clockwise).

For program linkages which consist of 20 programs, **END** cannot be set after the 20th program.

10. Press the $\overline{\text{START}}$ button. **STO B**, for example, is displayed.
11. Press the $\overline{\text{START}}$ button to save the program linkage.
As confirmation, **Multi program store ..** will be shown briefly.

19.4.3 Calling up program linkage

- Press the **ⓂPROG** button repeatedly until the parameter **RCL A...Z** is displayed.
- Use the adjusting knob **⊖** to set the program place you want.
- Press the **ⓂSTART** key. Briefly, **Multi program recall ...** is displayed as confirmation. The centrifugation data of the first program of the program linkage as well as the total runtime of the program linkage are displayed.

19.5 Automatic temporary storage

The program place 0 serves as temporary storage for centrifugation data of the last centrifugation run that took place. No programs can be saved at this program place.

After every start to a centrifugation run, the centrifugation data that is used for the run is automatically saved at program place "0" and can be accessed there.

20 Centrifugation

The centrifuge may only be operated in an automated system with loading and unloading by a robot. The centrifuge is controlled by the automated system via the interface. A centrifugation run can only be started via the interface.

The manufacturer of the automated system is responsible for the safety of the overall system including the centrifuge integrated in the system.

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

The centrifuge may only be controlled via the interface (key position "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

A centrifugation run can be stopped at any time by pushing the key **ⓂOPEN/STOP**.

During a centrifugation run, parameters can be selected and modified (see chapter, "Change settings during centrifugation run").

One can toggle between the RPM and RCF display at any time via the interface or with the **ⓂRPM** and **ⓂRCF** keys. When working with program links, it is not possible to switch with the **ⓂRPM** and **ⓂRCF** keys. When working with the RCF display, the centrifuging radius must be input.

If **Enter max cycles = <30000>** is displayed, first of all the maximum permitted number of running cycles specified on the rotor is entered, before the centrifugation run can be restarted (see chapter "Cycle Counter").

If the rotor has been changed, no centrifugation run takes place and this display is shown, e.g. **Rotor 4 Nmax= 4500 R=184 mm** (see chapter "Rotor Identification").

Operation errors and malfunctions will be shown (see Chapter "Faults").

20.1 Centrifugation with preselected time

- Set the centrifugation parameters or call up a program or a program linkage (see the chapters on "Enter centrifugation parameters", "Calling up programs" or "Program linkage").
- Start the centrifugation run via the interface. The LED in the **ⓂSTART** key blinks until the rotor has been imported, it is subsequently lit.
- After the time expires or if the centrifugation run is aborted via the interface or by pressing the **ⓂOPEN/STOP** key, the run-down is done with the selected run-down parameter. The run-down parameter is displayed e.g. **9.9**. The right-hand LED in the **ⓂOPEN/STOP** key is lit. Once the rotor has come to a standstill, the LED in the **ⓂSTART** key and the right LED in the **ⓂOPEN/STOP** key switch off and the left LED in the **ⓂOPEN/STOP** key lights up.

During the centrifugation run the rotor speed or the RCF value and the remaining time are displayed.

20.2 Continuous operation

- Set the minutes, seconds and hours to "0" or start a continual run program (see chapter "Enter centrifugation parameter" or "Calling up programs").
- Start the centrifugation run via the interface. The LED in the **START** key blinks until the rotor has been imported, it is subsequently lit. The time count starts from 00:00.
- End the centrifugation run via the interface or by pressing the **OPEN/STOP** key. Run-down is carried out with the selected run-down parameters. The run-down parameter is displayed e.g. 9. The right-hand LED in the **OPEN/STOP** key is lit. Once the rotor has come to a standstill, the LED in the **START** key and the right LED in the **OPEN/STOP** key switch off and the left LED in the **OPEN/STOP** key lights up.

During the centrifugation run the rotor speed or the RCF value and the elapsed time are displayed.

21 Change the settings during the centrifugation run

 During the centrifugation run centrifugation parameters can only be changed via the interface (key position "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

The runtime, speed, relative centrifugal force (RCF) as well as start-up and run-down parameters can be changed during the centrifugation run.

- Change the value of the desired parameter via the interface.
The changed setting is saved at program place "0" (see chapter "Automatic temporary storage"). The original program is not overwritten. The program place number is displayed in parentheses (). This means that the centrifugation data in the display no longer corresponds to the centrifugation data from the program place that has been saved.

22 Integral RCF

The integral RCF is a measure for the effect of sedimentation ($\int n^2 dt$). This numerical value is used to compare centrifugation runs.

22.1 Query integral RCF

 The integral RCF can be queried via the operating unit (key position "0") or via the interface (key position "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.
You can only query integral RCF if the integral RCF display is activated. See chapter "Activate or deactivate integral RCF display".
The integral RCF is not saved. After starting the next centrifugation run or after switching off the device, the integral RCF is deleted.
If the function, "**Timing begins at speed**" is selected, calculation of the integral RCF begins once the set speed has been reached.

- Press the **RCF** key as often as is necessary until the integral RCF is displayed, e.g. $\Sigma=4.8667e+05$ ($\Sigma=4.8667e+05 = 4,8667 \times 10^5 = 486670$).
- Press the **RCF** key. The centrifugation data are displayed again.
- If necessary, press the **RPM** key to switch over to the RPM display.

22.2 Activate or deactivate the integral RCF display

The integral RCF display is activated or deactivated as follows when the rotor is at rest:

This setting can only be made via the operating unit (key position "0" or "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

You can scroll backwards through the menu by pressing the **←** key.

The operation can be cancelled at any time by pressing the **↻** key. In this case, the settings are not stored.

- Keep the **Ⓜ** key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the **Ⓜ** key as often as necessary until the **->Settings** is displayed.
- Press the **Ⓢ** key. The **SOUND / BELL = off/on** is displayed.
- Press the **Ⓜ** key as often as is required until **RCF Integral = off/on** is displayed.
- Use the adjusting knob to set **Ⓢ** **off** or **on**.
off = integral RCF is deactivated,
on = integral RCF is activated.
- Press the **Ⓢ** key to save the setting.
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the **↻** key once to exit the "Settings" menu or press the **↻** key twice to exit the "Machine Menu".

23 Emergency Stop

An Emergency Stop can be executed via the operating unit (key position "0" or "LOCK 2") or via the interface (key position "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

- Press the **↻** key twice.

During a not-stop, run-out is carried out with braking level 9 (shortest run-out). Braking level **↵9** is displayed. If braking level 0 has been preselected, run-out is carried out with braking level **↵9d**. With braking level 9d, the run-out time is longer, technically speaking, than braking level 9.

24 Cycle counter

The period of use of the rotor is limited to 50000 running cycles (centrifugation runs).
 If these maximum permitted number of running cycles is reached the rotor must be immediately replaced with a new rotor for safety reasons.
 If this is ignored, it can lead to material damage and injury to personal.
 Status of the running cycles must be queried during each centrifugation run via the interface, see parameters 00635 and 00685 in interface description AH5680-01EN.

The settings can be made via the operating unit (key position "0" or "LOCK 2") or via the interface (key position "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

The centrifuge is equipped with a cycle counter, which counts the running cycles (centrifugation runs).

If a rotor is recognized the first time, the centrifugation run is canceled. After pressing any key, **Enter max cycles = <30000>** is displayed. The maximum permitted number of running cycles specified on the rotor has to be entered before the centrifugation run can be restarted (see chapter "After starting the first centrifugation run, enter the maximum permitted number of running cycles").

If the entered maximum permitted number of running cycles of the rotor is exceeded, the message *** MAX CYCLES PASSED *** is displayed after every start of a centrifugation run and the corresponding bit in the parameters 00635 and 00685 is set. The centrifugation run has to be restarted.

If *** MAX CYCLES PASSED *** is displayed and the corresponding bit in the parameters 00635 and 00685 is set, the rotor must be immediately replaced with a new rotor for safety reasons.

After replacing the rotor, the cycle counter has to be reset to "0" (see chapter "Reset cycle counter to "0" and enter the maximum permitted number of cycles").

24.1 After starting the first centrifugation run, enter the maximum permitted number of running cycles

- **Enter max cycles = <30000>** is displayed.
Use the knob \odot to set the maximum permitted number of running cycles specified on the rotor.
- Press the \odot START key to save the setting.
Briefly, **Store max cycles ...** is displayed as confirmation.

24.2 Reset cycle counter to "0" and enter the maximum permitted number of cycles

The cycle counter may only be reset to "0" if the rotor has been replaced with a new rotor first.

If the rotor is at standstill, this can be set as follows:

You can scroll backwards through the menu by pressing the \odot T°C key.
 The operation can be cancelled at any time by pressing the \odot OPEN / STOP key. In this case, the settings are not stored.

- Keep the \odot PROG key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the \odot PROG key as often as necessary until the **->Operating Time** is displayed.
- Press the \odot START key. The external operating hours are displayed, e.g. **OP Time ext = 0h25m**.
- Press the \odot PROG key as often as is required until the running cycles are displayed, e.g. **Cycles = 80001 of 80000**.
- Press the \odot RCF key. The number of running cycles is displayed in parentheses \acute{a} ñ , e.g. **Cycles = <80001> of 80000**.
- Turn the knob \odot to the left to reset the number of running cycles to "0".

If the running cycles are not reset to "0", press the \odot START key, which causes **Max cycles (= actual cycles)** to be displayed and the setting is not saved.

- Press the \odot RCF key. The maximum permitted number of running cycles is displayed in parentheses $\langle \rangle$, e.g. **Cycles = 0 of <80000>**.
- Use the knob \odot to set the maximum permitted number of running cycles specified on the rotor.

- Press the **START** key to save the settings.
In confirmation, **Store cycles ...** is displayed briefly. Next the running cycles are displayed, e.g. **Cycles = 0 of 80000**.
- Press the **OPEN/STOP** key twice to exit the "Operating Time" menu or press the **OPEN/STOP** key three times to exit the "Machine Menu".

24.3 Activating the cycle counter

The cycle counter is activated ex works.

You can scroll backwards through the menu by pressing the **TEMP** key.

The operation can be cancelled at any time by pressing the **OPEN/STOP** key. In this case, the settings are not stored.

The cycle counter can be activated/deactivated as follows when the rotor is at a standstill:

- Keep the **PROG** key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the **PROG** key as often as is required until the **->Operating Time** is displayed.
- Press the **START** key. The external operating hours are displayed, e.g. **OP Time ext = 0h25m**.
- Press the **PROG** key as often as is required until **Cycles = disabled** is displayed with the cycle counter deactivated.
If running cycles are displayed, e.g. **Cycles = 5120 of 80000**, the cycle counter is already activated.
- Press the **RCP** key as often as is required until the maximum permitted number of running cycles is displayed in parentheses $\langle \rangle$, e.g. **Cycles = 0 of (0)**.
- Use the knob \odot to set the maximum permitted number of running cycles specified on the rotor.
- Press the **START** key to save the setting.
In confirmation, **Store cycles ...** is displayed briefly. Next the running cycles are displayed, e.g. **Cycles = 0 of 80000**.
- Press the **OPEN/STOP** key twice to exit the "Operating Time" menu or press the **OPEN/STOP** key three times to exit the "Machine Menu".

25 Activate or deactivate the function, "Dual time mode"

The function "Dual time mode" can be activated or deactivated as follows when the rotor is at a standstill:

This setting can only be made via the operating unit (key position "0" or "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

You can scroll backwards through the menu by pressing the $\overline{\text{T/C}}$ key.

The operation can be cancelled at any time by pressing the $\overline{\text{OPEN/STOP}}$ key. In this case, the settings are not stored.

- Keep the $\overline{\text{PROG}}$ key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the $\overline{\text{PROG}}$ key as often as necessary until the **->Settings** is displayed.
- Press the $\overline{\text{START}}$ key. The **SOUND / BELL = off/on** is displayed.
- Press the $\overline{\text{PROG}}$ key as often as is required until **Dual time mode enabled/disabled** is displayed.
- Use the knob \odot to set **enabled** or **disabled**.
disabled = The function is deactivated,
enabled = The function is activated.
- Press the $\overline{\text{START}}$ key to save the setting.
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the $\overline{\text{OPEN/STOP}}$ key once to exit the "Settings" menu or press the $\overline{\text{OPEN/STOP}}$ key twice to exit the "Machine Menu".

26 Activating or deactivating start-up and run-down times

The start-up and run-down times can be activated or deactivated as follows with the rotor at a standstill:

This setting can only be made via the operating unit (key position "0" or "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

You can scroll backwards through the menu by pressing the $\overline{\text{T/C}}$ key.

The operation can be cancelled at any time by pressing the $\overline{\text{OPEN/STOP}}$ key. In this case, the settings are not stored.

- Keep the $\overline{\text{PROG}}$ key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the $\overline{\text{PROG}}$ key as often as necessary until the **->Settings** is displayed.
- Press the $\overline{\text{START}}$ key. The **SOUND / BELL = off/on** is displayed.
- Press the $\overline{\text{PROG}}$ key as often as is required until **Ramp Unit = Steps / Steps / Time** is displayed.
- Use the adjusting knob to set \odot **Steps** or **Steps / Time**.
Steps = start-up and run-down times deactivated,
Steps / Time = start-up and run-down times activated.
- Press the $\overline{\text{START}}$ key to save the setting.
Briefly, **Store Settings ..** is displayed as confirmation, followed by **-> Settings**.
- Press the $\overline{\text{OPEN/STOP}}$ key once to exit the "Settings" menu or press the $\overline{\text{OPEN/STOP}}$ key twice to exit the "Machine Menu".

27 Acoustic signal

The acoustic signal sounds:

- After a disturbance occurs, in 2 second intervals.
- After completion of a centrifugation run and rotor standstill in 30 second intervals.

The acoustic signal is stopped by opening the lid or pressing any key.

The acoustical signal can be activated or deactivated as follows when the rotor is at a standstill:

The acoustical signal is deactivated ex works.

This setting can only be made via the operating unit (key position "0" or "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

You can scroll backwards through the menu by pressing the $\overline{\text{TFC}}$ key.

The operation can be cancelled at any time by pressing the $\overline{\text{OPEN/STOP}}$ key. In this case, the settings are not stored.

- Keep the $\overline{\text{PROG}}$ key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the $\overline{\text{PROG}}$ key as often as necessary until the **->Settings** is displayed.
- Press the $\overline{\text{START}}$ key. The **SOUND / BELL = off/on** is displayed.
SOUND / BELL : Signal after the centrifugation run ends.
- Use the adjusting knob \odot to set **off** or **on**.
- Press the $\overline{\text{PROG}}$ key. The **SOUND / BELL error = off/on** is displayed.
SOUND / BELL error : Signal after a disturbance occurs.
- Use the adjusting knob \odot to set **off** or **on**.
- Press the $\overline{\text{START}}$ key to save the setting.
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the $\overline{\text{OPEN/STOP}}$ key once to exit the "Settings" menu or press the $\overline{\text{OPEN/STOP}}$ key twice to exit the "Machine Menu".

28 Displayed centrifugation data after switch-on.

After switch-on the centrifugation data from program 1, or from the last program that was used, is displayed.

If the rotor is at standstill, this can be set as follows:

This setting can only be made via the operating unit (key position "0" or "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.

You can scroll backwards through the menu by pressing the $\overline{\text{TFC}}$ key.

The operation can be cancelled at any time by pressing the $\overline{\text{OPEN/STOP}}$ key. In this case, the settings are not stored.

- Keep the $\overline{\text{PROG}}$ key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the $\overline{\text{PROG}}$ key as often as necessary until the **->Settings** is displayed.
- Press the $\overline{\text{START}}$ key. The **SOUND / BELL = off/on** is displayed.
- Press the $\overline{\text{PROG}}$ key as often as necessary until **Start program = Last/First** is displayed.
- Use the adjusting knob \odot to set **Last** or **First**.
Last = the program used last, First = program 1.
- Press the $\overline{\text{START}}$ key to save the setting.
Briefly, **Store Settings ...** is displayed as confirmation, followed by **-> Settings**.
- Press the $\overline{\text{OPEN/STOP}}$ key once to exit the "Settings" menu or press the $\overline{\text{OPEN/STOP}}$ key twice to exit the "Machine Menu".

29 Centrifuge address

Ex works, the address is] = 29. Address is set.

30 Query the operating hours, centrifugation runs and cycle counter

The operating hours are divided up into internal and external operating hours.

Internal operating hours: Total time the device was switched on.

External operating hours: Total time of previous centrifugation runs.

The query is performed as follows when the rotor is at a standstill:

 The query can be made via the operating unit (key position "0" or "LOCK 2") or via the interface (key position "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.
You can scroll backwards through the menu by pressing the $\overline{T/C}$ key.

- Keep the $\overline{\text{PROG}}$ key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the $\overline{\text{PROG}}$ key as often as necessary until the **-> Operating Time** is displayed.
- Press the $\overline{\text{START}}$ key. The external operating hours are displayed, e.g. **OP Time ext = 0h25m**.
- Press the $\overline{\text{PROG}}$ key. The internal operating hours are displayed, e.g. **OP Time int = 1h36m**.
- Press the $\overline{\text{PROG}}$ key. The number of all centrifugation runs is displayed, e.g. **Number of Starts = 10**.
- Press the $\overline{\text{PROG}}$ key. The number of running cycles (centrifugation runs) of the current rotor code since the last time the cycle counter was reset to "0" and the permitted number of running cycles are displayed, e.g. **CYCLES = 5120 of 30000**.
- Press the $\overline{\text{PROG}}$ key. The number of all running cycles (centrifugation runs) of the current rotor code is displayed, e.g. **Rotor cycles total = 37490**. This value cannot be set.
- Press the $\overline{\text{OPEN/STOP}}$ key twice to exit the "Operating Time" menu or press the $\overline{\text{OPEN/STOP}}$ key three times to exit the "Machine Menu".

31 Query the system information

The following system information can be queried:

- Centrifuge model
- Mains voltage
- Rotor information
- Program version of the centrifuge
- Program version of the frequency converter

The system information can be queried as follows when the rotor is at a standstill:

 The query can be made via the operating unit (key position "0" or "LOCK 2") or via the interface (key position "LOCK 2"). The functions of the various key positions are described in the "Key-operated switch" chapter.
You can scroll backwards through the menu by pressing the $\overline{T/C}$ key.

- Keep the $\overline{\text{PROG}}$ key held down for 8 seconds.
After 8 seconds, ***** Machine Menu ***** appears in the display.
- Press the $\overline{\text{PROG}}$ key as often as necessary until the **->Info** is displayed.
- Press the $\overline{\text{START}}$ key. The centrifuge model is displayed.
- Press the $\overline{\text{PROG}}$ key. The mains voltage is displayed, e.g. **Mains Voltage : 230 V**.
- Press the $\overline{\text{PROG}}$ key. The rotor code (Rotor), the maximum speed of the rotor (Nmax) and a centrifuging radius (R) of the last rotor detected by the rotor detector are displayed. e.g.
Rotor 4* : Nmax = 4500 R=184.
The last recognized rotor is marked with a star (*). Use the adjusting knob \odot to have the information of the rotors permitted in the centrifuge displayed.

 The required centrifuging radius must be set according to the used accessories; see chapter "Entering centrifugation parameters".

- Press the $\overline{\text{PROG}}$ key. The program version of the centrifuge is displayed, e.g. **SW version = V01.00**.
- Press the $\overline{\text{PROG}}$ key. The program version of the frequency converter is displayed e.g. **FC-SW-Version = 4**.
- Press the $\overline{\text{OPEN/STOP}}$ key twice to exit the "Info" menu or press the $\overline{\text{OPEN/STOP}}$ key three times to exit the "Machine Menu".

32 Immediate display of centrifugation data after switch-on.


- Switch on the power supply. (Switch setting I).
- The first time the display changes optically, (inverse display) press any key and keep it held down. The centrifugation data is displayed immediately.

When handing over the automated system to the user himself has the key switch in key-position "LOCK 2" are and the key must be withdrawn

33 Key-operated switch

 All the key-operated switches have the same keys.
 Only authorized skilled personnel may operate the device in the "TEACH" key setting.
 When handing over the automated system to the user the key-operated switch has to be in the "LOCK 2" setting and the key must be removed.
 The key is to be kept in a safe place so that it is protected from unauthorized access.

The following functions can be set with the key-operated switch when the rotor is at a standstill:

| Key setting  | Function |
|--|---|
| 0 | No status display. The centrifuge can be operated via the operating unit. No centrifugation run can be started via the operating unit, however. The centrifuge cannot be controlled via the interface. |
| LOCK 2 | LOCK 2 is displayed in the "↗" field. The centrifuge can only be controlled via the interface. Only settings in the "Machine Menu" can be performed via the operating unit. |
| TEACH | Display e.g. <div style="display: flex; justify-content: space-around; align-items: flex-start;"> <div style="text-align: center;"> <p>Teach</p>  </div> <div style="text-align: center;"> <p>Open</p>  </div> <div style="text-align: center;"> <p>=Close</p>  </div> </div> <p>= : Position of the hatch  : Open the hatch (keep key pressed)  : Close the hatch (keep key pressed)  : Set rotor position 1. See the chapter "Setting rotor position 1".</p> |

34 Setting rotor position 1



Only authorized skilled personnel may set rotor position 1.

If the rotor was changed and if this has another rotor code than the previously used rotor, a rotor detection procedure must be carried out first before rotor position 1 is set (see chapter "Rotor recognition").

Every time the rotor is installed, rotor position 1 must be reset afterwards. Rotor position 1 is the position where place 1 of the rotor must be located for loading and unloading. Rotor position 1 can be set as follows when the rotor is at a standstill:


- Look through the window in the lid to make sure that the rotor is at a standstill.
- Turn the key to the "TEACH" setting. For example, **Teach Open =Close** might be displayed.
- Keep the  key pressed until the hatch is completely open. **Teach =Open Close** is displayed.
- Press the  key. The saved rotor position 1 is displayed, e.g. **Teach Pos 1 : 1812 <Start=ok>**.
- Turn the rotor slowly in an anti-clockwise direction until you hear an acoustic signal. Then continue turning until position 1 of the rotor is in the desired loading and unloading position.. Hold the rotor in this position.
- Press the  button to save the setting.

As a confirmation, **store Pos 1 ...** is briefly displayed, and then e.g. **Goto 1/24 2826 0 2826 s**.

If "**no zero pulse !**" is briefly displayed, the rotor was turned to the loading/unloading position before the acoustic signal sounded. In this case, repeat the rotor position 1 setting.

If **Close the Lid !** is displayed briefly, the lid is open. In this case, close the lid.

- To test the positioning, set the desired rotor place with the rotary knob , and then press the  key. The rotor continues to turn until the set rotor place is in the loading/unloading position.
- Press the  key to exit the menu. **Teach =Open Close** is displayed.
- Keep the  key pressed until the hatch is completely closed. **Teach Open =Close** is displayed.
- Turn the key back to the "LOCK 2" position.

35 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left(\frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge container. For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

36 Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm³

When centrifuging with maximum revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.

The speed must be reduced for materials or mixtures of materials with a higher density.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density [kg/dm}^3]}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, density 1.6 kg/dm³

$$n_{\text{red}} = \sqrt{\frac{1.2 \text{ kg/dm}^3}{1.6 \text{ kg/dm}^3}} \times 4000 \text{ RPM} = 3464 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

37 Rotor recognition

Rotor recognition is carried out after each start of the centrifugation run.

If the rotor has been changed, the centrifugation run is canceled after rotor recognition. The rotor code (Rotor), the maximum speed of the rotor (Nmax) and a centrifuging radius (R) of the newly detected rotor are displayed, e.g.

Rotor 4 Nmax= 4500 R=184 mm.

If the maximum speed of the rotor being used is less than the set speed, the speed is limited to the rotor's maximum speed. In this case, the program place number is displayed in parentheses ().

- Start the centrifugation run via the interface.

If the cycle counter is activated, the number of rotor cycles (centrifugation runs) of the current rotor code is displayed briefly when the lid is opened, e.g. **CYCLES 5120 of 30000** (see chapter "Cycle Counter").

38 Maintenance and servicing

The device can be contaminated.

Pull the mains plug before cleaning.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Centrifuges, rotors and accessories must not be cleaned in rinsing machines.
- They may only be cleaned by hand and disinfected with liquids.
- The water temperature must be between 20 – 25°C.
- Only detergents/disinfectants may be used which:
 - have a pH between 5 - 8
 - do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

38.1 Centrifuge (housing, lid and centrifuging chamber)

38.1.1 Surface cleaning and care

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- Ingredients of suitable detergents: soap, anionic tensides, non-ionic tensides.
- After using detergents, remove the detergent residue by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- The centrifuging chamber is to be checked for damage once a year.

If damage is found which is relevant to safety, the centrifuge may no longer be put into operation. In this case, notify Customer Service of the automated system.

38.1.2 Surface disinfection

- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately.
- Ingredients of suitable disinfectants: ethanol, n-propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
- The surfaces must be dried immediately after disinfecting.

38.1.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for removing radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants: anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
- The surfaces must be dried directly after removing the radioactive contaminants.

38.2 Rotor and accessories

38.2.1 Cleaning and care

- To avoid corrosion and changes to the materials, the rotor and accessories have to be cleaned regularly with soap or a mild cleaning agent and a moist cloth. Cleaning is recommended at least once a week. Contaminants must be removed immediately.
- Ingredients of suitable detergents: soap, anionic tensides, non-ionic tensides.
- After using detergents, remove detergent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotor and accessories have to be dried immediately after cleaning.
- Check the rotor and accessories weekly for wear and corrosion damage.

The rotor and accessories must no longer be used if they show signs of wear or corrosion.

- Check the firm seating of the rotor on a weekly basis.

38.2.2 Disinfection

- If infectious material should get on the rotor or accessories, they must be appropriately disinfected.
- Ingredients of suitable disinfectants: ethanol, n-propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove disinfectant residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotor and accessories must be dried directly after disinfection.

38.2.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for the removal of radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants: anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove agent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotor and accessories must be dried directly after removing the radioactive contaminants.

38.2.4 Rotors and accessories with limited service lives

The period of use of the rotor is limited to 80000 running cycles (centrifugation runs). The maximum permissible number of run cycles can be seen on the rotor.

For safety reasons, the rotor may no longer be used when the maximum allowed number of running cycles (marked on it) has been reached.

The device is equipped with a cycle counter which counts the running cycles (centrifugation runs). For a description, see the "Cycle counter" chapter.

38.3 Autoclaving

The rotor and accessories must not be autoclaved.

38.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts, glass splinters and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splinters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

39 Faults

As long as the centrifuge is in operation, the parameter 00634 must be queried via the interface continuously whether a fault has occurred. In case of failure, the error number is deposited in parameter 00634, see interface description AH5680-01EN.

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

Perform a MAINS RESET:

- Switch off the mains switch (switch position "0").
- Wait at least 10 seconds and then switch on the mains switch again (switch position "1").


| Message / fault | | Cause | Remedy |
|------------------|-----------------------------------|--|---|
| No display | | No voltage. Overvoltage protection tripped out. | <ul style="list-style-type: none"> - Check supply voltage. - Mains switch ON. |
| TACHO - ERROR | 1, 2, 96 | Faulty speedometer. Motor, electronics defective. | <ul style="list-style-type: none"> - Open the cover. - Switch off the mains switch (switch position "0"). - Wait at least 10 seconds. - Turn the rotor vigorously by hand. - Switch on the mains switch again (switch position "1"). The rotor must turn during switch-on. |
| IMBALANCE | 3 * | The rotor is unevenly loaded. | <ul style="list-style-type: none"> - Open the lid or hatch. - Check the loading of the rotor, see chapter "Loading the rotor". - Repeat the centrifugation run. |
| CONTROL - ERROR | 4.1 - 4.5, 6 | Error in lid locking or lid closure. | <ul style="list-style-type: none"> - Perform a MAINS RESET. |
| CONTROL - ERROR | 4.6 - 4.9 | Error / defect hatch control | |
| N > MAX | 5.0, 5.1 | Rotation too fast | |
| N < MIN | 13 | Rotation too slow | |
| ROTORCODE | 10.1, 10.3 | Incorrect rotor coding | |
| MAINS INTERRUPT | 11 * | Power failure during the centrifugation run. (The centrifugation run was not finished.) | <ul style="list-style-type: none"> - Open lid. - Push button. - Repeat the centrifugation run if necessary. |
| VERSION-ERROR | 12 | Mismatch between electronic components Error / defect electronics | <ul style="list-style-type: none"> - Perform a MAINS RESET. |
| CONTROL-ERROR | 22, 25.1-25.4 | Error / defect electronics | |
| CRC ERROR | 27.1 | Error / defect electronics | |
| SER I/O - ERROR | 31, 34, 36 | Error / defect electronics | |
| POS-ERROR | 40-49 | Error / defect positioning or hatch drive | |
| LOCK_ERROR | 57.0, 57.1 | Error / defect key-operated switch | |
| FU / CCI - ERROR | 60, 61.2-61.20, 61.128-61,131, 62 | Error / defect electronics / motor | |
| FU / CCI - ERROR | 61.1 | Supply voltage is too low Error / defect electronics / motor | <ul style="list-style-type: none"> - Check supply voltage. - Perform a MAINS RESET. |



* Error number does not appear in the display.

| Message / fault | | Cause | Remedy |
|---------------------------------|---------|---|---|
| SENSOR-ERROR | 90 | Error / defect electronics | – Perform a MAINS RESET. |
| SENSOR-ERROR | 91 - 93 | Error / defect unbalance sensor | |
| ° C * - ERROR | 97, 98 | Error / defect electronics | |
| NO ROTOR OR ROTORCODE ERROR | --- | No rotor installed. | – Open lid. |
| N > ROTOR MAX | --- | Faulty speedometer. | – Install rotor. |
| | | Speed in the selected program greater than the maximum speed of the rotor. | – Check the set speed. – Reduce the set speed |
| N > ROTOR MAX in Prog : e.g. 3 | --- | The rotor has been changed. The maximum speed of the installed rotor is higher than that of the rotor that was previously used and this has not yet been detected by rotor recognition. | – Set a speed up to the maximum speed of the rotor that was previously used. Press the START key to run rotor recognition. See chapter, "Rotor Recognition". |
| | | There is a program in the displayed program location whose speed is higher than the maximum speed of the rotor. | – Check the set speed. – Reduce the set speed |
| | | The rotor has been changed. The maximum speed of the installed rotor is higher than that of the rotor that was previously used and this has not yet been detected by rotor recognition. | – Set a speed up to the maximum speed of the rotor that was previously used. Press the START key to run rotor recognition. See chapter, "Rotor Recognition". |
| Runtime 00:00 in Prog : e.g. 3 | --- | There is a continuous operation program in the displayed program location. | – Replace the continuous operation program in the program linkage with a program with preset time. |
| Empty Program | --- | There is no program linkage saved on the displayed program location. | – Call up a program linkage. |
| Ramp Unit Time in Prog : e.g. 3 | --- | A program with start-up and/or run-down time is located at the displayed program place. | – Replace the program in the program link with a program with start-up and braking stage. |
| Acc time > Run time | --- | The set start-up time is longer than the run time. | – Set a start-up time which is shorter than the run time. – Perform a MAINS RESET. |
| FC INIT ERROR | --- | Error / defect electronics | |
| FC VERSION ERROR | --- | Error / defect electronics | |
| FATAL EEPROM ERROR | 1 - 5 | Error / defect electronics | |
| WATCHDOG RESET | | Error / defect electronics | – |

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1) Sample temp. rise during maximum speed and 1 hour running time

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2)  Centrifugation of flammable materials is not allowed; application is therefore restricted to **non-flammable** solvents.