

## **ENVIRONMENT/WATER ANALYSIS**

Use of HettCube cooled incubators to determine the biochemical oxygen demand

The decomposition of organic material by micro-organisms is an important biochemical oxidation process in nature. The quantity of oxygen required is dependent on the mass of the substance decomposed.

# The five-day BOD (BOD<sub>5</sub>) as an indicator in the analysis of wastewater

The relationship between oxygen demand and the quantity of organic substance decomposed is used in the analysis of wastewater. The BOD5 of a wastewater sample is determined and the higher its value the greater the contamination of water with substances that can be biodegraded. The five-day BOD is used in most cases. This is the quantity of oxygen that is consumed by micro-organisms upon aerobic decomposition of organic substances over a period of 5 days under controlled conditions.

## Methods of determination

Different methods are used to establish the biochemical oxygen demand. They include self-monitoring methods and the methods given in international standards such as ISO 5815-1 and ISO 5815-2.

#### The BOD<sub>5</sub> is determined in the field for various purposes, for example

- Process monitoring and design within wastewater treatment plants
- Monitoring of the organic composition of wastewater discharges, e.g. from cider press houses
- Testing of surface waters by the environmental authorities
- Calculation of the fees to be paid by industrial companies and operators of wastewater treatment plants

#### Incubation conditions

Temperature	Duration
20 ± 1 °C	5 days

### Advantages of HettCube incubators

- Maximal validated usable space on a small footprint
- 4.3 inch touch display for intuitve operation
- Very homogeneous and stable temperature, as well as precise temperature control
- True "one-hand-operation" and flexible positioning of the shelves
- Minimal energy consumption of < 0.06 kW/h at 37 °C</li>



Fig. 1: Final clarifier of a sewage clarification plant 1)



Fig. 2: Karlsruher bottles, preferably used when testing according to the international standard ISO 5815-1 (dilution method).<sup>2)</sup>

Image by courtesy of the municipality of Tuttlingen
Image by courtesy of the Windaus Labortechnik
GmbH & Co. KG

## Hettich solution

Model	Cat. No.
HettCube 200	62000
HettCube 400	64000
HettCube 600	66000

Model	
without IVD	Cat. No.
HettCube 200	62001
HettCube 400	64001
HettCube 600	66001