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## HettCube

60; 120

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HettCube

60; 120

Translation of the original operating manual

Read the instructions prior to performing any task!

© 2023

Andreas Hettich GmbH & Co. KG

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D-78532 Tuttlingen, Germany

Telephone: +49 (0)7461 705-0


Fax: +49 (0)7461 705-1125

Email: [info@hettichlab.com](mailto:info@hettichlab.com), [service@hettichlab.com](mailto:service@hettichlab.com)

Internet: [www.hettichlab.com](http://www.hettichlab.com)

## Table of contents

<b>1</b>	<b>Introduction.</b>	<b>7</b>
1.1	Symbols.	7
1.2	Important labels on the packaging.	8
1.3	Personal protective equipment.	9
1.4	Personnel qualification.	9
1.5	Intended use.	9
1.6	Replacement parts/consumable material.	10
1.7	Scope of delivery.	10
1.8	Return shipments.	11
<b>2</b>	<b>Safety information.</b>	<b>12</b>
<b>3</b>	<b>Device description.</b>	<b>14</b>
3.1	Standard device.	14
3.2	Accessories.	15
<b>4</b>	<b>Transport and storage.</b>	<b>16</b>
<b>5</b>	<b>Commissioning.</b>	<b>18</b>
5.1	Unpacking the incubator.	18
5.2	Setting up, connecting and switching on the incubator.	19
5.2.1	Installing the incubator.	19
5.2.2	Incubator connection.	21
5.3	Inserting and removing the slide-in modules.	23
5.3.1	Standard slide-in modules.	23
5.4	Install glass door.	24
5.5	Changing the door swing.	24
5.5.1	Door.	24
5.5.2	Glass door.	25
<b>6</b>	<b>Operation.</b>	<b>26</b>
6.1	Operating elements.	27
6.2	Loading.	27
6.3	Door locking mechanism.	28
6.4	Standard check before every use.	28
6.5	Switching on the incubator.	28
6.6	Unit operation.	28
6.6.1	External memory (USB flash drive).	29
6.6.2	First boot.	30
6.6.3	🏠 Main screen.	30
6.6.4	🏠 Information panel.	30
6.6.5	Meanings of symbols.	33

6.6.6	Upper expandable and configurable menu. . . . .	36
6.6.7	Alarm bar. . . . .	36
6.7	Quick Program. . . . .	37
6.8	 Programs. . . . .	39
6.8.1	Creating / editing a program. . . . .	39
6.8.2	Segment editing. . . . .	40
6.8.3	Protection class. . . . .	44
6.8.4	Priority. . . . .	44
6.8.5	Repetition / Loop. . . . .	45
6.9	Starting the program. . . . .	45
6.9.1	The first option. . . . .	45
6.9.2	The second option. . . . .	46
6.10	Quick Change of parameters. . . . .	47
6.10.1	Quick Change of set temperature. . . . .	47
6.10.2	Quick Change of set time . . . . .	48
6.10.3	Quick change of fan speed and air flap position (degree of opening). . . . .	49
6.11	 Statistics. . . . .	50
6.12	 Data record. . . . .	51
6.13	 Event log. . . . .	52
6.14	 Info. . . . .	54
6.15	 User settings panel. . . . .	55
6.15.1	Unlocking the touchscreen. . . . .	56
6.16	 Time. . . . .	56
6.17	 Alarms. . . . .	57
6.17.1	STM function. . . . .	58
6.17.2	Mute option.. . . . .	60
6.18	 Network. . . . .	61
6.19	 Correction. . . . .	62
6.20	INTERFACE. . . . .	63
6.21	Temp. protection. . . . .	63
<b>7</b>	<b>Cleaning, disinfection and maintenance. . . . .</b>	<b>65</b>
7.1	Cleaning. . . . .	66
7.1.1	Cleaning the touchscreen. . . . .	67
7.2	Disinfection. . . . .	68
7.3	Removal of radioactive contaminants. . . . .	69
7.4	Autoclaving. . . . .	69
7.5	Maintenance. . . . .	69

<b>8</b>	<b>Troubleshooting. ....</b>	<b>71</b>
8.1	Troubleshooting. ....	71
8.2	Possible defects. ....	71
<b>9</b>	<b>Technical data. ....</b>	<b>73</b>
9.1	Incubator unit 60, 120. ....	73
9.2	Type plate. ....	76
9.3	Dimensions. ....	76
<b>10</b>	<b>Disposal. ....</b>	<b>78</b>
<b>11</b>	<b>Free and Open Source Software. ....</b>	<b>79</b>
<b>12</b>	<b>Index. ....</b>	<b>82</b>
<b>13</b>	<b>Appendix. ....</b>	<b>83</b>
A	List of standards. ....	84

# 1 Introduction

## 1.1 Symbols

### Signal words

Signal word	Meaning
<b>DANGER</b>	This combination of symbol and signal word indicates an immediate dangerous situation that will result in death or serious injury if it is not avoided.
<b>WARNING</b>	This combination of symbol and signal word indicates a possible dangerous situation that can result in death or serious injury if it is not avoided.
<b>CAUTION</b>	This combination of symbol and signal word indicates a possible dangerous situation that can result in minor injury if it is not avoided.
<b>NOTICE</b>	This combination of symbol and signal word indicates a possible dangerous situation that can result in material and environmental damage if it is not avoided.

### Warning categories

Warning signs	Type of danger
	Warning – danger zone.
	Warning – biological hazard.
	Warning – hand injuries.
	Warning – high-voltage.
	Warning – falling objects.

General symbols

- ➔ This listing symbol denotes descriptions of tasks that you must perform.
- This dot is for denoting lists.

Cross references are indicated as follows: ➔ Chapter 1.1 'Symbols' on page 7

Graphic symbols on the incubator

 *If warning signs/symbols on the incubator are no longer recognisable, the operator must replace them immediately.*

*The figures below indicate the positions of the warning signs and symbols found on the incubator.*

Rating plate

Rating plate displaying details of the unit's technical data.

1.2 Important labels on the packaging

 This way up.  
Shows the correct upright position of the package.

 Fragile, handle with care.  
Refers to a medical device that can break or be damaged if handled carelessly.

 keep dry  
The shipping packaging must be kept away from rain and kept in a dry environment.

 Stack limit.  
Maximum number of identical packages which may be stacked on the bottommost package, whereby "n" stands for the number of permissible packages. The bottommost package is not included in "n".

### 1.3 Personal protective equipment

The assessment of the personal protective equipment must take place on site according to the actual risks which are posed by the used substances, the procedures and the environmental conditions.

The operating company must create a hazard analysis according to the corresponding standards and directives and enact work instructions which ensure the safe handling of the incubator and accessories.

### 1.4 Personnel qualification

Repairs may only be performed by a person authorized by the manufacturer.

*Interventions and modifications to incubators by persons not authorized by the company Andreas Hettich GmbH & Co. KG will be at one's own risk and result in the voiding of all warranty claims, as well as the voiding of liability claims against Andreas Hettich GmbH & Co. KG.*

#### **Service technician**

A service technician is someone who was trained and authorized by Andreas Hettich GmbH & Co. KG to perform service work on the incubator.

#### **User**

The users intended for the commissioning, operation and use of the incubator include persons who work in laboratories who have read and understood the operating instructions.

The operating company must ensure that the incubator is only used as intended.

### 1.5 Intended use

#### **Intended use**

- This incubator, which is used in particular for the cultivation of micro-organisms (e.g. bacteria and fungi), is a laboratory device in the sense of the Low Voltage Directive 2014/35/EU.
- Temperature control and incubation tasks for materials and samples with equally high requirements are also possible.
- This is the sole intended use of the incubator.
- Intended use also includes observing all instructions in the Operating Manual and compliance with the required inspection and maintenance work.

#### **Non-intended use**

Any other use or use beyond the intended use is a non-intended use. Andreas Hettich GmbH & Co. KG will not be liable for any resulting damage.

- The cultivation of human sample material within the temperature range of the incubator for diagnostic purposes within the meaning of IVD Directive 98/79/EC within the scope of medical examinations of people with the aim of ruling out or specifying a pathological condition.
- The cultivation/temperature control of human material, but also of other substances or objects within the temperature range of the incubator in connection with therapeutic actions (within the meaning of the Medical Device Directive 93/42/EEC) on humans.
- Operation in an explosive, biologically, chemically or radioactively contaminated atmosphere.
- Loading the incubator with combustible or explosive materials or with materials that chemically react with one another in a highly exothermic manner.
- The presence of persons or animals inside the incubator.
- Using the incubator or its accessories as a climbing aid.
- Operation outdoors.
- Use on moving platforms such as ships, rolling stock or aircraft.
- Operation of devices in interior areas that generate heat or vibrations.

## 1.6 Replacement parts/consumable material

Only original spare parts and original accessories approved by the manufacturer may be used.

### WARNING

If no original spare parts or original accessories are used, the warranty as well as liability claims against the company Andreas Hettich GmbH & Co. KG are void.

## 1.7 Scope of delivery

- 1 mains cable
- 1 incubator
- 1 Operating Manual
- 2 keys
- 1 plug for the access port on the side wall
- 2 standard shelves
- 4 rails
- 1 inspection certificate

## 1.8 Return shipments

If it is necessary to send the incubator or its accessories back to Hettich,

- request the necessary returns documentation from Hettich or your responsible Hettich sales partner. This includes the returns number and a form for certifying the harmlessness of the article to be sent back.
- Based on the legal regulations, and primarily to protect the employees of any Hettich sales partner and Hettich, thoroughly clean and disinfect the incubator and its accessories. Then certify its harmlessness by signing the appropriate certificate.
- In the event of a missing certificate of harmlessness, Hettich reserves the right to invoice the customer for the costs of cleaning and disinfection.
- The incubator should only be shipped in the original packaging.

## 2 Safety information

- The Operating Manual must be read and observed before commissioning the incubator. Only persons who have read and understood the Operating Manual may operate the incubator.
- In addition to the Operating Manual and the binding accident prevention regulations, the recognised technical rules of sound and safe work practice must also be observed. The Operating Manual shall be supplemented with further instructions for accident prevention and environmental protection based on the existing national regulations in the country of use.
- The incubator has been built to the state-of-the-art engineering standards and is safe to operate. However, it may pose a risk to the user or third parties if it is not operated by personnel familiar with the incubator or if it is used improperly or for non-intended purposes.
- Unit operators must be thoroughly informed of possible health hazards that can arise from the sample material used and, if necessary, take appropriate measures to rule out such hazards.
- The incubator may only be operated if it is set up correctly. If it is incorrectly installed, hazards to persons and sample material cannot be ruled out.
- Do not use the floor of the unit interior, its shelves, drawers or door as something to stand on, or as a support.
- The rate of drying of the culture media during incubation is influenced by various factors, such as
  - Ambient conditions (e.g. relative humidity of the surrounding environment)
  - Handling and working methods (e.g. upstream processes, duration and frequency of times when the unit door is open).
  - Culture or incubation conditions (e.g. incubation time and temperature)
  - Set fan speed for the forced convection variant
  - Air flap position. These influences should be taken into account, especially with longer incubation durations.
- Sample material should not be placed outside the defined usable space. The specified temperature data apply to the defined usable space.
- There is a risk of an uncontrolled increase in the utilized space temperature if equipment with a heat load is operated in the uncooled incubator or if heat-generating reactions or processes are carried out.
- Repairs may only be carried out by persons authorised by the manufacturer.
- Only genuine spare parts and approved genuine accessories from Andreas Hettich GmbH & Co. KG may be used.
- Only the original power cable or a technical equivalent with the same technical data can be used.

- The safety and reliability of the incubator is only ensured if
  - the incubator is operated in accordance with the Operating Manual.
  - the electrical installation at the place of installation of the incubator complies with the pertinent requirements of the EN/IEC standards.
- The incubators have a ventilation air flap on the rear of the device. Hot air is expelled through the air flap when the device is in operation.
- The incubator has a temperature sensor. To check the readings of the temperature sensor, an external temperature measurement should be carried out.
- It is prohibited to:
  - store inflammable or volatile substances inside the unit – the unit is not provided with explosion protection,
  - touch electrically live parts of the unit,
  - operate the unit when your hands are wet,
  - place water containers on top of the unit,
  - climb onto the unit,
  - overload the shelves (the maximum load is specified in the technical data).

*The incubator is safe if used as intended and observing the descriptions and information given in this documentation.*

### CAUTION

The incubators have a ventilation air flap on the rear of the device. Do not insert your arm, fingers or any objects into the opening while the device is in operation. Such action can otherwise lead to personal injury or damage to equipment.

### 3 Device description

#### 3.1 Standard device


Fig. 1: Unit front side

- 1 Housing
- 2 Mains power switch
- 3 Device base
- 4 Door handle
- 5 Door lock
- 6 Touchscreen
- 7 USB port

Touchscreen including

- Temperature protection
- Tolerance range asymmetrically adjustable.
- Door alarm individually adjustable.

The incubator is equipped with a temperature monitor of protection class 3.1 in accordance with DIN12880:2007-05. The temperature monitor is used to protect the incubator (device protection), its surroundings and the sample material (sample protection) against impermissible temperature overshoots.


Fig. 2: Front of unit, door open

- 1 Temperature sensor
- 2 Inserting shelf
- 3 Interior made of high-quality stainless steel
- 4 Glass door (optional)
- 5 Glass door handle (optional)
- 6 Access port

For cable glands, use the access port on the side of the unit.

The incubator is fitted with a Ø 30-mm access port on the side of the unit.

Cables from external measuring systems can be fed into the interior through this opening.

*After routing the cables through the opening, the access port must be sealed on the unit's side to prevent temperature deviations in the interior. If the access port is not used, it must be sealed with a plug.*

*If external devices are used in incubators without cooling, overheating can quickly occur for which the incubator cannot compensate. The over-temperature can result in sample material being negatively affected. Long-term over-temperatures can damage the incubator.*

*When operating a thermal load inside the inner space, monitor the device*


Fig. 3: Rear of the unit

- 1 LAN connection
- 2 Air flap
- 3 Mains cable socket
- 4 Fuse

All devices have a key-operated lock. Two keys are included in the scope of supply, which can be found affixed to the rear of the device.

All devices have a sensor that detects whether the door is open or closed. When the door is open, the open door icon appears on the touchscreen. The number above the icon displays the counter which indicates how often the door has been opened.

Reset the counter by tapping the icon. The counter is also reset when the device is switched off.

**Free access to the incubator data**

*Due to the unsecured USB port on the front of the device, anyone can transfer data to a USB stick.*

The USB port on the front panel is only for transfer of data from the device's internal memory to an external data storage device or for service purposes

### 3.2 Accessories

#### Accessories

*Only genuine spare parts and approved genuine accessories from Andreas Hettich GmbH & Co. KG may be used.*

Brief description of accessories: Flexible equipment items with which the device can be retrofitted

Item number	description
60045	standard wire shelf for HettCube 60 <sup>1) 2)</sup>
60046	Perforated steel wire shelf for HettCube 60 <sup>1) 2)</sup>
60047	Reinforced shelf HettCube 60 <sup>1) 3)</sup>
60051	Internal glass door HettCube 60
60052	standard wire shelf for HettCube 120 <sup>1) 2)</sup>
60053	Perforated steel wire shelf for HettCube 120 <sup>1) 2)</sup>
60054	Reinforced shelf for HettCube 120 <sup>1) 3)</sup>
60055	Internal glass door for HettCube 120

<sup>1)</sup> All shelves are supplied with the appropriate support rails.

<sup>2)</sup> Maximum load: 25kg

<sup>3)</sup> Maximum load: 50kg

## 4 Transport and storage

### Storage conditions

#### DANGER

A danger is posed by penetrating liquids.

There is an electric shock hazard for personnel.

- If liquid has penetrated into the incubator, the incubator must not be put into operation.
- In such a case, contact a service technician of the manufacturer.

- The incubator can be stored in the original packaging.
- If the incubator is stored unpacked, then please remove all items from the chamber, disconnect the device from the mains, clean and dry the chamber and leave the door open to prevent unpleasant odours.
- The incubator may only be stored in dry rooms.
- The storage temperature must be between 0°C and +50°C.
- The humidity must be of non-condensing nature and between 10% and 70%.

### Receiving

#### NOTICE

When you receive the device, visually inspect it to assess its condition in the presence of the personnel delivering the goods, and ensure that the scope of supply is complete.

The courier company is responsible for any damage caused during transport.

### Transport

- Before transporting the incubator, make sure that all connection lines are removed.

#### DANGER

A danger is posed by penetrating liquids.

There is an electric shock hazard for personnel.

- If liquid has penetrated into the incubator, the incubator must not be put into operation.
- In such a case, contact a service technician of the manufacturer.

- The device may not be transported when it is loaded.
- During transport, note the weight of the incubator.
- When transporting with a transport aid (e.g. transport cart), this must be able to bear a weight at least 1.6 times the permissible transport weight of the incubator.
- Secure the incubator against tipping over and falling down during transport.

- The incubator may only be transported standing straight up.
- The incubator must not be lifted and transported by the door handle or the door.

## 5 Commissioning

- Ventilation openings may not be covered or blocked.
- Before commissioning, check to make sure the mains voltage is correct according to the nameplate and that the mains cable is connected properly.
- If necessary, additional devices and accessories can be mounted and/or connected.

Every time before using the incubator, make sure that

- the incubator is dry.

### 5.1 Unpacking the incubator

*When lifting the incubator, bear in mind the weight being lifted.*

*To avoid injury, only lift the incubator with the assistance of an appropriate number of helpers.*

*Please check the device and the outer packaging before unpacking. Report any transport damage immediately.*

1. → Remove the packaging straps.

2. →

*The packaging can be easily removed, even in rooms with a low ceiling. It is recommended that the room height be at least twice as high as the device.*

*If handled carefully, the packaging can be reused.*

Remove the packaging.

#### Unloading from the pallet

#### WARNING

**There is a risk of toppling when the incubator is taken off the transport pallet.**

There is a risk of the incubator toppling over. Persons could be injured if the incubator topples over.

- The instructions in this manual must be observed.

## 5.2 Setting up, connecting and switching on the incubator

### 5.2.1 Installing the incubator

#### NOTICE

If there is a cold-to-warm temperature difference, there is a danger that condensate will form on the electrical components of the incubator.

The forming condensate can cause a short circuit.


Allow the incubator to acclimate for at least three hours at the installation site before putting the device into operation.

- Please inspect the device for mechanical integrity, and ensure that all items in the delivery scope have been delivered by checking against your order.
- When choosing its place of installation, bear in mind the weight of the incubator when loaded, and its loading process; see → *Chapter 9 'Technical data' on page 73.*
- The installation location must not be in an area of direct sunlight or near to any sources of heat.
- The device must not be installed in the following places:
  - outdoors,
  - in damp places or places which can easily become flooded,
  - near flammable or volatile substances,
  - near acids or in corrosive environments.
- Ventilation openings must not be blocked. A clearance of 100 mm must be maintained from the ventilation slots and ventilation openings of the incubator.  
The room's ceiling height must be at least 300 mm higher than the height of the device.

1. → Position and level the incubator on a flat, non-combustible surface in a stable position.
2. → Adjust the height of the shelves and drawers as needed.

#### 5.2.1.1 Assembly stacking kit

Assembly fixing plate (lower incubator)


- 1 Screw
- 2 Fixing plate

1. → Unscrew the screws (1) in the lid of the lower incubator.
2. → Fasten the fixing plate (2) on the lid using the screws (1).

3. ▶ Position the lower incubator in the installation location.

Wall distance at least 100 mm but no more than 130 mm.  
Observe the instructions from chapter ▶ *Chapter 5.2 'Setting up, connecting and switching on the incubator'* on page 19.



**Assembly fixing kit (upper incubator)**

- 1 Fixing kit
- 2 Screw


1. ▶ Screw out the right and left upper screws on the rear wall.
2. ▶ Secure the fixing kit (1) on the right and left of the rear wall using the screws (2).

**Assembly of the upper incubator**

- 1 Upper incubator
- 2 Lower incubator
- 3 Device base
- 4 Fixing plate

1. ▶ Lift the incubator with an appropriate number of helpers and position on the lower incubator. The device feet (3) of the upper incubator (1) must be located in the fixing plates (4) of the lower incubator (2).

- 1 Socket screw
  - 2 Hole for fastening element
2. By loosening the socket screw (1) on the fixing kit, adjust the length at the wall distance.
  3. Retighten the socket screw (1).

**WARNING**
**Risk of injury due to falling loads**

Risk of injury due to insufficient pull-out resistance of the fastening elements.

- Select the fastening elements with a minimum pull-out resistance of 200 N.

4. Attach the fixing kit (2) to the wall using a suitable fastening element (e.g. rawl plug and screw).

*Alternatively another fastening option can be selected.*

*However, ensure that the fastening is implemented using 2 points at the top of the incubator and that the fixation elements each have a strength or pull-out resistance of 200 N.*

## 5.2.2 Incubator connection

### Power supply connection

**DANGER**

**A hazard is posed by a damaged mains supply line.**

There is an electric shock hazard for personnel.

- Check the incubator and the mains supply line for damage. If damage is discovered, the incubator must not be operated.
- The incubator may only be operated at a socket with a personal protection switch.

**DANGER**

A danger is posed by penetrating liquids.

There is an electric shock hazard for personnel.

- If liquid has penetrated into the incubator, the incubator must not be put into operation.
- The incubator may only be operated at a socket with a personal protection switch.
- In such a case, contact a service technician of the manufacturer.

**NOTICE**

There is a danger of damaging the incubator due to an incorrect input voltage.

- Before connecting the incubator, check the mains voltage.


Fig. 4: Power supply connection

- 1 Mains cable socket

1. →

*Ensure the mains power supplies the correct operating voltage as indicated on the rating plate!*

Connect the supplied mains power cable to the power supply socket on the rear of the incubator.

2. →

Plug the other end into the main power outlet.

*The mains cable must be freely accessible at all times so that the incubator can be disconnected from the mains.*

**Other connections**


Fig. 5: Other connections

- 1 LAN connection
- 2 Mains cable socket

- The installation should be protected by a 16A fuse and equipped with a residual current device.

## 5.3 Inserting and removing the slide-in modules

### CAUTION

Risk of injury due to sharp edges.

There may be sharp edges on the support rail.

- Wear cut-resistant gloves when inserting and removing the rails.

*When changing the rails of the shelves, the rails on the left and right sides must be at the same height.*

*Only change the configuration of the interior when the device is switched off and is not loaded.*

### 5.3.1 Standard slide-in modules

#### Inserting

#### *Note on the shelves*

*All shelves are supplied with the appropriate support rails.*

*Only use the slide-in units with the associated support rails.*


Fig. 6: Fitting a support rail in place


Fig. 7: Engaging the front end

The standard shelves are held by support rails.

1. ➤ Fit the support rail into the snap-in rail at the desired height at a 45 degree angle.
2. ➤ Engage the support rail in a downwards direction.
3. ➤ Fit the support rail on the opposite side.
4. ➤ Insert the shelves into the support rails.
5. ➤ Removing the shelves and support rails is done in the reverse order.

### CAUTION

The shelves are not secured with an end stop, and can fall down if pulled out completely. There is a risk of injury due to the shelf and load falling down.

## 5.4 Install glass door


Fig. 8: Installation of glass door (optional)

A glass door can be ordered and installed as an accessory.

Two people are required for the changeover.

1. ➤ Unscrew the existing screws on the desired hinge side.
2. ➤ Screw with two of these screws on the lower hinge (2) and with two of these screws on the upper hinge (2) of the door.  
Pay attention to the orientation of the glass door when screwing it on.
3. ➤ After attaching, check the function of the glass door (screws tightened enough, door is attracted by the magnet when it is closed.)

## 5.5 Changing the door swing

### 5.5.1 Door


Fig. 9: Changing the door swing

- 1 Upper hinge
- 2 Lower hinge
- 3 Space for lower hinge
- 4 Space for upper hinge

Two people are required for the changeover.

1. ➤ Remove the two screws (4) at the provided location of the upper hinge.
2. ➤ Remove the three screws from the lower hinge (2). The door must be securely held by a second person.
3. ➤ Remove the door.
4. ➤ Turn the door through 180 degrees (so that the bottom is now on top).

5. ➤ Remove the two screws from the upper hinge (1).
6. ➤ Remove the upper hinge (1).
7. ➤ Fit the lower hinge with 3 screws on the other side of the device (3). Loosely tighten the screws.
8. ➤ Fit the upper hinge (1) into the door.
9. ➤ Hang the door on the lower hinge. The door must be securely held by a second person.
10. ➤ Fit the upper hinge (4) with 2 screws.
11. ➤ Adjust the door so that the door seal makes contact. Only then tighten the screws of the lower hinge.

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Two people are required for the changeover.

1. ➤ Remove the two screws (2) from each of the hinge positions provided.
2. ➤ Remove the two screws from the lower hinge (1) and the two screws from the upper hinge (2) of the door. The door must be securely held by a second person.
3. ➤ Turn the door through 180 degrees (so that the bottom is now on top).
4. ➤ Screw in the two screws on the lower hinge (2) and the two screws on the upper hinge (2) of the door.
5. ➤ Screw two screws into each of the places of the previous hinges (1).
6. ➤ After attaching, check the function of the glass door (screws tightened enough, door is attracted by the magnet when it is closed.)

## 6 Operation

**DANGER**

A danger is posed by penetrating liquids.

There is an electric shock hazard for personnel.

- If liquid has penetrated the incubator during operation, do not touch the incubator and switch it off on the building side (e.g. via the fuse).
- In such a case, contact a service technician of the manufacturer.

**WARNING**

There is a risk of incorrect temperature control of the samples.

There is a risk of incorrect temperature control of microbiological samples, which can lead to misinterpretation of the test results or destruction of the samples.

- Arrange for maintenance and calibration to be carried out → *Chapter 7.5 'Maintenance' on page 69.*
- Note the information in this manual concerning the effect of heat generation from auxiliary equipment operated in the incubator.
- Bear in mind the temperature change profile.
- Do not place microbiological samples on the incubator floor. Use the shelves.
- Position the microbiological sample material within the usable space.
- Comply with the permissible ambient conditions.
- Prevent the incubator from being exposed to direct sunlight.
- Close the incubator door when not inserting or removing samples.
- After a power failure, check to ensure that the samples are still suitable for use.

**CAUTION**

A hazard is posed due to the improper pulling-out or pushing-in of the various slide-in modules or drawers. Avoid strong acceleration.

- Carefully pull out or push in slide-in modules, drawers or applications.

## 6.1 Operating elements


Fig. 11: Front-side controls

- 1 Door lock
- 2 Switch [Incubator On/Off]
- 3 USB port
- 4 Touchscreen

## 6.2 Loading

The maximum weight load per standard shelf is 25 kg.

Always load each level evenly.


Load the incubator so that the sample material is only inside the usable space and so that there is sufficient air circulation inside the incubator.

The following rules must be observed in order to ensure proper air circulation and stable conditions for storing the samples in the unit chamber:

- the max height of the samples should not exceed 1/3 of the space between the shelves;
- approx. 1/3 of the width and depth of the shelf should remain empty, while the distances between the samples as well as between the samples and the wall should be approximately equal.

The figure shows an example placement of samples in the chamber.

Following the above rules will provide optimum parameters for temperature fluctuation and variation.

## 6.3 Door locking mechanism

The incubator can be locked to prevent unauthorised persons from opening the door. Individual locking cylinders are used. If the keys are lost, the lock can be changed by the manufacturer's service department.

*Before locking, ensure you are completely certain that there are no persons inside the incubator.*

1.  Insert the key in the lock.
2.  Turn the key to the left to lock the incubator.
3.  Turn the key to the right to unlock the incubator.

The outer door is equipped with an open door sensor. If the door is not closed properly, an audible alarm is triggered. The door alarm delay can be set as follows: 30s, 1 min, 2 min, 5 min or 10 min.

## 6.4 Standard check before every use

Before each use of the incubator, ensure that

- the mains voltage corresponds to the specification on the rating plate;
- the incubator is level and sits securely on its feet;
- the door closes tightly;
- there are no outstanding alarm or error messages;
- the ventilation slots of the incubator are free from obstructions;
- the distances ( *Chapter 5.2.1 'Installing the incubator' on page 19*) are complied with;
- the incubator exhibits no obvious damage.
  - This applies to the incubator itself as well as to all accessories and optional equipment.

## 6.5 Switching on the incubator

-  Switch the incubator on by actuating the unit switch.
  - The control unit starts up.  
The start screen is displayed.

## 6.6 Unit operation

This symbol means that a given window can be moved in the direction shown in the picture.

### 6.6.1 External memory (USB flash drive)

The external memory (USB flash drive) can be used to copy data record, event log and service log from the unit's internal memory. Before first use, the USB flash drive must first be formatted with the FAT32 file system. To do this, the USB flash drive must be plugged into the USB port next to the display on the front of the device. After a few seconds, the message "USB flashdrive connected" is displayed on the screen. The USB flash drive must likewise be plugged into the USB port on the front side in order to transfer data. The following steps are then necessary:

#### *Recommendation*

*The manufacturer recommends the following USB sticks:*

- Kingston DataTraveler 100 64Gb
- SanDisk Ultra fit 64Gb
- SanDisk Cruzer Ultra 64Gb
- Goodram UUN2 8Gb

1. Go to the main menu .
2. Go to the data record .
3. Hit  and choose desired file type: \*.csv, or \*.plx.
4. Hit . The data have now been copied.

*After copying the data to the USB flash drive, before removing the flash drive from the USB port, it should first be unmounted by pressing the icon  in the top drop-down list. If the USB flash drive is not unmounted prior to connecting it to the computer, a message about USB flash drive damage may be displayed with a repair proposal, although the USB flash drive has actually not been damaged.*

*The only purpose of the USB port is to connect a flash memory – a USB flash drive or a card reader with a memory card. Other devices (e.g. an external hard disk) must not be connected without having first consulted the manufacturer. The USB port could otherwise be damaged.*

Data saved in the \*.csv file can be opened in a spreadsheet. Data saved as \*.plx can be opened in the LabDesk software program (accessory). This program allows, among other things, for data preview in the form of a table or a graph. You can also generate a report with statistics on a selected range of data.

### 6.6.2 First boot


Fig. 12: Downloading files

During first boot-up, the screen displays information about saving the "Download" folder to the USB flash drive. To do this, insert the USB flash drive and wait a moment until the hardware has been detected, then press  .

If you press  , you quit downloading the folder. The window will appear again during the next boot. You can tick "Don't show again" so that the window will no longer be displayed when the device is switched on. You can always download the "Download" folder in the Info submenu.

### 6.6.3 Main screen


Fig. 13: Main screen

After switching on the device, the main screen appears.

It contains information about the device status. After starting the program, additional information appears on the screen.


Fig. 14: Main screen – Running program

- 1 Main menu
- 2 Set temperature (temperature adjustment range 20°C - 65°C)
- 3 Segment duration time
- 4 Quick change of fan speed
- 5 Start/Stop the program
- 6 Current date / time
- 7 Active control panel
- 8 Program status
- 9 Information/control panel to be moved sideways
- 10 Alarm bar
- 11 Current temperature in the chamber

### 6.6.4 Information panel

There are two different windows in the information panel. Switching between them is done by swiping left or right.

Fig. 15: Information panel

The ●○ icon displays information on which window is active.

### 6.6.4.1 Alarms window

Fig. 16: Alarms window - icons

The icon  displayed on the second page of the information window allows you to go to the alarms window.

In the alarms window there's a list showing currently active alarms or the alarms that have occurred but have not yet been confirmed. When an alarm is active, the alarm bar is red and the alarm event is displayed in the list showing the status "active". When the alarm event stops, the state changes to "inactive".

- "Delete" button – confirms and removes the alarm from the list (only inactive alarms can be deleted).
- "Confirm" button – confirms an alarm.
- "Details" button – displays a preview of all instances of selected alarm.


Fig. 17: Alarms window

Fig. 18: Alarm details

### 6.6.4.2 Status window

Fig. 19: Status window - icons

The  icon displayed on the second page of the information window allows you to go to the status window.

The status of the device is also indicated by the description.


Fig. 20: Status description


program name	Name of the currently running program
current segment	Currently running segment / total number of segments in the program
priority	Program with time priority or parameter priority
time set	Set time of the currently running segment
time elapsed	Elapsed time since the segment started

time remaining	Remaining time until the end of the segment
current loop	Currently performing cycle / total number of cycles to perform
protection class	Information on the protection temperature of a running or completed program
over temperature under temperature	The protection parameters can be set in the program parameters
upper alarm temp lower alarm temp	Information on set alarms, separately for upper and lower temperature limit values. Setting the alarms

### 6.6.5 Meanings of symbols

	This icon enables you to go to the main screen.
	Automatic return to the home screen. Factory setting: disabled.
	This icon enables you to go to the main menu.
	Automatic screen lock. Factory setting: disabled.
	Unmounting the USB flash drive before removing it from the USB socket.
	Internal light is switched on. The icon appears in devices with an internal light (OWW). Light switched off – icon inactive. By clicking on the icon on the status window, you can switch the light On or Off manually (applies to devices with a viewing window on the door).

<p>4</p>  <p>4</p> 	<p>Closed door, open door. The number above the icon is a counter showing how often the door has been opened. Press the icon to cancel the counter. The counter is also cancelled by turning the unit off.</p>
	<p>For units with forced air convection: Fan icon. Rotating icon shows that the fan is running. Quick setting function – by clicking on the icon you can change the fan speed (without editing the program) within the range 0% to 100% in incubators sized up to 120 l, and from 10% to 100% in incubator sizes from 180 l - 1000 l.</p>
	<p>Ramp status: Chamber is currently heating up.</p>
	<p>Set temperature is reached.</p>
<p>2022.09.01 07:08</p> 	<p>The program will start on the given date / at the given time. Start delay activated.</p>
	<p>Icon is visible only when the chamber is heating up.</p>
	<p>Available when the program is running. Clicking the icon allows you to quickly change the set temperature (Quick Change function).</p>
	<p>When the program is running, click the icon to quickly change the time of program duration (Quick Change function). Indicates the time that has elapsed from the program start.</p>
	<p>Countdown of the time remaining to the end of the program.</p>
 	<p>The arrow icon allows navigation between segments, program parameters, and the summary overview.</p>
	<p>Starting the selected program. In the list of programs - The program is running.</p>
	<p>Stopping the program.</p>

	Adding a new program to the program list. The user can create up to 5 programs.
	Editing the selected program from the list. In the program list, a new program has been created but not yet approved.
	Removing selected program from the list.
	Goes to the alarms window.
	Going to the status window with information about the program parameters.
	Going to the menu to create, edit, delete and start programs.
	Cancelling adding or editing of the program. Cancel changes
	Editing individual program segments (the program can have 6 segments at maximum).
	Immediate start of the program selected from the program list.
	Delayed start of the program from the list of programs. The program starts according to the set date and time.
	Calling up the quick program function
	Turning off the alarm tone. Critical alarms (i.e. damage to the temperature sensor, temperature protection, etc.) continue emitting the alarm tone.
	When the program is running, click the icon to quickly change the fan speed (Quick Change function).
	In the program editing/creating mode: the degree of air-flap opening. In the operating mode (program started) clicking on the icon allows you to quickly change the air flap position (degree of opening) without having to edit the program.
<b>STM</b>	<p>The active Smart Temperature Monitor (STM) function informs the user of any problems incurred with reaching or maintaining the set temperature.</p> <ul style="list-style-type: none"> <li>■ White text - option enabled, the program is stopped</li> <li>■ Blue text - option enabled, the program is running</li> <li>■ Red text - warning about problems with reaching / maintaining the temperature</li> </ul>

### 6.6.6 Upper expandable and configurable menu

When the program is running, in the upper part of the main screen there is a bar menu with parameter icons (temperature, time, fan speed (only for incubators with forced convection) and air flap position (degree of opening)), which can be quickly changed (via the Quick Change function).

Depending on the model, you will see icons for all parameters that can be changed quickly by swiping your finger downwards. Among the options available in the bar you will see the following icons:

- Unmount USB flash drive.
- Mute option. Critical alarms (i.e. damage to the temperature sensor, temperature protection, etc.) continue emitting an alarm tone.
- Quick Change of:
  - program duration time,
  - set temperature,
  - fan speed and air flap position (degree of opening).

Fig. 21: Upper menu when the program is running

Fig. 22: Upper menu when the program is stopped

### 6.6.7 Alarm bar

The alarm bar provides quick visual information about the device status. The colour of the bar indicates the status of the device:

	turquoise bar – the device is working properly
	red bar and pulsing frame – active alarm

## 6.7 Quick Program

Quick program allows you to quickly start the program from the main screen position without having to enter the menu ☰.

The Quick program has several features that guarantee uninterrupted operation:

- You are unable to set the duration of the program – the time is always set to infinity;
- If the display fails, the program continues functioning;
- Once power supply has been restored after a power failure, the program resumes;
- The STOP button has been removed from the main window to prevent the program from being stopped inadvertently.





Fig. 23: Quick Program

Click the icon  in the main screen to go to the Quick Program. By clicking the appropriate icon you can set:

-  temperature
-  fan speed
-  air flap position (degree of opening)

Clicking the icon  starts the program in continuous mode (time set to infinity).

It has deliberately been made difficult to stop a Quick Program to prevent the program from being stopped inadvertently. To stop a program, you must go to the menu , click on the program window  and press and hold the STOP  button for five seconds. Over- and under-temperature protection is set automatically:

The over-temperature protection is the set temperature limit value of +5°C, while the under-temperature protection limit is the set temperature -5°C.

When the Quick Program is running, you can change the parameters (temperature, fan and air flap) by pressing the  or  icon. The next time you start the Quick Program, your previous settings will be saved.

After configuring the Quick Program, it appears in the list of programs. Quick Program is displayed at the top of the list by default.


Fig. 24: Quick program in program list

In Quick Program editing mode, you can change:

- the data recording interval setting,
- the protection class setting.

If the program is running, you can change the parameters (temperature, air flap position and fan speed (only for incubators with forced air circulation)) by pressing icon or . The next time you start the Quick Program, your previous settings will be saved.

## 6.8 Programs

In this window, you can run the selected program, add a new one, edit the program or delete it. The user can create 5 independent programs.

	Starting the selected program.
	Stopping the program.
	Adding a new program to the program list. The user can create up to 5 programs.
	Editing the selected program from the list.
	Removing the selected program from the list.

☰
programs: 4 / 5
🏠

name	segments	priority
<b>Program 1</b>	<b>1</b>	<b>parameters</b>
<b>Program 2</b>	<b>3</b>	<b>time</b>
<b>Program 3</b>	<b>1</b>	<b>parameters</b>
<b>Program 4</b>	<b>1</b>	<b>parameters</b>

Fig. 25: List of programs

Information on the number of created programs / the maximum number of programs that can be created is at the top of the screen (programs: 5/5).

### 6.8.1 Creating / editing a program

By pressing button or , a window with program parameters appears.

The program name is assigned automatically and cannot be changed.

In this window you can set:

- Number of segments – max. of 6 segments
- Interval – frequency at which data are saved in the data record (1 min, 2 min, 5 min, 10 min, 15 min, 30 min, 1 h)
- Protection class – for more information, see [Chapter 6.8.3 'Protection class'](#) on page 44.

- Protection temperature – temperature range for the protection class, for more information see .
- Priority – the priority of time or parameters; for more information see ➔ *Chapter 6.8.4 'Priority' on page 44.*
- Loop – the number of program repetitions; for more information see ➔ *Chapter 6.8.5 'Repetition / Loop' on page 45.*


Fig. 26: Program parameters

Cancelling the adding or editing of the program. Cancel changes.

Editing individual program segments (the program can have up to 6 segments at maximum).

Where multiple parameters are shown, the window can be scrolled up and down.

## 6.8.2 Segment editing

For each of the 5 programs, you can set time-temperature profiles with a maximum of 6 segments. This makes it possible to increase the incubation temperature of the samples in a step-wise fashion. This can protect the sample from thermal shock, for example. Example of program operation with programmed segments:

Program 1

Segment 1: temp. 30°C, time 2 hours (after reaching the temperature 30°C, it is maintained for 2 hours)

Segment 2: temp. 35°C, time 3 hours (after reaching the temperature 35°C, it is maintained for 3 hours)

Segment 3: temp. 40°C, time 3 hours (after reaching the temperature 40°C, it is maintained for 3 hours)

Segment 4: temp. 45°C, time 2 hours (after reaching the temperature 45°C, it is maintained for 2 hours)

Segment 5: temp. 50°C, time 2 hours (after reaching the temperature 50°C, it is maintained for 2 hours)

Segment 6: temp. 55°C, time 1 hour (after reaching the temperature 55°C, it is maintained for 1 hour)

Press button  and the first program segment is displayed.

In this window you can set:

- temperature – set temperature that the device should reach in this segment (must be at least 5°C less than the over-temperature protection value).
- time – the duration for which the set temperature ([d hh:mm]) in days, hours and minutes must be maintained. It is possible to select continuous operation ∞ in the final segment.
- ramp time – the time for reaching the set temperature ([d hh:mm]) in days, hours, minutes,
- fan - fan speed (power) in percent (only for incubators with forced air circulation),
- air flap – air flap position (degree of opening)
- fan ramp - fan speed (power) when the set temperature is reached (only for incubators with forced air circulation)
- air flap ramp - air flap position (degree of opening) upon reaching the set temperature.

The active value is highlighted in blue. The item highlighted in red means that the value is out of range and you should enter a different value, e.g. the temperature is above or below the operating range of the device or the protection temperature limit.

*By default, the fan speed (power) is set to 100%. Reducing the fan speed can lead to operating errors, e.g. increased temperature fluctuations and changes (only for incubators with forced air circulation).*

*Ramp time – setting a short time will not accelerate reaching the ramp, but the ramp will be reached in the shortest possible time depending on the set temperature, ambient conditions and the capabilities of the heating system in the device.*

The ramp parameters are factory-set in accordance with the manufacturer's instructions. If it is necessary to set individual parameters when reaching the segment temperature, activate the ramp edit field  and set your own values.

*Where multiple parameters are shown, the window can be scrolled up and down.*


Fig. 27: Editing program segments


Fig. 28: Editing program segments

Navigation between segments, program parameters and the summary overview is done by swiping sideways or by hitting the icon  .

*If you automatically return to the home screen while editing a program, the edited program will not be lost, rather it is saved as a draft (see below).*

After switching to program editing, a message appears about the option of further changing the program settings.


Fig. 29: Program segments - Info

### 6.8.2.1 Summary of segments

In the segments summary overview, all segments can be seen together with the set parameters:

- number of segments,
- temperature, duration time, target time for reaching the temperature of a given segment,
- fan speed,
- air flap position (degree of opening)


Fig. 30: Summary overview of segments

-  Confirms and saves the changes.
-  Cancels the entered changes in the segments and goes to the program parameters.

If multiple parameters are shown, the window can be scrolled up and down.

### 6.8.3 Protection class

The device is equipped with sample protection – temperature protection which is performed on the basis of the temperature value measured by an independent temperature sensor called the safety sensor. The main aim of sample protection is to protect against uncontrolled rise or fall of temperature. When activated, the relay disconnects the heating circuit.


Fig. 31: Confirmation of protection alarm class 2.0

There are three protection classes for the incubator:

**Class 1.0 in accordance with DIN 12880** – This protection class is always active.

**Class 2.0 in accordance with DIN 12880** – the user can choose between class 2.0 and class 3.1 in the settings.

**Class 3.1 in accordance with DIN 12880** – the user can choose between class 2.0 and class 3.1 in the settings.

For further explanations, see chapter ➔ *Chapter 6.21 'Temp. protection' on page 63.*

### 6.8.4 Priority

Can be set in terms of:

Parameter:

**In the program without ramp** – The device starts the countdown of the segment time when the set temperature is reached.

**In Programs with ramp** – The device starts to count down the ramp time and then switches to the segment countdown once the set temperature is reached, irrespective of whether the ramp time has elapsed.

*It may occur that the device fails to reach the set temperature within the set time because the set time duration was too short. In such situations, the time for reaching the set temperature will be extended and the segment's time countdown will start when the set temperature is reached.*

Time:

**In Programs without ramp** – The device starts counting down the segment time when the program starts, irrespective of whether the temperature has been reached.

**In the program with ramp** – First, the device counts down the ramp time and after that period expires it proceeds to the countdown of the segment time. It does so irrespective of whether the temperature has been reached.

*It may occur that the device fails to reach the set temperature within the set time because the set time duration was too short. In such cases, countdown of the segment time will start before reaching the set temperature. Thus, the actual time of device operating in the set temperature will be shortened.*

### 6.8.5 Repetition / Loop

This option is available if the number of segments is 2 or more. When the program finishes the final segment, the device starts the program again from the first segment. You can specify if the program should be carried out once (Loop: 1) or multiple times (Loop: 2 to 255). To set the program to run continuously, check the option "∞". If the time of the final segment is set to infinite, it will only be treated as infinite in the final cycle. In other cycles it will be treated as 0.

#### Example:

Loop:3

Segment 1: temp. 30°C, time 2 h

Segment 2: temp. 40°C, time 2 h

Segment 3: temp. 60°C, time "∞"

The device will run Segments 1 and 2 three times and then move on to Segment 3, which will continue indefinitely. In incubators that do not have a cooling system, the temperature in the chamber decreases naturally.

## 6.9 Starting the program

The created program can be started in two ways:

### 6.9.1 The first option

Go to the main menu and press the icon "Programs" .

Then select the program you want to activate and press the Start button .


Fig. 32: Main menu


Fig. 33: Program management menu

If the program is running, the symbol ▶ appears next to the program name on the list. The symbol ✎ means that the program has been edited, but the changes have not been confirmed.


Fig. 34: List of programs with the selected status

### 6.9.2 The second option

- In the main screen, press the icon in the upper right corner ⏸.
- Select the program you want to start. You have two additional options to start the program:

	Immediate start of the program.
	Scheduled program start according to the set date and time.


Fig. 35: Main screen

PROGRAM		
name	segments	priority
Program 1	2	parameters
Program 2	1	parameters
Program 3	1	parameters
Program 4	1	parameters

Fig. 36: Selection of the program from the list

A delayed start of a backdated program is possible (up to 7 days prior). This is possible for the programs with time priority. Program segments that would last from the back date to the current date are skipped.

If the program is running, the symbol ▶ appears next to the program name on the list.

## 6.10 Quick Change of parameters

Although the ramp time has been included in the program, the Quick Change of parameters will take place immediately while the temperature is being reached.

### 6.10.1 Quick Change of set temperature

In order to quickly change the value of the set temperature of a running program, press the icon in the main screen. The value of the temperature should be selected by scrolling through list upward or downward. Click to confirm the change.

The temperature cannot be higher than the over-temperature protection limit - 2°C or lower than the under-temperature protection limit + 2°C.


Fig. 37: Quick Change of set temperature


Fig. 38: Quick Change of set temperature

### 6.10.2 Quick Change of set time

In order to quickly change the duration of a running program, press the icon in the main screen. Select the number of days, hours and minutes by scrolling the list up or down. Click to confirm the change. To set continuous operation, press  $\infty$ .

To change the time display, press:

- – to display the elapsed time
- – to display the remaining time

To change only the time display, you do not have to confirm the change by hitting .


Fig. 39: Quick Change of set time


Fig. 40: Quick Change of set time

### 6.10.3 Quick change of fan speed and air flap position (degree of opening)

To quickly change the fan speed (only for incubators with forced recirculation) and/or the air flap position (degree of opening), press the  icon on the main screen. The value should be selected by scrolling the list up or down. Press  to confirm the change.


Fig. 41: Quick Change of the fan speed of the ventilator and/or the air flap position (degree of opening)


Fig. 42: Quick Change of the fan speed of the ventilator and/or the air flap position (degree of opening)

## 6.11 Statistics

Go to the main menu  and press the  icon. This window shows the statistics for the running or finished program. Statistics are calculated separately for each segment. Data logging for calculation starts after 30 seconds from reaching the set temperature in the segment.

Further data are registered every minute.

The following information is available:

- set temperature [°C] – the set temperature in the segment,
- minimum temperature [°C] – the lowest recorded temperature,
- maximum temperature [°C] – the highest recorded temperature,
- average temperature [°C] – the average temperature recorded,
- segment – status of the segment:
  - in progress – currently running segment (data are being constantly updated),
  - finished – the segment has been completed,
  - interrupted – the user interrupted the segment before the set time had elapsed,
- segment 1/2 – the number of the currently monitored segment / number of the currently performed or completed segment. Navigating between the segments is done by swiping up or down.

*You cannot overview any segment that has not started yet.*


Fig. 43: Statistics

## 6.12 Data record

Go to the main menu  and press the  icon. The data recording window contains the following information:

- time and date of sample registration [date],
- temperature value measured by the main sensor in the chamber [temp.].

It is possible to save 10,000 data records for a maximum period of 6 months. When the memory is full, the oldest data are overwritten. The data appear in the table in the order that they were added, not in chronological order by date.

The most recently added record is at the top. Data are only saved when the program is running. The frequency of date-saving depends on the set program parameters.

*When a data record is opened, all data are downloaded. If the user interrupts the data download, press  to continue downloading the rest of the data.*


Fig. 44: Data records

	Press to continue downloading data
	Recording data to the USB flash drive. Both .csv files (separated by semicolons when opening with a spreadsheet) and .plx files (opened using the LabDesk application (accessory)) are available.  Before removing the USB flash drive from the USB port, it must be unmounted.
	Deleting data


Fig. 45: Progress bar

**Deletion of data**  
  
 Are you sure you want to delete data?  
 This operation can not be undone!  
  
✓ ✕

Fig. 46: Deleting data

### 6.13 ⓘ Event log

Go to the main menu ☰ and press the ⓘ icon. The window displays information about recorded events, alarms and faults.

☰	date	name	code	🏠
	2023.03.24 11:38	<b>i</b> Program updated	2.01.0.1.024	•
	2023.03.24 11:38	<b>i</b> Temperature Correction Changed	2.01.0.1.028	
	2023.03.24 11:38	<b>i</b> Program Start	2.01.0.1.007	↑
	2023.03.24 11:37	<b>i</b> Program Stop	2.01.0.1.008	↓
	2023.03.24 11:37	<b>i</b> Program Edit	2.01.0.1.010	
	2023.03.24 11:37	<b>i</b> Program Edit	2.01.0.1.010	
	2023.03.24 11:37	<b>i</b> Program updated	2.01.0.1.024	
	2023.03.24 11:37	<b>i</b> Program Edit	2.01.0.1.010	

Fig. 47: Event log

	<p>Recording of data to the USB flash drive. Both .csv files (separated by semicolons when opening with a spreadsheet) and .plx files (opened using the LabDesk application (accessory)) are available.</p> <p>Before removing the USB flash drive from the USB port, it must be unmounted.</p>
	<p>Deleting data</p>

The events in the event log are sorted chronologically. However, while it may occur that the event "Program restarted" will not be displayed according to the chronology, the date and time of the event will be correct. This is not an error.

*Before removing the USB flash drive from the USB port, it must be unmounted.*

Information signs in the event log:

	<p>Information event</p>
	<p>Alarm event</p>
	<p>Error</p>
	<p>Warning</p>

Possible events:

<p>Program Start</p>	<p>Starting the program</p>
<p>Program Stop</p>	<p>Stop the program.</p>
<p>Program Edit</p>	<p>Changing the program parameters</p>
<p>Program End</p>	<p>Program is completed</p>
<p>Device On</p>	<p>The device is switched on (at the main switch)</p>
<p>Device Off</p>	<p>The device is switched off (at the main switch)</p>
<p>Door opened</p>	<p>The unit door is open</p>
<p>Open door alarm start</p>	<p>Open door alarm has been activated</p>
<p>Door closed</p>	<p>The unit door is closed</p>
<p>Open door alarm stop</p>	<p>Open door alarm has been deactivated</p>

Program Restarted	The program has resumed after a power failure
Under-Protection Start	Under-temperature protection has been activated
Under-Protection Stop	Under-temperature protection has been deactivated
Upper temp. alarm Start	Over-temperature protection has been activated
Upper temp. alarm End	Over-temperature protection has been deactivated
Date/time change	Date/time has/have been changed
Lower temp. alarm Start	Activation of the alarm if the temperature falls below the set temperature
Lower temp. alarm End	Deactivation of the alarm for when the temperature falls below the set temperature
Upper temp. alarm Start	Activation of the alarm if the temperature exceeds the set temperature
Upper temp. alarm End	Deactivation of the alarm for when the temperature exceeds the set temperature
Program saved	New program has been saved
Program deleted	Program has been deleted
Program updated	Program has been updated
Time Zone Changed	In the time settings, the time zone has been changed
Temperature Correction Changed	The main sensor temperature correction has been changed
Emergency stop of the program	The program has been automatically stopped – a situation occurred in which the program could not be continued. PLEASE CONTACT SERVICE
Power Fail Start	Power failure / Unit fuse has blown.
Power Fail Stop	Power supply has been restored, return to maintain program parameters

## 6.14 Info

Go to the main menu  and press the  icon. The window contains the following information:

- Name of the device
- Temperature range of the device
- Serial number of the device


- Software version
- Manufacturer's address
- Manufacturer's website

Fig. 48: Info window (example)

Press the icon to save the "Download" folder (with the Operating Manual) on the USB flash drive. After you have inserted the USB flash drive into the USB port, wait a few seconds until the message "Flashdrive connected" appears on the display. Press the icon to write the service data on the USB flash drive – contact Customer Service for more information.

Press to go to the main screen.

## 6.15 User settings panel

Go to the main menu and press the icon. In this window you can:

	Change the name of the unit – by default, the device serial number is entered.
	Change the language in the unit's menu.
	Set the time after which the screen will be dimmed.
	Turn the sound On or Off. Critical alarms will continue emitting an alarm tone.
	Set the time after which the user is returned to the home screen. Available settings: off, 1 min, 3 min, 5 min, 10 min.




Fig. 49: User settings panel

### 6.15.1 Unlocking the touchscreen

When the automatic touchscreen lock is enabled, slide the blue circle into the white circle to unlock the screen.

### 6.16 Time

Go to the main menu and press the icon. In this window you can:

- Change of the date / system time

*If the date/system time is/are changed to a later date/time and compared with the data and events stored in the memory, they will remain in the register. If the date/system time is/are changed to an earlier date/time than the date/time stored in the memory, they will be transferred to the archive.*

After changing the date/system time, the device will be restarted.

- Change time zone – changing the time zone will not affect the date/system time in the data and events previously saved.

To change the date/system time, press  in the window. A window appears in which you can make changes.


Fig. 50: Time zone change

	Confirms changes and restarts the device
	Cancel the entered changes



Fig. 51: Date change

## 6.17 Alarms

Go to the main menu  and press the  icon. You can set the alarm parameters.

- lower alarm – An alarm is activated if the temperature drops below the value indicated in this field.
- upper alarm – An alarm is activated if the temperature exceeds the value indicated in this field.

The lower and upper alarms are only active after the set temperature has been reached.

If the alarm is active, the alarm is triggered after a temperature drop following opening of a door.

- delay temp alarm: The alarm is activated with a delay (of 1 min, 2 min, 5 min, 10 min, or 15 min) after the permitted temperature limit is exceeded.
- delay door alarm: The door alarm is be activated when the door remains open for the time selected by the user (30 s, 1 min, 2 min, 5 min, or 10 min).
- Turning the STM function On or Off (to read the description of this function, go to Section → Chapter 6.17.1 'STM function' on page 58)


Fig. 52: Alarms

	Confirms changes
	Cancel the entered changes

In the field "lower alarm temp" you can set a value anywhere from -0.5°C to -5°C, and in the field "upper alarm temp" a value anywhere from +0.5°C to +5°C.

### 6.17.1 STM function

The Smart Temperature Monitor (STM) function informs the user if any problems incurred with reaching or maintaining the set temperature. The user can activate/deactivate the function. If the STM function is ON, the STM icon appears on the screen next to the temperature of the main sensor.


Fig. 53: Option enabled, the program is stopped

The function status is signalled by the following colour coding:

- No STM on the display – option disabled
- Black text – option enabled, the program is stopped
- Blue text – option enabled (temperature monitoring), the program is running
- Red text – option enabled, warning about problems with reaching or maintaining the temperature


Fig. 54: Option enabled (temperature monitoring), the program is running




Fig. 55: Option enabled, warning about problems with reaching or maintaining the temperature

Possible causes:

1.  damaged heater,
2.  the sample inserted into the chamber absorbs / dissipates too much energy.
3.  Door does not close tightly, air duct is open.

If the colour was red before opening the door, the colour changes to turquoise after opening the door.

In case of activation of the function (detection of problems with reaching or maintaining the temperature):

- Warning 4.00.0.1.009 appears in the event log
- The colour of the STM inscription changes to red and remains red as long as the error persists occurs, the segment is changed or the program is turned off
- When the status changes from red to turquoise, an entry about the end of the function 4.00.0.1.010 appears in the event log

### 6.17.2 Mute option.

The icon  in the main screen in the upper menu allows the alarm tones (open door alarm, exceeding temperature range) to be temporarily switched off, e.g. to prevent the door alarm from sounding while scheduled loading of samples into the chamber is underway. While it is possible to turn off the sound for 5, 10 or 15 minutes, the alarm tones of critical alarms (e.g. damage to the temperature sensor, over- and under-temperature protection) will be still emitted.



Fig. 56: Setting the time of the mute function

## 6.18 Network

Go to the main menu  and press the  icon. You can change the LAN settings in this window (Figure 51):

- IP – the device's IP address
- Mask – an Ethernet network mask to which the device is connected
- Gate – the IP address of the server or router that manages the Ethernet network
- DNS – the IP address of the domain name system
- MAC – the address of the network card, read-only
- DHCP – you can select whether the server that allocates IP addresses is to run on the local network. You can then skip setting IP addresses, masks and gates

The following icons indicate the connection status:

	Device connected to the network
	Device disconnected from the network


Fig. 57: LAN settings

	Confirms changes
	Cancels the entered changes

### 6.19 +/- Correction

Go to the main menu and press the +/- icon. In this window you can correct the temperature value shown on the display by adding the correction value. The set correction value applies to the whole temperature range of the device. For example, if the average temperature displayed by the device is 40.0°C and the average temperature measured by an independent, external sensor is 40.5°C, the correction should be set to +0.5°C. The mean temperature should be calculated over a chosen period of time, such as 30 minutes. The available correction range is between -5°C and +5°C.

*The device has been calibrated by the manufacturer in accordance with the applicable standards. The temperature shown on the display corresponds very closely to the temperature near to the sensor in the chamber. It is not necessary to utilise the user's calibration to ensure correct operation of the device.*

*Users conduct temperature correction under their own responsibility, and must be aware of consequences of changing the manufacturer's settings. If the user calibrates the unit, the existing calibration certificate loses its validity.*


Fig. 58: User's correction

	Confirms changes
	Cancels the entered changes

## 6.20 INTERFACE

### MODBUS TCP

The device allows status monitoring using the MODBUS TCP communication interface.

Connection parameters:

- IP address: the same as that of the device
- Port: 502

Register INPUT REGISTERS

function READ\_INPUT\_REGISTERS (0x30001)

Address	Type	Multiplier	Description
0	int	10	Temperature from the main sensor
3	Bool	-	Open door
4	bit	-	b0 – door alarm b1 – upper temperature alarm b2 – lower temperature alarm b3 – over-protection b5 – main sensor error b7 – protection sensor error b8 – temperature sensors error b10 – hardware error b11 – MRW error

## 6.21 Temp. protection


The device is factory-fitted with a sample protection system providing temperature protection. If any of the elements responsible for maintaining the set temperature are damaged, or the user inadvertently initiates a temperature setting, the set protection function will activate.

### Temperature protection class

By default, the incubator units are equipped with protection class 2.0 and 3.1 in accordance with DIN 12880. The user can then program the temperature value of the upper protection limit. If the set temperature limit is exceeded, the heating system shuts down.

With class 3.1: when the temperature returns to the allowed range, the device will resume operation.

The figure below shows how this works.

- 1 Malfunction
- 2 Class 1.0
- 3 Class 2.0
- 4 Class 3.1

Set temperature 44°C.

**Class 1.0:**

Device continues heating for as long as a malfunction exists (oscillates about the maximum temperature).

**Class 2.0:**

Device stops heating and cools down to the ambient temperature.

The user programs the protection temperature, and as soon as it is exceeded, the heaters are switched off and the alarm triggers. As long as the temperature remains above the protection temperature limit, the alarm remains on. To lower the temperature more quickly, the door can be opened (for example). The alarm stops as soon as the temperature is below the protection temperature limit. To restart the device, the user must turn the device Off and then On again.

**Class 3.1:**

Device stops heating as soon as the set over-temperature is reached. The temperature oscillates about the set temperature.

The user programs the protection temperature and once it's exceeded (caused by a failure), the heaters are cut off. Once the temperature returns to within the allowed range, the device resumes operation.

The set temperature in the segment cannot be higher than the over-temperature protection minus 2°C, e.g. if the over-temperature protection limit value is 50°C, then the maximum temperature that can be set in the segment is 48°C.

The temperature protection value for protection class 3.1 ranges:

- from the set temperature of +2°C to the maximum operating temperature of the device of +10°C.

## 7 Cleaning, disinfection and maintenance

### DANGER

#### Risk of injury due to electric shock

- Disconnect the device from the electric power supply before performing any cleaning work!
- The rechargeable battery back-up serving the controller should also be turned off.

### DANGER

#### A danger is posed by penetrating liquids.

There is an electric shock hazard for personnel.

### WARNING

#### There is a risk of contamination of the microbiological sample material.

If cleaning and/or disinfection are not sufficient after contamination, there is a risk that microbiological sample material will be contaminated with foreign germs.

- Be conscientious and thorough when cleaning and disinfecting.
- Follow the instructions of the cleaning agent and disinfectant manufacturers.
- Observe the range of effects of the used agent.

Observe the following rules when cleaning and disinfecting the unit:

- The incubator must **not** be cleaned in any equipment cleaning machines.
- The incubator must **not** be cleaned under running water, in water baths, or using compressed air.
- The units may only be cleaned manually and using liquid disinfection agents.
- To prevent any corrosion due to use of detergents and disinfectants, it is essential to follow the special application instructions provided by the manufacturers of the detergent and disinfectant.
- The water temperature must be between 20 and 25°C.
- Only detergents and disinfectants may be used that
  - are in the pH range 5-8,
  - do not contain caustic alkalis, peroxides, chlorine compounds, acids or alkalis.

*The outside surface of the incubator housing is powder-coated.*

*Disconnect the mains power plug before cleaning.*

*The rails and slide-in modules can be removed from the unit's interior to facilitate cleaning.*

*The operator must check with the manufacturer that the intended procedure will not damage the incubator before using any cleaning or decontamination procedures other than those recommended by the manufacturer.*

*Stainless steel products are manufactured with rust-proof steel. When used under standard laboratory conditions, they do not rust. However it is possible that stains (which may look like rust) to form on the steel surface (due e.g. to the kind of samples incubated in the chamber). In such cases, we recommend using a cleaning solution (to clean the stains) suitable for this particular application, such as Pelox.*

*When cleaning stainless steel products with a specialised cleaning solution, care should be taken to observe the suggestions and recommendations given in the instructions for use (or in the safety data sheet) of the cleaning solution.*

## 7.1 Cleaning

### Surface cleaning and care

#### **! WARNING**

The incubator may be contaminated.

There is a risk of contamination by bacteria, viruses or other pathogens. Based on its own hazard analysis, the operator must define whether and in what form personal protective equipment must be worn.

- ➔ ■ Clean the housing and the interior of the incubator at least once a week, if necessary using soap or a mild detergent and a damp cloth. Cleaning serves firstly to promote hygiene, and secondly to prevent corrosion due to adhering impurities.
- Ingredients of suitable detergents: Soap, anionic surfactants, non-ionic surfactants.
- After using detergents, remove any detergent residues by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- The interior must be checked annually for damage.
- Electrical parts should not come into contact with water or detergent.
- Clean **the touchscreen** using a soft cloth or cleansing foam for cleaning touchscreens.
- **The USB port** can be cleaned with a vacuum cleaner to prevent accumulation of dirt inside the port.

### Interior cleaning

#### ! WARNING

The incubator may be contaminated.

There is a risk of contamination by bacteria, viruses or other pathogens. Based on its own hazard analysis, the operator must define whether and in what form personal protective equipment must be worn.

1. → Before cleaning the interior of the device, empty the inner chamber.
2. → Open the door of the device and, if necessary, wait until the chamber has cooled down. Then remove the shelves and start cleaning the device.
3. → Only water or water with mild detergent should be used.
4. → Allow the device to dry completely after cleaning.
5. → During cleaning you should make sure not to damage the temperature sensors located inside the chamber.

#### 7.1.1 Cleaning the touchscreen

The touchscreen is exposed to dirt, so it must be cleaned regularly. To clean the touchscreen, use a clean, dry microfibre cloth. Such cloths are made of a very fine material that collects dirt well.

*Before using the cloth, make sure the surface is free of crumbs and other particles which during cleaning can act like sandpaper and scratch the surface of the screen.*

If stains cannot be removed by dry wiping, the cloth can be lightly dampened with water.

*Do not use paper towels to clean the screen as it may cause microdamage.*

Before cleaning, lock the screen by pressing  in the top drop-down list.


Fig. 59: Locking the screen

The touchscreen can now be cleaned.

To unlock the touchscreen, slide the blue circle into the white circle.

## 7.2 Disinfection

- ➔ ■ If infectious material gets into the interior, disinfect it immediately.
- Ingredients of suitable disinfectants: Ethanol, n-propanol, isopropyl alcohol, ethyl hexanol, corrosion inhibitors.
- After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
- The surfaces must be dried immediately after disinfecting.

### ⚠ WARNING

**Risk due to insufficient disinfection after contamination!**

There is a risk of infection for the user. Observe your laboratory rules (hygiene plan, etc.) for handling biological agents.

### 7.3 Removal of radioactive contaminants

- ■ The agent must be specifically labelled as being an agent for removing radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants: anionic tensides, non-ionic tensides.
- After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
- The surfaces must be dried directly after removing the radioactive contaminants.

### 7.4 Autoclaving

*Autoclaving accelerates the ageing of plastics. It can also cause changes to the colour of plastics.*

The shelves, support rails and the trays in the interior can be autoclaved at a temperature of 121°C / 250°F (for 20 min).

These parts must be removed from the unit interior before autoclaving. No statement can be made about the resulting degree of sterility.

### 7.5 Maintenance

**WARNING**

**There is a risk that the sample temperature control will be faulty.**

There is a risk that the temperature control of micro-biological samples is not correct, which could result in a misinterpretation of the test results or the destruction of the samples.

- Have maintenance and calibration done.

To ensure a long service life and error-free operation of the incubator, the incubator should be serviced at regular intervals by the manufacturer or by personnel authorised by the manufacturer.

The manufacturer offers maintenance agreements for this purpose. The operator is obliged to ensure that maintenance is carried out by an authorised service technician.

Interval	Maintenance work	Personnel
Before each use	Before each use, ensure that the incubator is in proper technical condition and that no restrictions are to be expected dependent on the desired application.	User
Annual maintenance	Annual maintenance - Inspection by one of the manufacturer's service technicians.	Service technician

Interval	Maintenance work	Personnel
Annual maintenance	Display on the screen - Inspection by one of the manufacturer's service technicians.	Service technician

## 8 Troubleshooting

*If the error cannot be remedied using the troubleshooting table, notify the service technician. Please specify the type of device and the serial number. Both numbers can be found on the name plate of the incubator.*

*If a malfunction occurs, there is an optical and acoustic alarm.*

*After a malfunction (power failure, door open too long, etc.), the operating company must assess whether the samples are still usable. He must therefore have the corresponding qualification.*

### 8.1 Troubleshooting

Before you contact Service Department:

1. ➤ Make sure that your operation of the incubator complies with the unit's Operating Manual.
2. ➤ Restart the unit to ascertain whether the unit is indeed not functioning properly. If it is still not functioning properly, disconnect the unit again from the mains and repeat the operation after one hour.

### 8.2 Possible defects

Malfunction	What needs to be checked?	What to you need to do?
The unit will not switch on.	Check if the unit is correctly plugged into the mains power source.	Correctly connect the unit to the mains socket-outlet.
	Check if the circuit-breaker has tripped.	Press the circuit breaker on the back of the device.
	Check the voltage at the mains socket-outlet.	Connect the device to a different socket-outlet, preferably on a different electrical circuit. Have the electrical installation checked by a licensed electrician.
	Check if the power cable is damaged.	Replace the power cable.
The unit is not heating up	Check if the unit door is properly closed.	Clean the seal.
	Check if the fan is switched on (only for incubators with forced air circulation).	Set the fan operation in the program.

Malfunction	What needs to be checked?	What to you need to do?
	Check if the ambient temperature is within the permissible value range specified in the technical data table.	Adjust the ambient temperature to the value specified in this Operating Manual.
	Check that thermal protection class 3.1 or class 2 has not been triggered.	Check the set values. Contact Customer Service.
The unit functions but is making too much noise.	Check that the device is not touching any other objects.	Remove any adjacent objects
The door has dropped or is skewed	Check if the door is properly levelled.	Level the unit using the adjustable feet. If this does not help, contact Customer Service.

## 9 Technical data

Technical data are specified within a tolerance of  $\pm 5\%$ , while the working capacity of the chamber is always smaller. All of the technical data below apply to standard incubator units (without optional accessories).

### 9.1 Incubator unit 60, 120

	HettCube 60 (natural convection)	HettCube 60 (forced convection)	HettCube 120 (natural convection)	HettCube 120 (forced convection)
Order no.	69601-10	69601-20	61201-10	61201-20
Mains voltage (+/- 10%)	220-240 V			
Mains frequency	50-60 Hz	50 Hz	50-60 Hz	50 Hz
power consumption	500 VA			
Incubator temperature range	5 K above ambient to +65°C	> 7.5 K above ambient to +65 °C	5 K above ambient to +65°C	> 6.3 K above ambient to +65 °C
Adjustable temperature	20°C to 65 °C			
External dimensions (without access port) W x D x H in mm	590 x 620 x 705	590 x 620 x 705	660 x 710 x 850	660 x 710 x 850
Internal dimensions W x D x H in mm	400 x 360 x 390	400 x 360 x 390	460 x 450 x 540	460 x 450 x 540
Internal volume in l	56	56	112	112
Footprint in m <sup>2</sup>	0.36	0.36	0.46	0.46
Weight in kg	62	62	82	82
Number of shelves (included in scope of supply)	2	2	2	2
Max. number of shelves	5	5	7	7
Temporal temperature deviation (+37°C)	+/- 0.7 K	+/- 0.1 K	+/- 0.4 K	+/- 0.2 K
Spatial temperature deviation (+37°C)	+/- 0.8 K	+/- 0.3 K	+/- 0.8 K	+/- 0.3 K
Maximum load capacity (kg)	40 kg	40 kg	60 kg	60 kg

	HettCube 60 (natural convection)	HettCube 60 (forced convection)	HettCube 120 (natural convection)	HettCube 120 (forced convection)
Ambient conditions (EN/IEC 61010-1)				
Installation site	indoors only			
Altitude	up to 2000 m above sea level			
Ambient temperature	10 °C to 28 °C			
Humidity	Maximum relative air humidity 75 %, non-condensing.			
Overvoltage category (IEC 60364-4-443)	II			
Pollution level	2			
Device protection class	I			
Not suitable for use in potentially explosive atmospheres.				
EMC				
Emitted interference, interference immunity	EN / IEC 61326-2-6, Class B;			
IP rating (EN 60529)	IP 20			
	HettCube 60 (natural convection)	HettCube 60 (forced convection)	HettCube 120 (natural convection)	HettCube 120 (forced convection)
Order no. (100-127 V)	69601-11	69601-21	61201-11	61201-21
Mains voltage (+- 10%)	100-127 V			
Mains frequency	50-60 Hz	60 Hz	50-60 Hz	60 Hz
power consumption	620 VA			
Incubator temperature range	5 K above ambient to +65°C	> 7.5 K above ambient to +65 °C	5 K above ambient to +65°C	> 6.3 K above ambient to +65 °C
Adjustable temperature	20°C to 65 °C			
External dimensions (without access port) W x D x H in mm	590 x 620 x 705	590 x 620 x 705	660 x 710 x 850	660 x 710 x 850

	HettCube 60 (natural convection)	HettCube 60 (forced convection)	HettCube 120 (natural convection)	HettCube 120 (forced convection)
Internal dimensions W x D x H in mm	400 x 360 x 390	400 x 360 x 390	460 x 450 x 540	460 x 450 x 540
Internal volume in l	56	56	112	112
Footprint in m <sup>2</sup>	0.36	0.36	0.46	0.46
Weight in kg	62	62	82	82
Number of shelves (included in scope of supply)	2	2	2	2
Max. number of shelves	5	5	7	7
Temporal temper- ature deviation (+37°C)	+/- 0.7 K	+/- 0.1 K	+/- 0.4 K	+/- 0.2 K
Spatial temperature deviation (+37°C)	+/- 0.8 K	+/- 0.3 K	+/- 0.8 K	+/- 0.3 K
Maximum load capacity (kg)	40 kg	40 kg	60 kg	60 kg
Ambient conditions (EN/IEC 61010-1)				
Installation site	indoors only			
Altitude	up to 2000 m above sea level			
Ambient tempera- ture	10 °C to 28 °C			
Humidity	Maximum relative air humidity 75 %, non-condensing.			
Overvoltage cat- egory (IEC 60364-4-443)	II			
Pollution level	2			
Device protection class	I			
Not suitable for use in potentially explosive atmospheres.				
EMC				
Emitted interfer- ence, interference immunity	EN / IEC 61326-2-6, Class B;			
IP rating (EN 60529)	IP 20			

## 9.2 Type plate


Fig. 60: Type plate

- 1 Manufacturer logo
- 2 Product designation
- 3 Type number
- 4 Serial number
- 5 possibly EAC marking, CE marking
- 6 Year of manufacture
- 7 Symbol: Observe the operating instructions.
- 8 Manufacturer address
- 9 Temperature range / nominal temperature
- 10 Power values
- 11 Voltage range
- 12 Mains frequency

## 9.3 Dimensions

### Dimensions


Fig. 61: Dimensions

	HettCube 60	HettCube 120
A	590 mm	660 mm
B	710 mm	850 mm
C	620 mm	710 mm

Tab. 1: Dimension (without 50-mm mains plug)

	HettCube 60	HettCube 120
D	400 mm	460 mm
E	390 mm	540 mm
F	360 mm	450 mm

Tab. 2: Interior dimensions

## Access port dimensions


Fig. 62: Access port dimensions

Ø30 mm	HettCube 60	HettCube 120
G	approx. 345 mm	approx. 310 mm
H	approx. 275 mm	approx. 270 mm
I	approx. 20 mm	approx. 20 mm

## 10 Disposal

Before its disposal, the incubator must be decontaminated and cleaned to protect people, the environment and materials.

All applicable pertinent legal regulations must be observed when disposing of the incubator.

According to Directive 2012/19/EU , laboratory incubator devices may no longer be disposed of with household waste.

Incubators without a cooling function belong to group 4 (large devices) and are assigned to the business-to-business sector.

Incubators with a cooling function belong to group 1 (heat exchangers) and are assigned to the business-to-business sector.

The crossed-through waste-bin symbol indicates that the incubator must not be disposed of with household waste.

Regulations governing disposal of such devices may differ in individual countries. If necessary, please contact your supplier to learn more.

## 11 Free and Open Source Software

### Open Source Licensed Used

Our SMART unit controllers use the following components:

Component	License
jfoenix 9.0.10	MIT
io.dataFx dataFx 8.0.7	BSD 3-Clause
io.dataFx flow 8.0.1	BSD 3-Clause
fontawesomefx 8.2	Apache 2.0
joda-time 2.10.10	Apache 2.0
org.jfree jfreechart 1.5.0	GNU LGPL
org.jfree jfreechart-fx 1.0.1	GNU LGPL
com.github.kenglxn.qrgen.javase 2.6.0	Apache 2.0
org.apache.maven.plugins 2.5.2	Apache 2.0
org.springframework 2.3.10.RELEASE	Apache 2.0
com.fasterxml.jackson 2.12.3	Apache 2.0
com.google.guava guava 22.0	Apache 2.0
org.mockito mockito-core 3.9.0	MIT
joda-time joda-time 2.10.10	Apache 2.0
org.scream3r jssc 2.8.0	GNU LGPL 3
com.intelligt.modbus jlibmodbus 1.2.9.1	Apache 2.0
commons-io 2.8.0	Apache 2.0
org.slf4j slf4j-ext 1.7.25	MIT
org.powermock	Apache 2.0
com.sun.mail javax.mail 1.6.1	CDDL/GPLv2+CE
net.samuelcampos usbdrivedetector 2.0.4	Apache 2.0
net.lingala.zip4j zip4j 1.3.2	Apache 2.0
de.undercouch bson4jackson 2.7.0	Apache 2.0
org.postgresql postgresql 42.2.20	BSD-2-Clause
flask 1.2.2	BSD License (BSD-3-Clause)

Component	License
pyserial 3.5	BSD License
unicorn 20.1.0	MIT License (MIT)
itsdangerous 2.0.1	BSD License (BSD-3-Clause)
ntplib 0.4.0	MIT License (MIT)
requests 2.26.0	Apache Software License (Apache 2.0)
librouteros 3.2.0	GNU General Public License v2 (GPLv2) (GNU GPLv2)
jinja2 3.0.3	BSD License (BSD-3-Clause)
warkzeug 2.0.3	BSD License (BSD-3-Clause)
markupsafe 2.1.1	BSD License (BSD-3-Clause)
certify	Mozilla Public License 2.0 (MPL 2.0) (MPL-2.0)
charset-normalizer 2.0.12	MIT License (MIT)
click 8.1.3	BSD License (BSD-3-Clause)
colorama 0.4.5	BSD License (BSD)
idna 3.3	BSD License (BSD-3-Clause)
setuptools 65.3.0	MIT License
urllib3 1.26.12	MIT License (MIT)
psycopy2	GNU Library or Lesser General Public License (LGPL) (LGPL with exceptions)
Bellsoft Liberica JDK8	GPL2 + <a href="https://bell-sw.com/assets/files/LibericaJDK-8-9-10-licenses.pdf">https://bell-sw.com/assets/files/LibericaJDK-8-9-10-licenses.pdf</a>
Bellsoft Liberica	JDK15 GPL2 + <a href="https://bell-sw.com/assets/files/LibericaJDK-11-0-11-3-d-party%20licenses.pdf">https://bell-sw.com/assets/files/LibericaJDK-11-0-11-3-d-party%20licenses.pdf</a>
PostgreSQL	PostgreSQL License, a liberal Open Source license, similar to the BSD or MIT licenses

#### Licences URL:

Component	License
Apache 2.0	<a href="http://www.apache.org/licenses/LICENSE-2.0">http://www.apache.org/licenses/LICENSE-2.0</a>
BSD-2-Clause	<a href="https://jdbc.postgresql.org/about/license.html">https://jdbc.postgresql.org/about/license.html</a>
BSD 3-Clause	<a href="https://opensource.org/licenses/bsd-license.php">https://opensource.org/licenses/bsd-license.php</a>
MIT	<a href="https://github.com/mockito/mockito/blob/release/3.x/LICENSE">https://github.com/mockito/mockito/blob/release/3.x/LICENSE</a>
GNU LGPL	<a href="http://www.gnu.org/licenses/lgpl.txt">http://www.gnu.org/licenses/lgpl.txt</a>
GNU LGPL 3	<a href="http://www.gnu.org/licenses/lgpl-3.0.txt">http://www.gnu.org/licenses/lgpl-3.0.txt</a>

Component	License
CDDL/GPLv2+CE	<a href="https://javaee.github.io/javamail/LICENSE">https://javaee.github.io/javamail/LICENSE</a>
Mozilla Public License 2.0 (MPL 2.0) (MPL-2.0)	<a href="https://www.mozilla.org/en-US/MPL/2.0/">https://www.mozilla.org/en-US/MPL/2.0/</a>
OSI Approved (Eclipse Public License v2.0 / Eclipse Distribution License v1.0)	<a href="https://www.eclipse.org/legal/epl-2.0/">https://www.eclipse.org/legal/epl-2.0/</a>
GPL2	<a href="https://openjdk.org/legal/gplv2+ce.html">https://openjdk.org/legal/gplv2+ce.html</a>
PostgreSQL License	<a href="https://opensource.org/licenses/postgresql">https://opensource.org/licenses/postgresql</a>

## 12 Index

<b>A</b>		<b>S</b>	
Access port dimensions. . . . .	77	Scope of delivery. . . . .	10
Accessories. . . . .	15	Signal words. . . . .	7
Autoclaving. . . . .	69	Standard check. . . . .	28
<b>C</b>		Standard shelves	
Cleaning. . . . .	66, 67	inserting. . . . .	23
Commissioning. . . . .	18	Standard unit. . . . .	14
Connection		Start screen. . . . .	28
Power supply. . . . .	21	Storage conditions. . . . .	16
Connections. . . . .	27	Switching the unit on. . . . .	28
Controls. . . . .	27	Symbols. . . . .	8
<b>D</b>		<b>T</b>	
Dimensions. . . . .	76	Troubleshooting. . . . .	71
Disinfection. . . . .	68	Type plate. . . . .	76
Disposal. . . . .	78	<b>U</b>	
Door lock. . . . .	28	Unit description. . . . .	14
<b>G</b>		Unit switch. . . . .	28
Graphic symbols on the incubator. . . . .	8	Unpacking the incubator. . . . .	18
<b>I</b>		<b>W</b>	
Incubator assemblies. . . . .	27	Warning symbols. . . . .	7
Installing the incubator. . . . .	19		
Assembly stacking kit. . . . .	19		
Intended use. . . . .	9		
<b>L</b>			
Labels			
on the packaging. . . . .	8		
Loading. . . . .	27		
<b>M</b>			
Maintenance. . . . .	65, 69		
<b>N</b>			
Non-intended use. . . . .	9		
<b>O</b>			
Operation. . . . .	26		
Original replacement parts. . . . .	10		
Other connections. . . . .	22		
<b>P</b>			
Personal protective equipment. . . . .	9		
Personnel qualification. . . . .	9		
Personnel qualifications. . . . .	9		
Protective equipment. . . . .	9		
<b>R</b>			
Radioactive contaminants			
Removal. . . . .	69		
Replacement parts. . . . .	10		
Return. . . . .	11		

## 13 Appendix

## A List of standards

### Standards and regulations applicable to this device

The device is a high-technology product. It is subject to extensive testing and certification procedures in accordance with the latest version of the following standards and regulations:

### Electrical and mechanical safety for construction and final testing:

Standard series: IEC 61010 (corresponds to the standard series DIN EN 61010)

- IEC 61010-1 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1: General requirements" (pollution degree 2, overvoltage category II)
- IEC 61010-2-010 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials"

### Electromagnetic compatibility:

EN 61326-1 "Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements"

### Risk management:

DIN EN ISO 12100 "Safety of machinery — General principles for design — Risk assessment and risk reduction"

### Restriction of Hazardous Substances (RoHS II):

EN 50581 "Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances"

### Certified quality management system in accordance with

- ISO 9001 "Quality management systems - Requirements"

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