

ENVIRONMENT/WATER ANALYSIS

The use of HettCube incubators to detect the presence of Pseudomonas aeruginosa in drinking water

Pseudomonas aeruginosa is an opportunistic pathogen that is ubiquitous in environments containing sufficient moisture. The micro-organism thrives in a wide range of conditions, from 9 °C to 42 °C, and has a high natural resistance to antibiotics. It does not generally pose a danger to healthy people who ingest the organism through drinking water, but it can cause inflammation if it comes into contact with skin lesions. The danger lies in the risk of nosocomial infections. The risk is highest for immunosuppressed patients and those in intensive care. It is therefore necessary for clinics and other facilities caring for people to make certain that their water is free of P. aeruginosa. It must therefore be tested every 6 months in accordance with a guideline of the Robert Koch Institute.

Methods of detection

The presence of P. aeruginosa is determined by the method given in EN ISO 16266 through membrane filtration and incubation of the filter on Cetrimide agar. The colour change and fluorescence under UV light indicate its presence. Presumptive colonies are cultured further on agar and tested to establish if they form ammonia from aceta-mide, are oxidase-positive or fluoresce in King B agar.

Importance of the test for the presence of Pseudomonas aeruginosa in practice

Drinking water is tested at regular intervals to protect consumers. One such test is for the presence of P. aeruginosa. The nature and frequency of the tests is laid down in national regulations such as the "Trinkwasserverordnung" [Drinking Water Act] in Germany and the Safe Drinking Water Act in the USA.

Incubation conditions in accordance with EN ISO 16266:2008

	Temperature	Duration
Agar and broth	36 ± 1 °C	40 - 48 h

Advantages of HettCube incubators

- Maximal validated usable space on a small footprint
- 4.3 inch touch display for intuitve operation
- Very homogeneous and stable temperature, as well as precise temperature control
- True "one-hand-operation" and flexible positioning of the shelves
- Minimal energy consumption of < 0.06 kW/h at 37 °C
- Low noise level of \leq 44 dB(A)
- Optimierte Beladungskapazität durch einzigartiges Zubehör und Optionen



Fig. 1 *: Detection of P. aeruginosa on Cetrimide agar



Fig. 2 *: Demonstration of its presence using the oxidase test (left) and acetamide broth (right)

Image by courtesy of the CVUA Karlsruhe

Hettich solution

HettCube 200 62000 HettCube 400 64000 HettCube 600 66000 Model without IVD Cat. No.	Model	Cat. No.
HettCube 600 66000 Model Image: Comparison of the second	HettCube 200	62000
Model	HettCube 400	64000
	HettCube 600	66000
	Model without IVD	Cat. No.

without IVD	Gal. NO.
HettCube 200	62001
HettCube 400	64001
HettCube 600	66001