UNIVERSAL 320 / 320 R

Inhalt des Dokuments / content of the document

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Operating instructions (EN)

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Operating instructions

UNIVERSAL 320 / 320 R

Translation of the original operating instructions

©2022 - All rights reserved

Andreas Hettich GmbH

Föhrenstraße 12

D-78532 Tuttlingen, Germany

Telephone: +49 (0)7461/705-0

Fax: +49 (0)7461/705-1125

Email: info@hettichlab.com, service@hettichlab.com

Internet: www.hettichlab.com

Table of contents

1	About this document.	5
1.1	Use of this document.	5
1.2	Gender reference.	5
1.3	Symbols and labels in this document.	5
2	Safety.	5
2.1	Intended use.	5
2.2	Personnel requirements.	6
2.3	Operator's responsibility.	7
2.4	Safety instructions.	7
3	Device overview.	9
3.1	Technical data.	9
3.2	Certifications and logos.	13
3.3	European registration.	13
3.4	Important labels on the packaging.	14
3.5	Important labels on the device.	14
3.6	Operating and indicator elements.	15
3.6.1	Control.	15
3.6.2	Indicator elements.	15
3.6.3	Controls.	16
3.7	Original spare parts.	16
3.8	Scope of delivery.	16
3.9	Returns.	17
4	Transport and storage.	17
4.1	Transport and storage conditions.	17
4.2	Fastening the transport lock.	18
5	Commissioning.	19
5.1	Unpacking the centrifuge.	19
5.2	Removing the transport lock.	19
5.3	Setting up and connecting the centrifuge.	20
5.4	Switching the centrifuge on and off.. . . .	22
6	Operation	22
6.1	Opening and closing the lid.	22
6.2	Removing and installing the rotor.	23
6.3	Inserting and removing buckets.	25
6.4	Inserting and removing adapters.	26
6.5	Loading.	26
6.6	Opening and closing the biosafety system.	28
6.6.1	Explanation.	28
6.6.2	Lid with screw cap and hole	29
6.6.3	Lid with bracket and spring-type lock.	29
6.6.4	Lid with screw cap.	30

6.7	Centrifugation.	30
6.7.1	Centrifugation in continuous operation.	30
6.7.2	Centrifugation with time preselection.	31
6.7.3	Short-term centrifugation.	31
6.8	Quick stop function.	31
7	Software operation.	32
7.1	Centrifugation parameters.	32
7.1.1	Relative centrifugal force, RCF.	32
7.1.2	Centrifugation of substances or mixtures of substances with a density higher than 1.2 kg/dm ³	32
7.2	Programming.	33
7.2.1	Write protection for programs.	33
7.2.2	Opening or loading programs.	33
7.2.3	Entering or changing programs.	33
7.3	Rotor detection.	34
7.4	Cooling (for centrifuges with cooling).	34
7.4.1	Instructions, cooling.	34
7.4.2	Standby cooling.	34
7.4.3	Precooling the rotor.	34
7.5	Machine Menu.	35
7.5.1	Querying system information.	35
7.5.2	Querying operating hours.	35
7.5.3	Audible signal.	36
7.5.3.1	General.	36
7.5.3.2	Setting an audible signal.	36
8	Cleaning and care.	36
8.1	Overview table.	36
8.2	Instructions for cleaning and disinfection.	37
8.3	Cleaning.	37
8.4	Disinfection.	38
8.5	Maintenance.	39
9	Troubleshooting.	40
9.1	Fault description.	40
9.2	Performing a MAINS RESET.	42
9.3	Emergency release.	42
10	Disposal.	42
10.1	General instructions.	42
11	Index.	44

1 About this document

1.1 Use of this document

- Read this document carefully and in full before commissioning the device for the first time.
Observe other enclosed instruction sheets where necessary.
- This document is part of the device and must be kept within easy reach.
- This document must be included if the device is passed on to a third party.
- The current version of the document in the available languages can be found on the manufacturer's website: ➡ <https://www.hettichlab.com/de/download-center/>

1.2 Gender reference

The employed masculine or feminine language form is to facilitate reading. In the spirit of equal treatment, corresponding terms apply in principle to all genders and do not imply any valuation.

1.3 Symbols and labels in this document

General symbols

The following markers are used in this document to highlight instructions, results, listings, references and other elements:

Marker	Explanation
1.  2.  3.  ... 	Step-by-step instructions
	Results of action steps
	References to sections of the document and other applicable documents
■ ... ■ ...	Listings without a fixed order
[Buttons]	Controls (for example: buttons, switches)
'Indicator'	Indicator elements (for example: signal lights, screen elements)

2 Safety

2.1 Intended use

Intended use

The centrifuge **UNIVERSAL 320 / 320 R** is an in vitro diagnostic medical device according to the In Vitro Diagnostic Medical Devices Regulation (EU) 2017/746.

The device is used to separate samples of human origin into their constituent parts for subsequent further processing. The user can set each of the variable physical parameters within the limits set by the device.

The centrifuge may only be used by qualified personnel in closed laboratories. The centrifuge is only intended for the use referred to above. Intended use also includes observing all instructions in the user manual and compliance with inspection and maintenance. Any other use or use beyond this scope is considered improper. Andreas Hettich GmbH shall not be liable for any damage arising from such non-compliant use.

Non-intended use

- The centrifuge is not suitable for use in explosive or radioactive, or biologically or chemically-contaminated atmospheres.
- The user must take appropriate actions when centrifuging hazardous substances or mixtures of substances that are toxic, radioactive or contaminated with pathogenic microorganisms.
The manufacturer generally recommends using only centrifuge tubes with special screw caps designed for use with hazardous substances. Use sealable centrifuge tubes with a biosafety system for materials of risk groups 3 and 4.
- The manufacturer does not recommend centrifugation of flammable or explosive materials.
- The manufacturer does not recommend centrifugation of materials that react chemically with one another with high activation energy.

Foreseeable misuse

The manufacturer recommends using only accessories approved that it has approved for the intended purpose.

Only operate the centrifuge under supervision.

2.2 Personnel requirements

Required qualifications

The user has read the Operating Manual in full and familiarised themselves with the device.

NOTICE

Damage to the device by unauthorised personnel

- Tampering with and modifications to devices by unauthorised persons are at the operating organisation's own risk and will result in the loss of all warranty and liability claims.

Trained user

The user has been educated and trained in laboratory work and is able to carry out the work assigned to them, and to recognise and prevent potential hazards independently.

Personal protective equipment

Lack of personal protective equipment or unsuitable personal protective equipment increases the risk of impaired health and injury.

- Only use personal protective equipment that is in proper condition.
- Only use personal protective equipment that is adapted to the person (correct size, for example).
- Observe instructions on other protective equipment for specific activities.

2.3 Operator's responsibility

Follow the instructions in this document for proper and safe use of the device.

Keep the user manual for future reference.

Provide information

- Following the instructions in this document will help:
 - To avoid dangerous situations.
 - To minimise repair costs and downtime.
 - To increase the reliability and service life of the device.
- The operator is responsible for compliance with company regulations, standards and national laws.
- Note and keep the revision of the document separate from the document. If lost, the document can be replaced in the correct revision.
- Keep the user manual available at the place where the device is used.
- Pass the user manual on to the buyer when the device is sold.

Personnel training

Lack of knowledge when working with the device may result in serious injury or death.

- Instruct personnel on their tasks and the associated risks in accordance with the instruction.

2.4 Safety instructions

Reporting of events and incidents

In the event of incidents or notifiable events involving the device or its accessories, these must be reported to the manufacturer and, where applicable, to the competent authority where the user and/or the patient is registered.

Manufacturer:

Andreas Hettich GmbH

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78532 Tuttlingen, Germany

Phone: +49 7461 705 0

E-mail: info@hettichlab.com

DANGER

Risk of contamination

Inadequate cleaning or failure to observe the cleaning instructions can lead to contamination risks.

- Observe national and local regulations on safety and accident prevention.
- Observe cleaning instructions.
- Wear personal protective equipment when cleaning the device.
- Comply with laboratory regulations (e.g. TRBAs, the German Protection against Infection Act, hygiene plan) for handling biological agents.

DANGER

Risk of fire and explosion

Risk of accidents, injuries or damage to property due to fire or explosion.

- Observe regulations and directives for handling chemicals and hazardous substances.
- Do not use corrosive chemicals.
- Do not use dangerous chemicals.
- Do not use corrosive extraction agents.
- Do not use strong acids.

WARNING

Risk of injury

Inadequate or late maintenance can result in injuries.

- Observe maintenance intervals.
- Check the device for visible damage or defects.
If any visible damage or defects are present, immediately remove the device from service and inform a service technician.

WARNING

Risk of electric shock

Liquids that penetrate the device may cause electric shocks.

- The device must be protected from external contact with liquids.
- Do not pour any liquids inside the device.
- The original transport packaging must be used when transporting the device.

WARNING

Contamination with hazardous substances and substance mixtures

Observe the following actions for substances and substance mixtures that are toxic, radioactive and/or contaminated with pathogenic microorganisms:

- As a rule, use only centrifuge tubes with special screw caps for hazardous substances.
- Use sealable centrifuge tubes with a biosafety system for materials of risk groups 3 and 4.
- If no biosafety system is used, the device is not micro-biologically tight in the sense of standard EN / IEC 61010-2-020.
- Contact the manufacturer if necessary.

WARNING
Risk of injury and device damage due to the rotor

A loose rotor can cause injuries and device damage.

- Check that the rotor is firmly seated.
- Follow maintenance intervals.

CAUTION
Risk of injury

Long hair and clothing can get caught in the rotor during manual movement.

- Tie long hair back.
- Do not allow garments to hang in the centrifuging chamber.

NOTICE
Damage

Incorrect voltage or frequency.

- Only operate the device according to the specifications on the rating plate.
- Compliance with the instructions for use.

NOTICE
Damage

Aborting the program early can cause damage to the device and samples.

- Do not switch off, perform an emergency release or pull out the mains plug.

3 Device overview

3.1 Technical data

Manufacturer	Andreas Hettich GmbH D-78532 Tuttlingen	
Model	UNIVERSAL 320	
Type	1401	1401-01
Mains voltage ($\pm 10\%$)	200-240 V 1~	100-127 V 1~
Mains frequency	50-60 Hz	50-60 Hz
power consumption	400 VA	400 VA
Power consumption	2.0 A	4.0 A
max. capacity	4 x 200 ml	
max. permissible density	1.2 kg/dm ³	

max. speed (RPM)	16000		
max. acceleration (RCF)	24900		
max. kinetic energy	9800 Nm		
Obligation to perform checks (DGUV Rules 100-500) (valid only in Germany)	No		
Ambient conditions (EN / IEC 61010-1):			
Installation site	indoors only		
Altitude	up to 2000 m above sea level		
Ambient temperature	2 °C to 35 °C		
Humidity	maximum relative humidity 80% for temperatures up to 31 °C, decreasing linearly to 50% relative humidity at 40 °C.		
Overvoltage category (IEC 60364-4-443)	II		
Pollution level	2		
Device protection class	I not suitable for use in potentially explosive atmospheres.		
EMC:			
Emitted EM interference, EM interference immunity	EN / IEC 61326-1 Class B	FCC Class B	
Noise level (rotor-dependent)	≤68 dB(A)		
Dimensions:			
Width	401 mm		
Depth	529 mm		
Altitude	346 mm		
Weight	approx. 31 kg		
Manufacturer	Andreas Hettich GmbH D-78532 Tuttlingen		
Model	UNIVERSAL 320 R		
Type	1406		1406-01
Mains voltage (±10%)	200-240 V 1~	240 V 1~	115-127 V 1~

Mains frequency	50 Hz	60 Hz	60 Hz
power consumption	800 VA		950 VA
Power consumption	4.0 A		8.0 A
Refrigerant	R452A		
max. capacity	4 x 200 ml		
max. permissible density	1.2 kg/dm³		
max. speed (RPM)	16000		
max. acceleration (RCF)	24900		
max. kinetic energy	9800 Nm		
Obligation to perform checks (DGUV Rules 100-500) (valid only in Germany)	No		
Ambient conditions (EN / IEC 61010-1):			
Installation site	indoors only		
Altitude	up to 2000 m above sea level		
Ambient temperature	5 °C to 35 °C		
Humidity	maximum relative humidity 80% for temperatures up to 31 °C, decreasing linearly to 50% relative humidity at 40 °C.		
Overvoltage category (IEC 60364-4-443)	II		
Pollution level	2		
Device protection class	I not suitable for use in potentially explosive atmospheres.		
EMC:			
Emitted EM interference, EM interference immunity	EN / IEC 61326-1 Class B		FCC Class B
Noise level (rotor-dependent)	≤64 dB(A)		
Dimensions:			
Width	407 mm		
Depth	698 mm		

Altitude	346 mm
Weight	approx. 52 kg

Rating plate


Fig. 1: Rating plate


- 1 Item number
- 2 Serial number
- 3 Revision
- 4 Equipment number
- 5 Data matrix code
- 6 any labelling indicating whether medical device or in vitro diagnostic medical device
- 7 Global Trade Item Number (GTIN)
- 8 Date of manufacture
- 9 Serial number
- 10 any EAC mark, CE mark
- 11 Country of manufacture
- 12 Date of manufacture
- 13 Mains frequency
- 14 Maximum kinetic energy
- 15 Maximum permissible density
- 16 Manufacturer's address
- 17 any Coolant circuit pressure
- 18 any Coolant capacity
- 19 any Coolant type
- 20 Revs per minute
- 21 Performance values
- 22 Mains voltage
- 23 any Device designation
- 24 Manufacturer's logo

3.2 Certifications and logos

Certifications

 ISO 9001 Certified Quality Management System www.tuvsud.com/ms-cert	ISO 9001 Quality management system in accordance with ISO 9001
 ISO 14001 Certified Environmental Management System www.tuvsud.com/ms-cert	ISO 14001 Environmental management in accordance with ISO 14001
	EN ISO 13485 Quality management in accordance with ISO 13485

Logos

	Made in Germany Device was developed and made in Germany.
---	--

3.3 European registration

Device conformity

Device conformity according to EU directives.

Single Registration Number

SRN: DE-MF-000010680

Basic-UDI-DI

Basic-UDI-DI	Device assignment
040506740100139R	UNIVERSAL 320 / 320 R (in vitro diagnostic medical device)

3.4 Important labels on the packaging

TOP

This is the correct upright position of the shipping container for transport and/or storage.

FRAGILE GOODS

The contents of the shipping container are fragile, so it must be handled with care.

PROTECT FROM MOISTURE

Protect the transport packaging from moisture and keep it in a dry environment.

TEMPERATURE LIMITATION

The shipping container must be stored, transported and handled within the indicated temperature range (-20 °C to +60 °C).

HUMIDITY LIMITATION

The shipping container must be stored, transported and handled within the indicated air humidity range (10 % to 80 %, non-condensing).

STACK LIMITATION BASED ON QUANTITY

Maximum number of identical packages that may be stacked on the lowest package, "n" standing for the number of packages allowed. The lowest package is not included in "n".

3.5 Important labels on the device

The signs on the device must not be removed or covered, or have anything pasted over them.

Attention, general danger area.

Ensure you read the instructions for commissioning and operation and observe the safety instructions before using the device.

Biohazard warning.

Direction of rotation of the rotor.

The orientation of the arrow indicates the rotor's direction of rotation.

Symbol for the separate collection of electrical and electronic equipment, in accordance with Directive 2012/19/EU (WEEE).

Use in European Union countries, Norway and Switzerland.

3.6 Operating and indicator elements

3.6.1 Control


Fig. 2: Control (device with cooling)


Fig. 3: Control (device without cooling)

3.6.2 Indicator elements

Fig. 4: [START/IMPULS] button

- The button lights up during the centrifugation run for as long as the rotor is not yet at a standstill.

Fig. 5: [STOP/OPEN] button

- The right side of the button lights up when the centrifuge is in ramp-down. The rotor has not yet stopped.
- The left side of the button lights up when the rotor is stationary.
- The light on the left side of the button goes out when the lid is unlocked.

3.6.3 Controls

Fig. 6: [Rotary knob]

- Setting the individual parameters.
Turning anticlockwise decreases the value.
Turning clockwise increases the value.

Fig. 7: [Mains switch]

- Switch the device on and off.

Fig. 8: [Cooling] button

- Start the centrifugation run to pre-cool the rotor (only for centrifuges with cooling).
- The precooling speed is adjustable. The default value is 10,000 RPM.

Fig. 9: [RCF] button

- Toggle between RCF indicator and RPM indicator.
- Relative centrifugal force, RCF.
The RCF is displayed in brackets > <.
- Speed, RPM.

Fig. 10: [SELECT] button

- Selecting the individual parameters.
- Scroll forward in the menus.

Fig. 11: [START/IMPULSE] button

- Start centrifugation run.
- Short-term centrifugation. The centrifugation run takes place as long as the button is being pressed.
- Save entries and changes.

Fig. 12: [STOP/OPEN] button

- End the centrifugation run.
The rotor coasts to a stop with the preselected ramp-down parameter.
- Pressing the button twice triggers the quick stop function.
- Unlock the lid.
- Exit parameter input and the menus.

3.7 Original spare parts

Use only original spare parts from the manufacturer and approved accessories.

3.8 Scope of delivery

The following accessories are supplied with the centrifuge:

- 1 hex key (SW5 x 100)
- 1 grease for the trunnions

- 1 power cable
- 1 user manual
- 1 instruction sheet, transport lock

Rotors and the corresponding accessories are supplied depending on the order.

3.9 Returns

An original Return Material Authorisation (RMA) form from the manufacturer must always be requested for a return. Secure and reliable acceptance and booking in of the goods with the manufacturer is not possible without an original RMA form from the manufacturer. The Return Material Authorisation (RMA) form contains a Declaration of No Objection (UBE), which must be completed in full and enclosed with the return.

If the device and/or accessories are returned to the manufacturer, the complete return shipment must be cleaned and decontaminated by the sender. If returns are not cleaned and/or decontaminated or are insufficiently cleaned and/or decontaminated, this will be performed by the manufacturer and charged to the sender.

The original transport locks must be attached for return shipment, see ➔ *Chapter 4 'Transport and storage' on page 17*. The device must be shipped in its original packaging.

4 Transport and storage

4.1 Transport and storage conditions

Transport conditions

NOTICE

Damage

The device may be damaged if it is not secured during transport.

- Secure the transport locks before transport.
- Observe the transport instructions.

NOTICE

Danger due to condensation in the event of temperature differences

Moisture may damage electrical components.

- Ensure that all surfaces are dry prior to commissioning or maintenance.
- If the temperature changes, wait until the device or component has become acclimatised.
- Prevent moisture from penetrating into sensitive components.
- If moisture forms, switch off the device immediately and allow it to dry properly.

- Before transporting, fasten the transport lock and disconnect the device from the mains socket.
- The transport temperature must be between -20 °C and +60 °C.
- Humidity must not be condensing. Humidity must be between 10% and 80%.

- Be aware of the weight of the device.
- When transporting using a transport aid (e.g., a pallet truck), the transport aid must be able to carry at least 1.6 times the transport weight of the device.
- Secure the device to prevent it tipping over and falling down during transport.
- Never transport the device sideways or upside down.

Storage conditions

- The device must be stored in the original packaging.
- Only store the device in dry rooms.
- The storage temperature must be between -20 °C and +60 °C.
- Humidity must not be condensing. Humidity must be between 10% and 80%.

4.2 Fastening the transport lock

Personnel:

- Trained user


1. ➤ Open the lid.

2. ➤ For UNIVERSAL 320 R:

Check the bellows (3) underneath the motor cover for correct seating.

The bellows (3) must be placed over the edge of the motor cover (1) and over the edge of the centrifuging chamber (2).

3. ➤ Close the lid.

- 1 Motor cover
- 2 Edge of the centrifuging chamber
- 3 Bellows


Fig. 13: Transport lock

- 1 Transport lock
- 2 Screw
- 3 Spacer sleeve

4. ➤ Lay the device down on its right-hand side.

5. ➤ Insert 3 transport locks (1).

6. ➤ Screw in 3 screws (2) with spacer sleeves (3).

5 Commissioning

5.1 Unpacking the centrifuge

**CAUTION**

Danger of crushing due to parts falling out of the transport packaging.

- Keep the device balanced during the unpacking process.
- Only open the packaging at the points provided for this purpose.

**CAUTION**

Risk of injury from lifting heavy loads.

- Provide an adequate number of helpers.
- Note the weight. See ➔ *Chapter 3.1 'Technical data' on page 9.*

**NOTICE**

Damage to the device due to improper lifting.

- Do not lift the centrifuge by the control panel or the control panel holder.

Personnel:

- Trained user

1. ➤ Unscrew the screws on the lid of the wooden packaging and store them.
2. ➤ Remove the lid.
3. ➤ Unscrew the screws on the side panels of the wooden packaging and store them.
4. ➤ Remove the side panels.
5. ➤ Remove padding and strips.
6. ➤ Remove the device and accessories by lifting them up out of the box.
7. ➤ Place the device on a stable and level surface.

5.2 Removing the transport lock

Personnel:

- Trained user

The lid is closed.



Fig. 14: Transport lock

- 1 Transport lock
- 2 Screw
- 3 Spacer sleeve

1. Lay the device down on its right-hand side.
2. Unscrew 3 screws (2) with 3 spacer sleeves (3).
3. Remove 3 transport locks (1).
4. Keep the screws, spacer sleeves and transport locks in a safe place.
5. Place the device upright.
6. Open the lid.
7. For UNIVERSAL 320 R:

Check the bellows (3) underneath the motor cover for correct seating.

The bellows (3) must be placed over the edge of the motor cover (1) and over the edge of the centrifuging chamber (2).

- 1 Motor cover
- 2 Edge of the centrifuging chamber
- 3 Bellows

5.3 Setting up and connecting the centrifuge

Setting up the centrifuge

WARNING

Risk of injury

Due to insufficient clearance from the centrifuge.

- As per EN / IEC 61010-2-020, no persons, hazardous materials or objects may be present within a safety zone of 300 mm around the centrifuge during a centrifugation run.
- Maintain a distance of 300 mm from the ventilation slots and ventilation openings of the centrifuge.
- The centrifuge's ventilation openings must never be blocked.

CAUTION
Risk of crushing and damage

Changes in position caused by vibration may cause the device to fall.

- Place the device on a stable and level surface.
- Select an installation area that is capable of bearing the weight of the device.
- Observe national and local regulations on safety and accident prevention.

NOTICE
Damage

Any deviation from the external temperature conditions will result in damage to the samples and device.

- Comply with the maximum and minimum permissible ambient temperatures.
- Avoid placing the device next to heat sources.
- Do not place the device in direct sunlight.
- Protect the device from frost.
- Maintain the required clearance around the device.

Personnel:

- Trained user

1. ➤ Place the device on a stable and level surface.
2. ➤ Maintain a distance of 300 mm around the device.
3. ➤ Comply with the ambient conditions in the technical data (➔ Chapter 3.1 'Technical data' on page 9).

Connecting the centrifuge

NOTICE
Damage

Property damage caused by unauthorised personnel.

- Do not allow persons without the appropriate authorisation to perform any work on or make changes to devices.
- Only authorised personnel may carry out maintenance and repairs.
- Obtain the manufacturer's approval or guidance before carrying out any work on the device.

**NOTICE****Danger due to condensation in the event of temperature differences**

Moisture may damage electrical components.

- Ensure that all surfaces are dry prior to commissioning or maintenance.
- If the temperature changes, wait until the device or component has become acclimatised.
- Prevent moisture from penetrating into sensitive components.
- If moisture forms, switch off the device immediately and allow it to dry properly.

Personnel:

- Trained user

1. ➤ A type B residual current circuit breaker must be used if the device is additionally protected with a residual current circuit breaker in the building installation.

When using a different type, the residual current circuit breaker may either not switch off the unit if there is a fault on the unit, or it may switch off the unit even though there is no fault on the unit.

2. ➤ Check whether the mains voltage and mains frequency match the specification on the rating plate.
3. ➤ Connect the device to a standard mains socket using the mains cable.

5.4 Switching the centrifuge on and off.

Switching the centrifuge on

Personnel:

- Trained user

- Set the mains switch to */I/*.

- ➡ The buttons flash, depending on the centrifuge type.

The following indicators appear one after the other, depending on the centrifuge type:

- the centrifuge model and program version
- When the lid is closed: 'OPEN OEFFNEN' indicator
- When the lid is open: The last centrifugation data used.

Switching off the centrifuge

The rotor is stationary.

- Set the mains switch to */0/*.

6 Operation

6.1 Opening and closing the lid

Opening the lid

Personnel:

- Trained user

The centrifuge is switched on

The rotor is stationary.

- ➔ Press the *[STOP/OPEN]* button.
 - ➡ The lid unlocks by means of a motor.
- The light on the left side of the *[STOP/OPEN]* button goes out.

Closing the lid

CAUTION

Crushing hazard when closing the lid.

Danger of fingers getting crushed when the closing motor pulls the lid against the seal.

- No parts of the operator's body should be in the hazard zone of the lid when closing the lid.
- To close the lid, press on the lid from above.

NOTICE

Damage to the device caused by the lid slamming.

- Close the lid slowly.
- Do not slam the lid.

*When the left side of the *[STOP/OPEN]* button flashes, press the *[STOP/OPEN]* button so that the motorised lid lock assumes the home position (open).*

Personnel:

- Trained user
- ➔ Close the lid and press the front edge of the lid down gently.
 - ➡ The lid locks using a motor.
- The left side of the *[STOP/OPEN]* button lights up.

6.2 Removing and installing the rotor

Removing the rotor with a clamping nut


Fig. 15: Rotor installation and removal

- 1 Groove
- 2 Motor shaft
- 3 Driver
- 4 Hole

Personnel:

- Trained user
- 1. ➔ Open the lid.
- 2. ➔ Loosen the rotor clamping nut using the supplied spanner.
 - ➡ After passing the working point for lifting the rotor, the rotor detaches from the cone of the motor shaft (2).
- 3. ➔ Turn the clamping nut until the rotor can be lifted off the motor shaft.
- 4. ➔ Remove the rotor.

Installing the rotor with a clamping nut

Personnel:

- Trained user

The lid is open.

1. ➤ Clean the motor shaft (2) and rotor hole (4).
2. ➤ Lightly grease the motor shaft (2), see ➔ Chapter 8.2 'Instructions for cleaning and disinfection' on page 37.
3. ➤ Place the rotor vertically on the motor shaft (2).
The driver (3) of the motor shaft must be in the groove (1) of the rotor. The orientation of the groove is marked on the rotor.
4. ➤ Hand-tighten the rotor clamping nut using the supplied spanner.
5. ➤ Check that the rotor is firmly seated.

Removing the rotor without a clamping nut

Removing the rotor

Personnel:

- Trained user

- Hold the rotor by the rotary handle (1) of the lid and lift it off the hub (2).



Fig. 16: Rotor installation and removal

- 1 Rotary handle
- 2 Hub

Removing the hub

1. ➤ Open the lid.
2. ➤ Unscrew the clamping nut.
➔ After passing the working point for lifting the rotor, the hub (1) detaches from the cone of the motor (2).
3. ➤ Remove the hub.

Fig. 17: Hub installation and removal

- 1 Hub
- 2 Motor shaft
- 3 Driver
- 4 Groove

Installing the rotor without a clamping nut

Installing the hub

Personnel:

- Trained user

1. ➤ Open the lid.
2. ➤ Clean the motor shaft (2) and rotor hole.
3. ➤ Lightly grease the motor shaft (2), see ➡ Chapter 8.2 'Instructions for cleaning and disinfection' on page 37.
4. ➤ Place the hub (1) vertically on the motor shaft (2).
The driver (3) of the motor shaft must be in the groove (4) of the hub.
Check that the hub is firmly seated.
5. ➤ Hand-tighten the clamping nut of the hub using the supplied hex key.
6. ➤ Check that the hub is firmly seated.

Installing the rotor

1. ➤ Clean the hub (2).
2. ➤ Lift the rotor by the rotary handle and place it vertically on the hub (2).
3. ➤ Push the rotor down as far as it will go.

6.3 Inserting and removing buckets

Inserting buckets

NOTICE



Damage to the device due to imbalances caused by incorrect loading of the rotor.

- Load all swing-out rotor locations with the same buckets.

Buckets marked with the number of the rotor location may only be used there.

Buckets marked with a set number may only be used together.


1. ➤ Check that the rotor is firmly seated.
2. ➤ Grease the trunnions (3).
3. ➤ Insert the bucket (1) into the rotor from above. The trunnions (3) must be in the grooves (2).
4. ➤ Push the bucket (1) down as far as it will go.

Removing the bucket

- Pull the bucket (1) vertically upwards out of the rotor.

6.4 Inserting and removing adapters

Inserting

the adapter

➔ Insert the adapter (1) vertically into the bucket (2) from above.

removing

➔ Remove the adapter (1) vertically upwards out of the bucket (2).

Adapter with positioning pin


Fig. 18: Adapter with positioning pin

- 1 Inserting
- 2 Positioning pin
- 3 Inserting

the groove

➔ Insert the adapter (1) into the bucket
The positioning pin (2) must be in the groove (3) of the bucket.

removing

➔ Remove the adapter (1) vertically upwards out of the bucket.

6.5 Loading

Filling centrifuge tubes

WARNING

Risk of injury from contaminated sample material.

Contaminated sample material escapes from the sample tube during centrifugation.

- Use centrifuge tubes with special screw caps for hazardous substances.
- For risk group 3 and 4 materials, use a biosafety system in addition to the sealable centrifuge tubes (see WHO's 'Laboratory Biosafety Manual').

NOTICE

Damage to the device due to highly corrosive substances.

Highly corrosive substances may impair the mechanical strength of rotors, buckets and accessories.

- Do not centrifuge highly corrosive substances.

Standard glass centrifuge tubes can be loaded up to RCF 4000 (DIN 58970 part 2).

Personnel:

- Trained user


➔ Fill centrifuge tubes outside the centrifuge.

The maximum capacity of the centrifuge tubes specified by the manufacturer must not be exceeded.

With angle rotors, the centrifuge tubes must only be filled to the extent that no liquid can be ejected from the tubes during the centrifugation run.

It must be ensured that there is a uniform fill level in the tubes in order to keep the weight differences in the centrifuge tubes as low as possible.

Loading swing-out rotors

Personnel:

- Trained user

1. ➔ Check that the rotor is firmly seated.

2. ➔ The centrifuge tubes must be distributed symmetrically across all rotor locations.

The weight of the permissible filling capacity is indicated on each rotor. The weight must not be exceeded.

No liquid must be allowed to enter the buckets and the centrifuging chamber when loading the buckets and swinging them out during the centrifugation run.

For containers with rubber inserts, there must always be the same number of rubber inserts under the centrifuge tubes.


All rotor locations must be filled with the same buckets. Certain buckets are marked with the number of the rotor location. The buckets must only be inserted in the corresponding rotor location.

Buckets marked with a set number (e.g. S001/4) may only be used in the set.

Loading the angle rotors

Personnel:

- Trained user

1. ➤ Check that the rotor is firmly seated.
2. ➤ The centrifuge tubes must be distributed evenly over all locations on the rotor.

No liquid must be allowed to enter the rotor and the centrifuging chamber when loading the rotor.

With rotors, the centrifuge tubes must only be filled to the extent that no liquid can be ejected from the tubes during the centrifugation run.

The weight of the permissible filling capacity is indicated on each rotor. The weight must not be exceeded.

6.6 Opening and closing the biosafety system

6.6.1 Explanation

The user must take appropriate actions when centrifuging hazardous substances or mixtures of substances that are toxic, radioactive or contaminated with pathogenic microorganisms.

Centrifuge tubes with special screw caps for hazardous substances must always be used.

For materials of risk group 3 and 4, a biosafety system must be used in addition to the sealable centrifuge tubes (see the World Health Organisation's "Laboratory Biosafety Manual").

In a biosafety system, a bioseal (sealing ring) prevents droplets and aerosols from escaping.

If the bucket of a biosafety system is used without the lid, the sealing ring must be removed from the bucket to prevent damage to the sealing ring during the centrifugation run.

Damaged biosafety systems are no longer microbiologically tight.

If no biosafety system is used, a centrifuge is not microbiologically tight in the sense of the EN / IEC 61010-2-020 standard.

Storage of biosafety systems

Biosafety systems must only be stored with the lid open to avoid damage to the sealing rings during storage.

6.6.2 Lid with screw cap and hole


Fig. 19: Biosafety system

- 1 Rotary handle
- 2 Rotor
- 3 Lid
- 4 Key

Closing

1. ➤ Place the lid (3) centrally on the rotor (2).
2. ➤ Insert the supplied key (4) into the hole in the rotary handle (1).
3. ➤ Turn the lid (3) at the key (4) clockwise until it is tightly closed.

Opening

1. ➤ Insert the supplied key (4) into the hole in the rotary handle (1).
2. ➤ Turn the lid (3) at the key (4) anticlockwise until it is open.
3. ➤ Remove the lid (3) from the rotor (2).

6.6.3 Lid with bracket and spring-type lock


Fig. 20: Biosafety system

- 1 Bracket position "open"
- 2 Bracket opening operations
- 3 Bracket position "close"
- 4 Bracket carrying position

Closing

1. ➤ Swivel the bracket to the "open" position (1).
The arrows of the labelling must point downwards so that the word "open" is legible.
2. ➤ Place the lid centrally on the bucket.
The two pins of the lid must be in the two openings of the bracket (2).
3. ➤ Swivel the bracket to the "close" position (3).
The arrows of the labelling must point downwards so that the word "close" is legible.
The bracket must rest on the bucket so that the bucket can swing out during the centrifugation run.

4. ➤ For transport or when inserting and removing the bucket, swing the bucket into the carrying position (4) and hold the bucket by the bracket.
 - The tightness of the biosafety system is also guaranteed in the carrying position.

Do not rock the biosafety system back and forth during transport, otherwise leak-tightness is no longer ensured.

Opening

1. ➤ Swivel the bracket to the "open" position (1).
The arrows of the labelling must point downwards so that the word "open" is legible.
2. ➤ Remove the lid from the bucket.

6.6.4 Lid with screw cap


Fig. 21: Biosafety system

- 1 Lid
- 2 Bucket

Closing

1. ➤ Place the lid (1) centrally on the bucket (2).
2. ➤ Turn the lid (1) clockwise until it is tightly closed.

Opening

1. ➤ Turn the lid (1) anticlockwise until it is open.
2. ➤ Remove the lid (1) from the bucket (2).

6.7 Centrifugation

6.7.1 Centrifugation in continuous operation

Personnel:

- Trained user

1. ➤ Set minutes and seconds at '∞' or retrieve a continuously running program.
2. ➤ Press the [START/IMPULS] button.
 - The centrifugation run is started.
 - The [START/IMPULSE] button lights up during the centrifugation run.
 - The timing starts at '00:00'.
 - The rotor speed or the RCF value, the temperature in the centrifuging chamber (only for centrifuges with cooling) and the elapsed time are displayed during the centrifugation run.
3. ➤ Press the [STOP/OPEN] button to cancel the centrifugation run.
 - Ramp-down takes place with the set brake level. The brake level is displayed.
 - An audible signal sounds when the rotor comes to a standstill.
 - 'OPEN OEFFNEN' is displayed.

6.7.2 Centrifugation with time preselection

Personnel:

- Trained user

1. ➤ Set centrifugation parameters or retrieve a program.
2. ➤ Press the *[START/IMPULS]* button.
 - ➡ The centrifugation run is started.

The *[START]* button lights up during the centrifugation run.

The rotor speed or the RCF value, the temperature in the centrifuging chamber (only for centrifuges with cooling) and the remaining time are displayed during the centrifugation run.
3. ➤ Ramp-down takes place with the selected brake level after the time has elapsed or if the centrifugation run is cancelled.
 - ➡ The brake level is displayed.

An audible signal sounds when the rotor comes to a standstill.

'OPEN OEFFNEN' is displayed.

The right side of the *[STOP/OPEN]* button lights up when the centrifuge is in ramp-down.

The left side of the *[STOP/OPEN]* button lights up when the rotor is at a standstill.

The light on the *[START/IMPULS]* button and the right side of the *[STOP/OPEN]* button go out.

6.7.3 Short-term centrifugation

Personnel:

- Trained user

1. ➤ Press and hold the *[START/IMPULS]* button.
 - ➡ The *[START/IMPULS]* button lights up during the centrifugation run.

Timing starts at 00:00.

The rotor speed or the RCF value, the temperature in the centrifuging chamber (only for centrifuges with cooling) and the elapsed time are displayed during the centrifugation run.
2. ➤ Release the *[START/IMPULSE]* button to end the centrifugation run.
 - ➡ Ramp-down takes place with the set brake level. The brake level is displayed.

An audible signal sounds when the rotor comes to a standstill.

'OPEN OEFFNEN' is displayed.

6.8 Quick stop function

Personnel:

- Trained user

- Press the *[STOP/OPEN]* button twice.
 - ➡ Ramp-down with brake level "9" (shortest ramp-down time) is displayed and executed.

7 Software operation

7.1 Centrifugation parameters

7.1.1 Relative centrifugal force, RCF

The relative centrifugal force RCF is dependent on the speed and the centrifuging radius.

The relative centrifugal force RCF is stated as a multiple of the acceleration due to gravity (g).

The relative centrifugal force RCF is a dimensionless numerical value and is used to compare the separation and sedimentation performance.

$$RCF = \left(\frac{RPM}{1000} \right)^2 * r * 1,118$$

$$RPM = \sqrt{\frac{RCF}{r * 1,118}} * 1000$$

RCF = relative centrifugal force

RPM = Speed

r = centrifuging radius in mm = distance from the centre of the axis of rotation to the bottom of the centrifuge tube.

7.1.2 Centrifugation of substances or mixtures of substances with a density higher than 1.2 kg/dm³

The density of the substances or mixtures of substances must not exceed 1.2 kg/dm³ during centrifugation at maximum speed. The maximum permissible speed must be reduced for substances or substance mixtures with a higher density. The permissible speed can be calculated using the following formula:

$$\text{Reduced speed } (n_{red}) = \sqrt{\frac{1,2}{\text{Greater density [kg/dm}^3]}} * \text{maximum speed [RPM]}$$

For example: Maximum speed 4000 RPM, density 1.6 kg/dm³

$$n_{red} = \sqrt{\frac{1,2(\text{kg/dm}^3)}{1,6(\text{kg/dm}^3)}} * 4000 \text{ RPM} = 3464 \text{ RPM}$$

If, in exceptional cases, the maximum load indicated on the bucket is exceeded, the speed must also be reduced. The permissible speed can be calculated using the following formula:

$$\text{Reduced speed } (n_{red}) = \sqrt{\frac{\text{maximum load [g]}}{\text{actual load [g]}}} * \text{maximum speed [RPM]}$$

For example: Maximum speed 4000 RPM, maximum load 300 g, actual load 350 g

$$n_{red} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} * 4000 \text{ RPM} = 3703 \text{ RPM}$$

Please contact the manufacturer if you are not sure.

7.2 Programming

7.2.1 Write protection for programs

The programs can be protected to prevent unintentional changes.

When the rotor is at a standstill, write protection can be enabled or disabled as follows:

1.  Press and hold the *[SELECT]* button.
 - ➡ 'SOUND/BELL' is displayed after 8 seconds.
2.  Press the *[SELECT]* button.
 - ➡ 'LOCK' is displayed.
3.  Use *[Rotary knob]* to set 'OFF' or 'ON'.
 - OFF = The programs are not write-protected
 - ON = The programs are write-protected
4.  Press the *[START/IMPULS]* button.
 - ➡ The setting is stored.
 - If ON is set: '*** lock ***' is displayed briefly.
 - If OFF is set: '*** ok ***' is displayed briefly.

7.2.2 Opening or loading programs

1.  Use the *[SELECT]* button to select the 'PROG RCL' parameter.
2.  Use the *[Rotary knob]* to set the desired program location.
3.  Press the *[START/IMPULS]* button.
 - ➡ '*** ok ***' is displayed briefly.
 - The centrifugation data of the desired program location is displayed
4.  To check the parameters: Press the *[SELECT]* button several times.
5.  To exit the parameter indicator: Press the *[OPEN/STOP]* button or do not press any button for 8 seconds.

7.2.3 Entering or changing programs

1.  Retrieve program.
2.  If required: Press the *[RCF]* button to toggle between RPM and RCF indicator ('>' '<').
3.  If required: Press the *[SELECT]* button to select the desired parameter and set it with the *[Rotary knob]*.
 - The parameters t/min and t/sec must be set to 0 using the *[Rotary knob]* to set continuous operation. Continuous operation is shown in the indicator with '∞'.
4.  Use the *[SELECT]* button to select the 'PROG STO' parameter.
5.  Use the *[Rotary knob]* to set the desired program location.
6.  Press the *[START/IMPULS]* button.
 - ➡ Settings are stored in the desired program location.
 - '*** ok ***' is displayed briefly.
 - The settings are always stored in program location # if the *[START/IMPULS]* button is pressed without the 'PROG STO' parameter being selected.

7.3 Rotor detection

- Rotor detection is performed after starting a centrifugation run.
- If the rotor has been changed, the centrifugation run is cancelled after rotor detection. The rotor code (rot) is displayed.
- If the maximum speed of the rotor used is less than the set speed, the speed is limited to the maximum rotor speed.

7.4 Cooling (for centrifuges with cooling)

7.4.1 Instructions, cooling

The temperature setpoint can be adjusted from -20 °C to +40 °C.
The lowest achievable temperature is rotor dependent.

7.4.2 Standby cooling

After a centrifugation run, standby cooling is delayed and the display shows '*Deckel entriegelt*'.

The delay time is adjustable from 1 to 5 minutes, in 1-minute increments. It is preset to 1 minute.

- The rotor is stationary.
 - The lid is open
1.  Press and hold the *[Cooling]* button.
 - ➡ '*t/min = X*' is displayed after 8 seconds.
 2.  Use the *[Rotary knob]* to set the delay time.
 3.  Press the *[START/IMPULS]* button.
 - ➡ The setting is stored.
 - *** ok *** is displayed briefly.
 4.  Press the *[STOP/OPEN]* button twice or wait 8 seconds to exit the menu.

7.4.3 Precooling the rotor

Starting

The rotor is stationary.

1.  Press the *[Cooling]* button.
2.  Press the *[STOP/OPEN]* button.
 - ➡ Precooling of the rotor is terminated.
 - Ramp-down takes place with the selected brake level.
 - The brake level is displayed.

Set

The precooling speed is adjustable from 500 RPM up to the maximum rotor speed in increments of 10 RPM. It is preset to 10000 RPM.

- The rotor is stationary.
 - The lid is open.
1.  Press and hold the *[Cooling]* button.
 - ➡ '*t/min = X*' is displayed after 8 seconds.
 2.  Press the *[Cooling]* button.
 - ➡ Precooling speed '*RPM = XXXX*' is displayed.
 3.  Use the *[Rotary knob]* to set the precooling speed.

4. ➤ Press the *[START/IMPULS]* button.
 - The setting is stored.
 - *** ok **** is displayed briefly.
5. ➤ Press the *[STOP/OPEN]* button twice or wait 8 seconds to exit the menu.

7.5 Machine Menu

7.5.1 Querying system information

Parameter query

The rotor is stationary.

1. ➤ Press and hold the *[SELECT]* button for 8 seconds.
 - *'SOUND/BELL'* is displayed.
2. ➤ Press the *[SELECT]* button repeatedly until *'FU/CCI -S.'* is displayed.
Program version for the frequency inverter
3. ➤ Press the *[SELECT]* button repeatedly until *'HOURS'* is displayed.
Internal operating hours (the time during which the centrifuge was switched on)
4. ➤ Turn to the right with the *[Rotary knob]*.
 - *'STARTS'* is displayed.
 - Number of centrifugation runs
5. ➤ Turn to the right with the *[Rotary knob]*.
 - *'ROTORCHG1'* is displayed.
 - Internal operating hour of the last rotor change
6. ➤ Turn to the right with the *[Rotary knob]*.
 - *'ROTORCHG2'* is displayed.
 - Internal operating hour of the penultimate rotor change
7. ➤ Turn to the right with the *[Rotary knob]*.
 - *'OPhoursCHG'* is displayed.
 - Internal operating hour of the last operating hours change
8. ➤ Turn to the right with the *[Rotary knob]*.
 - *'IMBALCHG'* is displayed.
 - Internal operating hour of the last imbalance cut-off change
9. ➤ Turn to the right with the *[Rotary knob]*.
 - *'OffsetCHG'* is displayed.
 - Internal operating hour of the last offset adjustment
10. ➤ Press the *[STOP/OPEN]* button to exit the menu.

7.5.2 Querying operating hours

The rotor is stationary.

1. ➤ Press and hold the *[SELECT]* button.
 - *'SOUND/BELL'* is displayed after 8 seconds.
2. ➤ Press the *[SELECT]* button repeatedly until *'CONTROL.'* is displayed.
 - *'CONTROL.'* and the operating hours are displayed.
3. ➤ Press the *[STOP/OPEN]* button to exit the menu.

7.5.3 Audible signal

7.5.3.1 General

The audible signal sounds:

- after a problem occurs in the 2 s interval.
- after completion of the centrifugation run and rotor standstill in the 30 s interval.

Opening the lid or pressing any button stops the audible signal.

7.5.3.2 Setting an audible signal

1. ➔ Press and hold the *[SELECT]* button.
 - ➔ 'SOUND / BELL ON' or 'SOUND / BELL OFF' is displayed after 8 seconds.
2. ➔ Use *[Rotary knob]* to set 'OFF' or 'ON'.
 - OFF = audible signal disabled
 - ON = audible signal enabled
3. ➔ Press the *[START/IMPULS]* button.
 - ➔ The setting is stored.
 - **** ok **** is displayed briefly.

8 Cleaning and care

8.1 Overview table

Chap.	Task to execute	if required	daily	weekly	Annually	Page
8	Cleaning and care					36
8.3	Cleaning					37
8.3	Cleaning the device		X			37
8.3	Cleaning biosafety systems			X		38
8.3	Cleaning accessories			X		38
8.4	Disinfection					38
8.4	Disinfecting the device	X				38
8.4	Disinfecting the accessories	X				38
8.5	Maintenance					39
8.5	Greasing the rubber seal of the centrifuging chamber			X		39
8.5	Greasing the rubber seal of the biosafety system			X		39
8.5	Trunnion greasing			X		39

Chap.	Task to execute	if required	daily	weekly	Annually	Page
8.5	Checking the accessories			X		39
8.5	Checking the biosafety system			X		39
8.5	Centrifuging chamber damage inspection				X	39
8.5	Greasing the motor shaft				X	39
8.5	Accessories with a limited service life	X				40
8.5	Replacing centrifuge tubes	X				40

8.2 Instructions for cleaning and disinfection

DANGER

Risk of contamination

Inadequate cleaning or failure to observe the cleaning instructions can lead to contamination risks.

- Observe national and local regulations on safety and accident prevention.
- Observe cleaning instructions.
- Wear personal protective equipment when cleaning the device.
- Comply with laboratory regulations (e.g. TRBAs, the German Protection against Infection Act, hygiene plan) for handling biological agents.

- The device and its accessories must not be cleaned in dishwashers.
- Only perform hand cleaning and liquid disinfection.
- The water temperature must not exceed 25 °C.
- To prevent any corrosion due to use of detergents or disinfectants, it is essential to follow the special application instructions provided by the manufacturers of the detergent or disinfectant.

Disinfectant:

- Surface disinfectant (not disinfectant for hands or instruments)
- Ethanol as the sole active substance.
Do not use an ethanol-propanol mixture to disinfect the viewing window in the lid of the device.
- Concentration is not less than 30 %
- pH: 6 – 8
- Non-corrosive

8.3 Cleaning

Cleaning the device

1. ➤ Open the lid.
2. ➤ Switch off the device and disconnect it from the power supply.
3. ➤ Remove accessories.

4. ➤ Clean the centrifuge housing and the centrifuging chamber with soap or a mild detergent and a damp cloth.
5. ➤ Remove any detergent residues with a damp cloth after using detergents.
6. ➤ The surfaces must be dried immediately after cleaning.
7. ➤ Dry the centrifuging chamber with an absorbent cloth if condensation forms.

Cleaning biosafety systems

1. ➤ Clean the biosafety system using the detergent and a damp cloth.
2. ➤ Remove any detergent residues with a damp cloth after using detergents.
3. ➤ Dry the accessories immediately after cleaning using a lint-free cloth and oil-free compressed air. Dry all cavities completely using oil-free compressed air.

Cleaning accessories

1. ➤ Clean the accessories using the detergent and a damp cloth.
2. ➤ Remove any detergent residues with a damp cloth after using detergents.
3. ➤ Dry the accessories immediately after cleaning using a lint-free cloth and oil-free compressed air. Dry all cavities completely using oil-free compressed air.

8.4 Disinfection

Disinfection must always be preceded by cleaning of the components concerned.

See ➔ Chapter 8.3 'Cleaning' on page 37

Disinfectant concentration and application time according to the manufacturer's instructions.

Disinfecting the device

CAUTION

Risk of injury due to ingress of water or other liquids.

- Protect the device against external liquids.
- Do not disinfect the device using spray.

1. ➤ Open the lid.
2. ➤ Switch off the device and disconnect it from the power supply.
3. ➤ Remove accessories.
4. ➤ Clean the housing and centrifuging chamber using disinfectant.
5. ➤ Remove any disinfectant residues with a damp cloth after using disinfectants.
6. ➤ The surfaces must be dried immediately after cleaning.

Disinfecting the accessories

1. ➤ Disinfect the accessories using the disinfectant.
2. ➤ Wet all cavities with bubble-free disinfectant.

3. ➤ Remove the disinfectant residues or leave them to dry after using disinfectants.

Autoclaving

The following accessories may be autoclaved at 121 °C / 250 °F (20 min):

- Swing-out rotors
- Aluminium angle rotors
- Metal buckets
- Lid with bioseal
- Inserting

No statement can be made about the resulting degree of sterility.

The lids of the rotors and bucket must be removed before autoclaving.

Autoclaving accelerates the ageing of materials. It may cause changes in colour. After autoclaving, the rotors and accessories are to be visually inspected for damage and any damaged parts are to be replaced immediately.

The sealing ring in question is to be replaced if there are signs of cracking, embrittlement or wear. For lids with non-replaceable sealing rings, the whole lid must be replaced.

The sealing rings must be replaced after autoclaving to ensure the tightness of the biosafety systems.

8.5 Maintenance

Greasing the rubber seal of the centrifuging chamber

- Rub the sealing ring lightly with a rubber care product.

Greasing the rubber seal of the biosafety system

- Rub the sealing ring lightly with a rubber care product.

Trunnion greasing

1. ➤ Remove accessories.
2. ➤ Clean the trunnions.
3. ➤ Remove any detergent residues with a damp cloth after using detergents.
4. ➤ Grease the trunnions and suspension with Hettich Tubenfett 4051.
5. ➤ Excess grease in the centrifuging chamber must be removed.

Checking the accessories

1. ➤ The accessories shall be checked for wear and corrosion damage.
2. ➤ Check that the rotor is firmly seated.

Checking the biosafety system

1. ➤ Visually check all parts of the biosafety system for damage.
2. ➤ Check the correct installation position of the sealing ring(s) of the biosafety system.
3. ➤ Replace the damaged parts of the biosafety system.
4. ➤ Replace the sealing ring in question immediately if there are signs of cracking, embrittlement or wear. For lids with non-replaceable sealing rings, the whole lid must be replaced.

Centrifuging chamber damage inspection

- Check the centrifuging chamber for damage.

Greasing the motor shaft

1. ➤ Remove accessories.

2. ➤ Clean the motor shaft.
3. ➤ Remove any detergent residues with a damp cloth after using detergents.
4. ➤ Grease the motor shaft with Hettich Tubenfett 4051.
5. ➤ Excess grease in the centrifuging chamber must be removed.

Accessories with a limited service life

The use of certain accessories is time-limited. For safety reasons, the accessories must no longer be used when either the maximum number of permissible run cycles marked on them or the expiry date marked on them has been reached.

- The maximum permissible number of run cycles or the expiry date can be seen on the accessories.

Replacing centrifuge tubes

CAUTION

Risk of injury from broken glass.

Broken glass may cause glass splinters and contaminated liquids to be found inside the centrifuge.

- Wear cut-resistant gloves.
- Wear protective goggles and a face mask.

Broken parts of the tube, glass splinters and spilled centrifuge material must be removed completely in the event of leakage or if a centrifuge tube breaks. Glass splinters that are not removed will cause further glass breakage.

The rubber inserts and the plastic sleeves of the rotors must be replaced after a glass breakage.

Disinfection must be carried out if the material is infectious.

9 Troubleshooting

9.1 Fault description

Customer service must be notified if the fault cannot be rectified based on the fault table. State the centrifuge type and serial number. Both numbers can be seen on the type plate of the centrifuge.

* Error number does not appear on the display.

Fault description	Cause	Remedy
no display	No power. Overcurrent protection fuse has tripped.	■ Check the supply voltage. ■ Set the mains switch to <i>/I/</i>.
TACHO - ERROR 1, 2, 96	Tacho defective. Motor, electronics defective.	■ Open the lid. ■ Set the mains switch to <i>/O/</i>. ■ Wait at least 10 seconds. ■ Turn the rotor vigorously by hand. ■ Set the mains switch to <i>/I/</i>. The rotor must rotate while switching on.
IMBALANCE 3*	The rotor is unevenly loaded.	■ Open the lid. ■ Check the loading of the rotor. ■ Repeat the centrifugation run.

Fault description	Cause	Remedy
CONTROL - ERROR 4, 6	Lid lock error.	Perform a MAINS RESET.
CONTROL - ERROR 8	Lid lock error	- Open the lid. - Set the mains switch to [0]. - Wait at least 10 seconds. - Turn the rotor vigorously by hand. - Set the mains switch to [I]. The rotor must rotate while switching on.
N > MAX 5	Overspeed.	Perform a MAINS RESET.
N < MIN 13	Underspeed.	Perform a MAINS RESET.
MAINS INTERRUPT 11*	Loss of mains power during the centrifugation run. The centrifugation run was not completed.	- Open the lid. - Press the [START/IMPULS] button. - If required: Repeat the centrifugation run.
ROTORCODE 10.1, 10.2	Rotor coding error.	Open the lid.
CONTROL-ERROR 21, 22, 25, 27, 29	Error/defect in electronics.	Perform a MAINS RESET.
CONTROL-ERROR 23	Error/defect in control panel.	Perform a MAINS RESET.
SER I/O-ERROR 30, 31, 33, 36	Error/defect in electronics.	Perform a MAINS RESET.
°C * -ERROR 51-53, 55	Error/defect in electronics.	Perform a MAINS RESET.
FU/CCI-ERROR 60-64, 67, 68, 82-86	Error/defect in electronics/motor.	Perform a MAINS RESET.
SYNC-ERROR 90	Error/defect in electronics.	Perform a MAINS RESET.
SENSOR-ERROR 91-93	Error/defect in imbalance sensor.	Perform a MAINS RESET.
KEYBOARD-ERROR	Error/defect in control panel.	Perform a MAINS RESET.
NO ROTOR	No rotor installed.	Open the lid and install the rotor.
N > ROTOR MAX	Speed in the selected program greater than the maximum rotor speed.	Check and correct the speed.
N > ROTOR MAX	The rotor has been changed. The built-in rotor has a higher maximum speed than the previously used rotor, and it has not yet been detected by the rotor detection function.	Set a speed up to the maximum speed of the previously used rotor. Press the [START/IMPULS] button to perform rotor detection.
 The left half of the display lights up.	-	Notify customer service.

9.2 Performing a MAINS RESET

1. ➤ Set the mains switch to [0].
2. ➤ Wait 10 seconds.
3. ➤ Set the mains switch to [I].

9.3 Emergency release

The lid cannot be unlocked by the motor in the event of a power failure. Emergency unlocking by hand must be performed.

! WARNING

Risk of electric shock due to maintenance and servicing work on live device.

- Disconnect the device from the mains before carrying out repairs and maintenance.

WARNING

Danger of cutting and crushing due to moving rotor.

- Do not open the lid until the rotor has stopped.


Fig. 22: Emergency release

1 Hole

Personnel:

- Trained user

1. ➤ Look through the window in the lid to ensure that the rotor is stationary.
2. ➤ Insert the hex key horizontally into the hole (1) and turn clockwise until the lid opens.
3. ➤ Remove the hex key from the hole (1).
4. ➤ Check whether the left side of the [STOP/OPEN] button flashes when power is restored.

When the left side of the [STOP/OPEN] button flashes, press the [STOP/OPEN] button so that the motorised lid lock assumes the home position (open) again.

10 Disposal

10.1 General instructions

The device can be disposed of via the manufacturer.

A Return Material Authorisation (RMA) form must always be requested for a return.

If necessary, contact the Technical Service Department of the manufacturer:

- **Andreas Hettich GmbH**
- Föhrenstraße 12
- 78532 Tuttlingen, Germany
- Phone: +49 7461 705 1400
- E-Mail: service@hettichlab.com

! WARNING

Risk of soiling and contamination

Damage to the environment and health due to improper disposal.

- Observe national and local environmental protection and disposal regulations to ensure proper disposal or recycling. Separate metals, non-metals, composite materials and auxiliary materials by type and dispose of them in an eco-friendly manner.
- The device may only be disassembled and disposed of by trained and authorised service technicians.

The device is intended for the commercial sector ("Business to Business" - B2B).

According to Directive 2012/19/EU, the devices may no longer be disposed of with household waste.

The devices are assigned to the following groups according to the Stiftung Elektro-Altgeräte Register (EAR (German foundation under civil law)):

- Group 1 (heat exchangers)
- Group 4 (large devices)

The crossed-out wheellie bin symbol indicates that the device must not be disposed of with household waste. Regulations governing disposal of such devices may differ in individual countries. If necessary, contact the supplier.

Fig. 23: Household waste ban

11 Index

A		
Accessories.	16	
biosafety systems.	38	
disinfection.	38	
the rotor.	39	
with limited service life.	40	
Audible signal		
an audible signal.	36	
Autoclaving.	39	
C		
Care		
Intervals.	36	
Centrifugation		
in continuous operation.	30	
with higher substance density.	32	
with time preselection.	31	
Centrifuging chamber		
inspection.	39	
Certifications.	13	
Checking		
the biosafety system.	39	
Cleaning.	37	
biosafety systems.	38	
Cleaning and disinfection		
Instructions.	37	
Connecting the centrifuge.	21	
Continuous operation.	30	
D		
Device		
biosafety systems.	37	
disinfection.	38	
Disinfection.	38	
Disposal.	42	
E		
Error messages.	40	
F		
Fastening		
removal.	19	
the transport lock.	18	
Filling.	26	
Foreseeable misuse.	6	
G		
General safety instructions.	7	
Greasing		
the rubber seal.	39	
I		
Intended use.	5	
L		
Labels		
on the device.	14	
on the packaging.	14	
Lid		
the lid.	22, 23	
Loading.	26	
Logos.	13	
M		
MAINS RESET.	42	
Maintenance.	39	
Intervals.	36	
N		
Non-intended use.	6	
O		
Opening		
for programs.	33	
programs.	33	
Operator's responsibility.	7	
Original spare parts.	16	
P		
Personal protective equipment.	6	
Personnel qualifications.	6	
Personnel training.	7	
Protective equipment.	6	
Q		
Querying		
system information.	35	
R		
Rating plate.	12	
Relative centrifugal force		
RCF.	32	
Removing		
the rotor.	23, 27	
Replacing		
centrifuge tubes.	40	
Return.	17	
Rotor		
the rotor.	27	
Rotor detection.	34	
S		
Safety instructions.	7	
Scope of delivery.	16	
Setting up the centrifuge.	20	
Short-term centrifugation.	31	
Spare parts.	16	
Storage conditions.	18	
Switching off.	22	
Switching the unit on.	22	
Symbols.	5	
T		
Transport condition.	17	
Troubleshooting.	40	
Trunnion		
greasing.	39	

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1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)



2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)



3) Observe the tube manufacturer's instructions.

4) can not be centrifugated when plug is attached, Scal. 10µl-300µl, 15ml, 30ml

5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

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





- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.
- 5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".



1324	1490 + 1492								
Ausschwingrotor 4-fach / Swing out rotor 4-times  	 mit Bioabdichtung / with bio-containment ⁵⁾								
	1383					1384	1396	1457	
									
									
	Kapazität / capacity	ml	2,6 – 2,9	4,9	4,5 - 5	1 -5	4 –7	50	85
Maße / dimensions Ø x L	mm	13 x 65	13 x 90	11 x 92	13 x 75	13 x 100	29 x 115	38 x 106	8 x 66
Anzahl p. Rotor / number p. rotor		20	20	20	20	20	4	4	28
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF ³⁾		3192	3192	3192	3129	3129	3328	3260	3215
Radius / radius	mm	141	141	141	141	141	147	144	142
 9 (97%)	sec	27							
 9	sec	30							
Temperatur / temperature	°C ¹⁾	- 6							
Probenerwärmung/Sample temp. rise	K ²⁾	10							

1324	1490 + 1492								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°	 mit Bioabdichtung / with bio-containment 5)								
	1459	4416	4417	6311	6318	1356	0761		
									
									
	Kapazität / capacity	ml	4 – 5,5	7,5 – 8,2	50	30	12	50	15
Maße / dimensions Ø x L	mm	15 x 75	15 x 92	29 x 107	26 x 95	17 x 100	29 x 115	17 x 120	44 x 100
Anzahl p. Rotor / number p. rotor		16	16	4	4	4	4	12	4
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF 3)		3215	3215	3283	3056	3328	3328	3328	3192
Radius / radius	mm	142	142	145	135	147	147	147	141
 9 (97%)	sec	27							
 9	sec	30							
Temperatur / temperature	°C 1)	- 6							
Probenerwärmung/Sample temp. rise	K 2)	10							

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






- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.
- 5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 7) not possible to close the lid



1324	1398								
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°									
	1482A + 2 x 0716				1482A				
	 +  (2 x)								
									
	Kapazität / capacity	ml	2,6 – 3,4	4 – 5,5	9 – 10	10	12	4 - 7	5 – 10
Maße / dimensions $\varnothing \times L$	mm	13 x 65	15 x 75	16 x 92	15 x 102	17 x 100	16 x 75	16 x 100	14 x 100
Anzahl p. Rotor / number p. rotor		16	16	16	16	16	16	16	16
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	³⁾	2875	2875	3192	3192	3192	3034	3034	3192
Radius / radius	mm	127	127	141	141	141	134	134	141
 9 (97%)	sec	27							
 9	sec	30							
Temperatur / temperature	°C ¹⁾	- 6							
Probenerwärmung/Sample temp. rise	K ²⁾	10							

1324	1398							
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°								
	1482A	1483A	1484	1484				
			 ohne / without E2110-A					
			Falcon ®					
Kapazität / capacity	ml	15	15	50	50			
Maße / dimensions $\varnothing \times L$	mm	17 x 100	17 x 120	29 x 115	29 x 115			
Anzahl p. Rotor / number p. rotor		16	16	4	4			
Drehzahl / speed	RPM	4500	4500	4500	4500			
RZB / RCF	³⁾	3192	3305	3260	3260			
Radius / radius	mm	141	146	144	144			
 9 (97%)	sec	27						
 9	sec	30						
Temperatur / temperature	°C ¹⁾	- 6						
Probenerwärmung/Sample temp. rise	K ²⁾	10						

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.

1494			1427 + 1421														
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$																	
	1357		5229 ⁷⁾		5229		1326		5230		5230		5231		5231 ⁷⁾		
																	
	Rhe- sus																
							 						 				
Kapazität / capacity	ml	1	0,4	4,5 – 5		2,7 - 3		6	4	7		4,5 – 5		15		5 - 10	
Maße / dimensions $\varnothing \times L$	mm	6 x 45		11 x 92		11x 66		12×82	12×60	12 x 100		11 x 92		17 x 100		16 x 100	
Anzahl p. Rotor / number p. rotor		120		48		48		48		48		48		24		24	
Drehzahl / speed	RPM	5000		5000		5000		5000		5000		5000		5000		5000	
RZB / RCF ³⁾		4081		4053		4053		4053		3941		3941		3941		3941	
Radius / radius	mm	146		145		145		145		141		141		141		141	
 9 (97%)	sec	30															
 9	sec	32															
Temperatur / temperature	°C ¹⁾	- 7															
Probenerwärmung/Sample temp. rise	K ²⁾	+ 15															

1494	1427 + 1421								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°									
	1731	1732				1732		5237	5279
									
									
Kapazität / capacity	ml	25	5	1 – 5	2,6 – 2,9	4,9	4 – 7	9	4 – 5,5
Maße / dimensions ∅ x L	mm	25 x 90	12/13x75	13 x 75	13 x 65	13 x 90	13 x 100	14 x 100	15 x 75
Anzahl p. Rotor / number p. rotor		8	32	32	32	32	32	24	20
Drehzahl / speed	RPM	5000	5000	5000	5000	5000	5000	5000	5000
RZB / RCF	³⁾	3969	4025	4025	4025	4025	4025	3941	4109
Radius / radius	mm	142	144	144	144	144	144	141	147
 9 (97%)	sec	30							
 9	sec	32							
Temperatur / temperature	°C ¹⁾	- 7							
Probenerwärmung/Sample temp. rise	K ²⁾	15							

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








- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.
- 7) not possible to close the lid
- 6) Remove the inserts




1494	1425								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°									
	1431				1436	1437	1458	1435	1439
									
									
	Kapazität / capacity	ml	9	15	7,5 - 8,2	50	100	1,1 – 1,4	25
Maße / dimensions $\varnothing \times L$	mm	14 x 100	17 x 100	15 x 92	34 x 100	44 x 100	8 x 66	24 x 100	16 x 92
Anzahl p. Rotor / number p. rotor		28	28	28	4	4	36	8	16
Drehzahl / speed	RPM	5000	5000	5000	5000	5000	5000	5000	5000
RZB / RCF	³⁾	3913	3913	3913	3913	3801	4025	3913	3913
Radius / radius	mm	140	140	140	140	136	144	140	140
 9 (97%)	sec	30							
 9	sec	32							
Temperatur / temperature	°C ¹⁾	- 10							
Probenerwärmung/Sample temp. rise	K ²⁾	12							

1494	1425								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°									
	1438								
									
									
Kapazität / capacity	ml	2,6 – 2,9	2,7 – 3	4,9	4,5 – 5	1 – 5	4 – 7	5	4
Maße / dimensions $\varnothing$ x L	mm	13 x 65	11 x 66	13 x 90	11 x 92	13 x 75	13 x 100	12/13x75	12 x 60
Anzahl p. Rotor / number p. rotor		28	28	28	28	28	28	28	28
Drehzahl / speed	RPM	5000							
RZB / RCF	³⁾	3913							
Radius / radius	mm	140							
 9 (97%)	sec	30							
 9	sec	32							
Temperatur / temperature	°C ¹⁾	- 10							
Probenerwärmung/Sample temp. rise	K ²⁾	12							

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















- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.
- 4) can not be centrifugated when plug is attached, Scal. 10µl-300µl, 15ml, 30ml
- 5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".


1554	1560 + 1561 / 1565																
 $\angle 90^\circ$ max. Laufzyklen / max. cycles 100 000	 +  / 																
	mit Bioabdichtung / with bio-containment ⁵⁾					-----											
	max. Laufzyklen / max. cycles: 80000																
	max. Beladung / max. load: 290 g																
	1571		1593		1589			1588									
																	
																	
																	
Kapazität / capacity	ml	1,5		2,0		5		5		6		7		9		15	
Maße / dimensions $\varnothing \times L$	mm	11 x 38				17 x 59		12 x 75		12 x 82		12 x 100		14 x 100		17 x 100	
Anzahl p. Rotor / number p. rotor		56				16		28				20					
Drehzahl / speed	RPM	4500															
RZB / RCF	³⁾	3328/2332				3328		3215									
Radius / radius	mm	147/103				147		142									
 9 (97%)	sec	28															
 9	sec	31															
Temperatur / temperature	$^\circ\text{C}$ ¹⁾	- 8															
Probenerwärmung/Sample temp. rise	K ²⁾	10															

1554	1560 + 1561 / 1565								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90° max. Laufzyklen / max. cycles 100 000									
	mit Bioabdichtung / with bio-containment 5)				-----				
	max. Laufzyklen / max. cycles: 80000								
	max. Beladung / max. load: 290 g								
	1572	1573	1574	1575	1576	1594		---	
									
									
Kapazität / capacity	ml	15	25	50	85	100	125	150	200
Maße / dimensions ∅ x L	mm	18 x 102	24 x 100	34 x 100	38 x 101	44 x 100	51 x 100	51 x 116	56 x 112
Anzahl p. Rotor / number p. rotor		16	4						
Drehzahl / speed	RPM	4500							
RZB / RCF	3)	3260	3056	3124	3260	3124	3328		
Radius / radius	mm	144	135	138	144	138	147		
 9 (97%)	sec	28							
 9	sec	31							
Temperatur / temperature	°C 1)	- 8							
Probenerwärmung/Sample temp. rise	K 2)	10							

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










- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.
- 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 7) not possible to close the lid



1554	1560 + 1561 / 1565								
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$ max. Laufzyklen / max. cycles 100 000	 +  /  mit Bioabdichtung / with bio-containment 5) ----- max. Laufzyklen / max. cycles: 80000 max. Beladung / max. load: 290 g								
	1578	1579	1581	1582		1583	1584	1585	
									
									
	Kapazität / capacity	ml	30	50	12	25	30	50	10
Maße / dimensions $\varnothing$ x L	mm	25 x 110	30 x 115	17 x 100	25 x 90	25 x 110	29 x 115	16 x 80	26 x 95
Anzahl p. Rotor / number p. rotor		4		16		4		20	
Drehzahl / speed	RPM	4500							
RZB / RCF	³⁾	3328		3260	3328		3215		3260
Radius / radius	mm	147		144	147		142		144
 9 (97%)	sec	28							
 9	sec	31							
Temperatur / temperature	$^{\circ}\text{C}$ ¹⁾	- 8							
Probenerwärmung/Sample temp. rise	K ²⁾	10							

1554	1560 + 1561 / 1565					
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90° max. Laufzyklen / max. cycles 100 000	 +  /  mit Bioabdichtung / with bio-containment 5) max. Laufzyklen / max. cycles: 80000 max. Beladung / max. load: 290 g					
	1586	1575	1587			
						
			0534 4)			
						
Kapazität / capacity	ml	50	85	94	30	
Maße / dimensions Ø x L	mm	29 x 107	38 x 106	38 x 110 16)	44 x 105	
Anzahl p. Rotor / number p. rotor		4				
Drehzahl / speed	RPM	4500				
RZB / RCF	3)	3260		3192		
Radius / radius	mm	144		141		
 9 (97%)	sec					28
 9	sec					31
Temperatur / temperature	°C 1)					- 8
Probenerwärmung/Sample temp. rise	K 2)					10

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


- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.
- 4) can not be centrifugated when plug is attached, Skal. 10µl-300µl, 15ml, 30ml
- 5) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 16) Dimensions with lid

1554	1559														
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$ max. Laufzyklen / max. cycles 100 000	 max. Laufzyklen / max. cycles: 80000 max. Beladung / max. load: 200 g														
	1486			1482A			1482A + 4 x 0715			1486			1482A		
															
															
	Kapazität / capacity	ml	5	6	7	9	15	4 - 7	4 - 5,5	2,6 - 3,4	2,7 - 3	4,5-5	4,9	9-10	10
Maße / dimensions $\varnothing \times L$	mm	12 x 75	12 x 82	12 x 100	14 x 100	17 x 100	16 x 75	15 x 75	13 x 65	11 x 66	11 x 92	13 x 90	16 x 92	15 x 102	
Anzahl p. Rotor / number p. rotor		20			16		16		20		20		16		
Drehzahl / speed	RPM	4500													
RZB / RCF	³⁾	3215			3260		3215		3215			3260			
Radius / radius	mm	142			144		142		142			144			
 9 (97%)	sec	28													
 9	sec	31													
Temperatur / temperature	°C ¹⁾	- 8													
Probenerwärmung/Sample temp. rise	K ²⁾	9													

1554	1559																							
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$ max. Laufzyklen / max. cycles 100 000																								
	max. Laufzyklen / max. cycles: 80000																							
	max. Beladung / max. load: 200 g																							
	1486		1482A		1488		1487		1483A		1484		1482A		1484									
																								
																								
	Kapazität / capacity		ml		1 - 5		4 - 7		4 - 7		8,5-10		8		12		15		50		12		50	
	Maße / dimensions $\varnothing \times L$		mm		13 x 75		13 x 100		16 x 75		16x 100		16 x 125		17 x 102		17 x 120		29 x 115		17 x 100		29 x 115	
	Anzahl p. Rotor / number p. rotor				20				16				12				16		4		16		4	
Drehzahl / speed		RPM		4500																				
RZB / RCF		³⁾		3215		3260		3328		3147		3351		3305		3260		3305						
Radius / radius		mm		142		144		147		139		148		146		144		146						
 9 (97%)		sec		28																				
 9		sec		31																				
Temperatur / temperature		°C ¹⁾		- 8																				
Probenerwärmung/Sample temp. rise		K ²⁾		9																				

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




- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
3) Observe the tube manufacturer's instructions.

1554		1563					
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$ max. Laufzyklen / max. cycles 100 000		 max. Laufzyklen / max. cycles: 80000 max. Beladung / max. load: 160					
		1592	1592 + E2109	E2110-A			
							
							
Kapazität / capacity	ml	12	15	50	50		
Maße / dimensions $\varnothing \times L$	mm	17 x 100	17 x 120	29 x 115	29 x 115		
Anzahl p. Rotor / number p. rotor		8					
Drehzahl / speed	RPM	4500					
RZB / RCF	³⁾	3260					
Radius / radius	mm	144					
 9 (97%)	sec	28					
 9	sec	31					
Temperatur / temperature	°C ¹⁾	- 8					
Probenerwärmung/Sample temp. rise	K ²⁾	10					

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






- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.
- 5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 6) Remove the inserts

1624	1481 + 1492									
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°	 mit Bioabdichtung / with bio-containment ⁵⁾ 1383 									
										
	Kapazität / capacity	ml	6	7	4,9	4,5 - 5,0	2,7 - 3,0	2,6 - 2,9	1 - 5	5
	Maße / dimensions $\varnothing \times L$	mm	12 x 82	12 x 100	13 x 90	11 x 92	11 x 66	13 x 65	13 x 75	12 x 75
	Anzahl p. Rotor / number p. rotor		20	20	20	20	20	20	20	20
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000	
RZB / RCF	³⁾	2558	2558	2558	2558	2558	2558	2558	2558	
Radius / radius	mm	143	143	143	143	143	143	143	143	
 9 (97%)	sec	20								
 9	sec	25								
Temperatur / temperature	°C ¹⁾	- 15								
Probenerwärmung/Sample temp. rise	K ²⁾	7								

1624	1481 + 1492								
Ausschwingrotor 4-fach / Swing out rotor 4-times   90°	 mit Bioabdichtung / with bio-containment ⁵⁾								
	1383	1384	1396	1457	1343		1363	1365	
									
									
Kapazität / capacity	ml	4 – 7,0	50	85	1,1 – 1,4	3	4	25	30
Maße / dimensions $\varnothing \times L$	mm	13 x 100	29 x 115	38 x 106	8 x 66	10 x 60	10 x 88	25 x 90	25 x 110
Anzahl p. Rotor / number p. rotor		20	4	4	28	36	36	4	4
Drehzahl / speed	RPM	4000	4000	4000	4000	4000	4000	4000	4000
RZB / RCF	³⁾	2558	2665	2612	2576	2630	2630	2343	2665
Radius / radius	mm	143	149	146	144	147	147	131	149
 9 (97%)	sec	20							
 9	sec	25							
Temperatur / temperature	°C ¹⁾	- 15							
Probenerwärmung/Sample temp. rise	K ²⁾	7							

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









CP Kulturplatte /
Culture plate

DWP Deep Well Platte /
Deep well plate


MS Micronic System /
Micronic system




1418	1467								
Winkelrotor 8-fach / Angle rotor 8-times  ↙ 45°									
	1054-A + 0701	1054-A						0716	
									
									
Kapazität / capacity	ml	4	5	1,1 – 1,4	2,7 - 3	2,6 – 2,9	1 - 5	5	9-10
Maße / dimensions Ø x L	mm	12 x 60	12 x 75	8 x 66	11 x 66	13 x 65	13 x 75	13 x 75	16 x 92
Anzahl p. Rotor / number p. rotor		32	32	32	32	32	32	32	32
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	³⁾	2694	2762	2762	2762	2762	2762	2762	3215
Radius / radius	mm	119	122	122	122	122	122	122	142
 (97%)	sec	30							
 1	sec	31							
Temperatur / temperature	°C ¹⁾	- 5							
Probenerwärmung/Sample temp. rise	K ²⁾	20							

1418	1467					1468			
Winkelrotor 8-fach / Angle rotor 8-times  									
	0716					E2109			
									
									
Kapazität / capacity	ml	15	5 - 10	12	12	4 – 7	15	50	50
Maße / dimensions Ø x L	mm	17 x 100	13 x 100	17 x 102	17 x 100	13 x 100	17 x 120	29 x 107	29 x 115
Anzahl p. Rotor / number p. rotor		32	32	32	32	32	32	8	8
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	³⁾	3215	3215	3215	3305	3283	3147	3147	3147
Radius / radius	mm	142	142	142	146	145	139	139	139
 9 (97%)	sec	30							
 9	sec	31							
Temperatur / temperature	°C ¹⁾	- 5							
Probenerwärmung/Sample temp. rise	K ²⁾	20							

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





- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.
- 5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 8) recommended for high-speed centrifugation
- 9) Load only each second position of the rotor
- 14) With running time longer than 20 min., the sample temp. rise up will be > 20°K (only with cooling centrifuges)

1553								
		8) 2031	2023		2024		0788	
Winkelrotor 30-fach / Angle rotor 30-times  45° mit Bioabdichtung / with bio-containment ⁵⁾	---							
								
							9)	
Kapazität / capacity	ml	2,0	1,5	0,8	0,5	0,4	0,2	0,5
Maße / dimensions Ø x L	mm	11 x 38	11 x 38	8 x 45	8 x 30	6 x 45	6 x 18	10,7 x 46
Anzahl p. Rotor / number p. rotor		30	30	30	30	30	30	15
Drehzahl / speed	RPM	14150	14150	14150	14150	14150	14150	14150
RZB / RCF ³⁾		21713	21713	21713	21713	21713	21713	20818
Radius / radius	mm	97	97	97	97	97	97	93
 9 (97%)	sec	35						
 9	sec	32						
Temperatur / temperature	°C ¹⁾	- 1						
Probenerwärmung/Sample temp. rise	K ¹⁵⁾	20						

1555								
		2024	2023		8) 2031			
Ausschwingrotor 24-fach / Swing out rotor 24-times  90° mit Bioabdichtung / with bio-containment ⁵⁾	---					---		
								
Kapazität / capacity	ml	0,2	0,4	0,5	0,8	1,5	2,0	
Maße / dimensions Ø x L	mm	6 x 18	6 x 45	8 x 30	8 x 45	11 x 38	11 x 38	
Anzahl p. Rotor / number p. rotor		24						
Drehzahl / speed	RPM	13000						
RZB / RCF ³⁾		18327						
Radius / radius	mm	97						
 9 (97%)	sec	36						
 9	sec	31						
Temperatur / temperature	°C ¹⁾	3						
Probenerwärmung/Sample temp. rise	K ¹⁰⁾	20						

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)

- 3) Observe the tube manufacturer's instructions.



- 5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".





- 8) recommended for high-speed centrifugation

- 9) Load only each second position of the rotor

- 10) With running time longer than 10 min., the sample temp. rise up will be > 20°K (only with cooling centrifuges)



- 15) With running time longer than 15 min., the sample temp. rise up will be > 20°K (only with cooling centrifuges)


1556										
Winkelrotor 6-fach / Angle rotor 6-times  mit Bioabdichtung / with bio-containment ⁵⁾	1449		1477	1478						
										
										
Kapazität / capacity	ml	1,5	2,0	0,5	10	15	7,5 – 8,2	9 - 10	10	5 - 10
Maße / dimensions $\varnothing$ x L	mm	11 x 38	10,7x46	16 x 80	17 x 100	15 x 92	16 x 92	15 x 102	16 x 100	
Anzahl p. Rotor / number p. rotor		24	24	12	6					
Drehzahl / speed	RPM	9000								
RZB / RCF ³⁾		10324	10324	10414	10052					
Radius / radius	mm	114	114	115	111					
 9 (97%)	sec	36								
 9	sec	32								
Temperatur / temperature	°C ¹⁾	0								
Probenerwärmung/Sample temp. rise	K ²⁾	20								

1556										
Winkelrotor 6-fach / Angle rotor 6-times  35° mit Bioabdichtung / with bio-containment ⁵⁾			1466	1454	1447	1446			1463	
	---									
										
Kapazität / capacity	ml	85	94	15	50	30	50	85	50	50
Maße / dimensions Ø x L	mm	38 x 106	38 x 110 16)	17 x 120	29 x 115	26 x 95	29 x 107	38 x 101	35 x 105	34 x 100
Anzahl p. Rotor / number p. rotor		6								
Drehzahl / speed	RPM	9000								
RZB / RCF ³⁾		10595		10052		9690	10142	10595		10414
Radius / radius	mm	117		111		107	112	117		115
 9 (97%)	sec	36								
 9	sec	32								
Temperatur / temperature	°C ¹⁾	0								
Probenerwärmung/Sample temp. rise	K ²⁾	20								

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.
- 5) in conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".
- 16) Dimensions with lid

1613		1054-A						1054-A /0701		
Winkelrotor 12-fach / Angle rotor 12-times  ∠ 35°										
										
Kapazität / capacity	ml	5	1 – 5	6	2,6 – 2,9 13 x 65	2,7 - 3 11 x 66	1,1 – 1,4	4	8,5 - 10	8
Maße / dimensions ∅ x L	mm	12/13 x 75	13 x 75	12 x 82	13 x 65	11 x 66	8 x 66	12 x 60	16 x 100	16 x 125
Anzahl p. Rotor / number p. rotor		12	12	12	12		12	12	12	6
Drehzahl / speed	RPM	6000	6000	6000	6000		6000	6000	6000	6000
RZB / RCF	³⁾	3300	3300	3300	3300		3300	3260	4146	4146
Radius / radius	mm	82	82	82	82		82	81	103	103
 9 (97%)	sec	15								
 9	sec	15								
Temperatur / temperature	°C ¹⁾	- 16								
Probenerwärmung/Sample temp. rise	K ²⁾	4								

1613									
Winkelrotor 12-fach / Angle rotor 12-times  $\angle 35^\circ$									
									
Kapazität / capacity	ml	4,5 - 5	4,9	7,5 – 8,2	9 – 10	10			
Maße / dimensions $\varnothing \times L$	mm	11 x 92	13 x 90	15 x 92	16 x 92	15 x 102			
Anzahl p. Rotor / number p. rotor		12	12	12	12	12			
Drehzahl / speed	RPM	6000	6000	6000	6000	6000			
RZB / RCF	³⁾	4146	4146	4146	4146	4146			
Radius / radius	mm	103	103	103	103	103			
 9 (97%)	sec	15							
 9	sec	15							
Temperatur / temperature	°C ¹⁾	- 16							
Probenerwärmung/Sample temp. rise	K ²⁾	4							

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







10) With running time longer than 10 min., the sample temp. rise up will be > 20°K (only with cooling centrifuges)




1494	1452								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°									
	1662						1670		
							 11)		
	1663	1664	1665	1666	1667	1668	1663	1664	
									
Kapazität / capacity	ml	1	2	4	8	3 x 2	4 x 1	1	2
Maße / dimensions $\varnothing$ / A	mm ²	6,2 / 30	8,7 / 60	12,4 / 120	17,5 / 240	8,7 / 60	6,2 / 30	6,2 / 30	8,7 / 60
Anzahl p. Rotor / number p. rotor		4	4	4	4	4	4	4	4
Filterkarten / filter cards		1675	1675	1675	1676	1677	1678	1692	1692
Drehzahl / speed	RPM	5000							
RZB / RCF	³⁾	2879							
Radius / radius	mm	103							
 9 (97%)	sec	30							
 9	sec	32							
Temperatur / temperature	°C ¹⁾	- 10							
Probenerwärmung/Sample temp. rise	K ²⁾	12							

1494	1452							
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$								
	1670							
	 ¹¹⁾							
	1665	1666	1667	1668				
								
Kapazität / capacity	ml	4	8	3 x 2	4 x 1			
Maße / dimensions $\varnothing$ / A	mm ²	12,4 / 120	17,5 / 240	8,7 / 60	6,2 / 30			
Anzahl p. Rotor / number p. rotor		4	4	4	4			
Filterkarten / filter cards		1692	1691	1694	1693			
Drehzahl / speed	RPM	5000	5000	5000	5000			
RZB / RCF	³⁾	2879	2879	2879	2879			
Radius / radius	mm	103	103	103	103			
 9 (97%)	sec	30						
 9	sec	32						
Temperatur / temperature	°C ¹⁾	- 10						
Probenerwärmung/Sample temp. rise	K ²⁾	12						

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
- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
- 3) Observe the tube manufacturer's instructions.
- 11) Object slide will not stand RCF values exceeding 1100




1624	1660 + 1661								
Ausschwingrotor 4-fach / Swing out rotor 4-times  ∠ 90°									
	1662						1670		
							 11)		
	1663	1664	1665	1666	1667	1668	1663	1664	
									
Kapazität / capacity	ml	1	2	4	8	3 x 2	4 x 1	1	2
Maße / dimensions Ø / A	mm ²	6,2 / 30	8,7 / 60	12,4 / 120	17,5 / 240	8,7 / 60	6,2 / 30	6,2 / 30	8,7 / 60
Anzahl p. Rotor / number p. rotor		4	4	4	4	4	4	4	4
Filterkarten / filter cards		1675	1675	1675	1676	1677	1678	1692	1692
Drehzahl / speed	RPM	4000							
RZB / RCF	³⁾	1646							
Radius / radius	mm	92							
 9 (97%)	sec	20							
 9	sec	25							
Temperatur / temperature	°C ¹⁾	- 16							
Probenerwärmung/Sample temp. rise	K ²⁾	3							

1624	1660 + 1661					1680			
Ausschwingrotor 4-fach / Swing out rotor 4-times  $\angle 90^\circ$									
	1670					1662			
	 11)								
	1665	1666	1667	1668		1671	1672	1673	
									
Kapazität / capacity	ml	4	8	3 x 2	4 x 1		[1] 0,5	[1] 0,5	[1] 0,5
Maße / dimensions $\varnothing$ / A	mm ²	12,4 / 120	17,5 / 240	8,7 / 60	6,2 / 30		6,2 / 30	8,7 / 60	12,4 / 120
Anzahl p. Rotor / number p. rotor		4	4	4	4		4	4	4
Filterkarten / filter cards		1692	1691	1694	1693		[1] 1696	[1] 1697	[1] 1698
Drehzahl / speed	RPM	4000	4000	4000	4000		4000	4000	4000
RZB / RCF	³⁾	1646	1646	1646	1646		1467	1467	1467
Radius / radius	mm	92	92	92	92		82	82	82
 9 (97%)	sec	20							
 9	sec	25							
 0	sec	390							
Temperatur / temperature	°C ¹⁾	- 16							
Probenerwärmung/Sample temp. rise	K ²⁾	3							

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
2) Sample temp. rise during maximum speed and 1 hour running time (only with centrifuges without cooling)
3) Observe the tube manufacturer's instructions.
11) Object slide will not stand RCF values exceeding 1100
[1] One-step method

1626	1660 + 1661								
Ausschwingrotor 6-fach / Swing out rotor 6-times  $\angle 90^\circ$	 + 								
	1662						1670		
							 11)		
	1663	1664	1665	1666	1667	1668	1663	1664	
									
Kapazität / capacity	ml	1	2	4	8	3 x 2	4 x 1	1	2
Maße / dimensions $\varnothing$ / A	mm ²	6,2 / 30	8,7 / 60	12,4 / 120	17,5 / 240	8,7 / 60	6,2 / 30	6,2 / 30	8,7 / 60
Anzahl p. Rotor / number p. rotor		6	6	6	6	6	6	6	6
Filterkarten / filter cards		1675	1675	1675	1676	1677	1678	1692	1692
Drehzahl / speed	RPM	4000							
RZB / RCF	3)	2039							
Radius / radius	mm	114							
 9 (97%)	sec	20							
 9	sec	22							
Temperatur / temperature	°C 1)	- 16							
Probenerwärmung/Sample temp. rise	K 2)	7							

1626	1660 + 1661				1680			
Ausschwingrotor 6-fach / Swing out rotor 6-times  $\angle 90^\circ$								
	1670				1662			
	 ¹¹⁾							
	1665	1666	1667	1668		1671	1672	1673
								
Kapazität / capacity	ml	4	8	3 x 2	4 x 1		[1] 0,5	[1] 0,5
Maße / dimensions $\varnothing$ / A	mm ²	12,4 / 120	17,5 / 240	8,7 / 60	6,2 / 30		6,2 / 30	8,7 / 60
Anzahl p. Rotor / number p. rotor		6	6	6	6		6	6
Filterkarten / filter cards		1692	1691	1694	1693		[1] 1696	[1] 1697
Drehzahl / speed	RPM	4000	4000	4000	4000		4000	4000
RZB / RCF	³⁾	2039	2039	2039	2039		1842	1842
Radius / radius	mm	114	114	114	114		103	103
 9 (97%)	sec	20						
 9	sec	22						
 0	sec	330						
Temperatur / temperature	°C ¹⁾	- 16						
Probenerwärmung/Sample temp. rise	K ²⁾	7						

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