

HETTWIPES

CENTRIFUGE AND EQUIPMENT WIPES

Keep It Clean

- Safely clean and disinfect all your Hettich products
- Will not dry, crack or damage your equipment
- No alcohol, citrus, ammonia, bleach, phenol
- No harsh chemicals
- Part of Hettich's preventative maintenance schedule for cleaning and disinfecting your equipment.

INGREDIENTS:

N-Alkyl (60% C14, 30% C15, 5% C12, 5% C18)	
Dimethyl benzyl ammonium chloride	0.105%
N-Alkyl (68% C12, 32% C14)	
Dimethyl ethylenbenzyl ammonium chloride	0.105%

ORDERING INFORMATION:

Hettich HettWipes (Case of 6)	HW185-6
-------------------------------	----------------

Approved for use with Hettich equipment only, not approved for use on other lab equipment. Not for infectious diseases. Please check the MSDS for chemical compatibility (available at: hettweb.com/hettwipes).



LAB ASSISTANT

Give your HettWipes personality. Comes fully outfitted with lab coat and eye stickers.



hettweb.com/hettwipes

CONTACT YOUR LOCAL HETTICH REPRESENTATIVE TO LEARN MORE: