

The Hettich logo is displayed in a white, italicized, sans-serif font against a black rectangular background. A thin red line is visible at the bottom edge of the black background.

Conquer multi-cellular samples
one layer at a time



The Hettich Cyto-System:
fast single layer slide preparation

OVER 100 YEARS OF QUALITY, SAFETY & RELIABILITY
For additional information on all of our products please visit www.hettweb.com

Hettich Cyto-System

Slide preparation is the optimal way of viewing cells in a single layer. The **Hettich Cyto-System** offers a dependable method to create a monolayer of cells in one quick easy step. From 0.2 mL to 8 mL, single funnel or up to 24 sample locations per spin. Optimize your lab with a quick rotor exchange to spin round bottom, conical/ skirted or other flat bottom tubes up to 100 mL, all on **one** centrifuge!

CYTO-SYSTEM OPTIONAL FEATURES

Cyto-rotor	Fits standard centrifuges
	4, 8 and 12 place
	Open, closed or quick-release versions
Slides	Coated and uncoated, etched location rings
Funnels	Single, double & large funnels
Filter Cards	Attached & unattached filter cards
Eco-friendly	EXCLUSIVE - Reusable funnel system
Market Comparison	Accessories fit Shandon Cytospin®
Sedimentation	Wet or dry samples
Max Capacity	12 slides (up to 2 sample locations per slide)



HETTICH CLIP SYSTEM

Flexibility for low and high throughput labs. Several rotor options with 4, 8 and 12 places. Closed rotors for improved containment and quick release rotors for easy transport.

Ideal for: Urinalysis, Abdominal & Spinal fluids, Cysts, Biopsies, Ascites



CLIP SYSTEM ACCESSORIES

Single and double funnels with attached or unattached filter cards plus sealing cap. Corresponding glass slides with etched location ring.

Ideal for: Pap tests, Urinalysis, Pleura-pericardium fluid, Fine Needle Aspiration



HETTICH BIO-SEAL SYSTEM

Certified bio-containment carrier with clear sealing lid. Optimize lab throughput with multiple funnel options in 2, 3, or 4 funnels per slide. Perfect funnel for every application from 0.2ml - 8.0ml.

Ideal for: CSF, Viral infections, Bronchial Secretions

Cyto-centrifuge with quick rotor exchange

Applications

CSF, Pap tests, Urinalysis, Cysts, Cultures, Ascities, Bronchial secretions, Viral Infections, Abdominal Fluids, ISH, FISH, Pleura-pericardium fluid

Laboratories

Cytology, Hematology, Oncology, Blood Labs, Pediatrics, Microbiology, Pathology, Gynecology, Neurology, Virology, Internal Medicine, Pediatrics, Immunohistochemistry

Shandon Cytospin® is a registered trademark of Thermo Fisher Scientific Inc



Hettich Centrifuges distributed by:

Copyright © 2008 Hettich Instruments, LP. All rights reserved. For the most up to date information on Hettich please visit www.hettweb.com.

♻️ Printed in the USA using soy based inks on 100% post consumer waste recycled paper. Paper is produced with 100% renewable electric energy & manufactured carbon neutral.

PS-CYTOLOGY-E RevA 09/08

For ordering information please visit:
www.hettweb.com/cytology
or contact your local representative