


TOP
S
P
O



HANDZENTRIFUGE
Manual Centrifuge

HANDZENTRIFUGE



Maintenance-free, bench-mounted manually operated centrifuge

This convenient centrifuge can be easily clamped to a bench or table, and operated without electricity. Its self-lubricating gear makes it virtually maintenance-free.

It can be operated with two different rotors.

With rotor **1014** conical 15 ml tubes (Cat. No. 0515, non-graduated or No. 0516, graduated) can be accelerated to 1,300 RCF.

Rotor **1025** also accommodates a variety of other commercial tubes (also with round bottom) up to 15 ml. The tubes may have a max. length of 125 mm and a max. diameter of 17 mm.



illustrated with rotor 1014
↘ 90°

illustrated with rotor 1025
↘ 50°

ROTORS	4-place	4-place
Max. capacity	4 x 15 ml ^{*)}	4 x 15 ml
Max. RPM (speed)	3,000 min ⁻¹	3,000 min ⁻¹
Max. RCF	1,300	1,076
Cat. No.	1014	1025

^{*)} conical only

HANDZENTRIFUGE

Manual centrifuge, without rotor	
Dimensions (H x W x D)	285 x 140 x 175 mm
Weight	approx. 0.9 kg
Cat. No.	1011