

**FREE &
EASY**



HAND CENTRIFUGE
Manual Centrifuge

Hettich
ZENTRIFUGEN

HAND CENTRIFUGE

Maintenance-free, bench-mounted, manual centrifuge

This convenient, manually operated centrifuge can be easily clamped to a bench or table - independently of an external power supply. Its self-lubricating gear makes it virtually maintenance-free.

It can be operated with two different rotors.

With rotor No. **1014** conical 15 ml tubes can be accelerated to 1,300 RCF.

Rotor No. **1025** also accommodates a variety of other commercial tubes (also with round bottom) up to 15 ml. Max. length of the tubes: 125 mm, max. diameter: 17 mm. Max. RCF: 1,076.



illustrated with rotor No. 1014

illustrated with rotor No. 1025

ROTORS

	4-place	4-place
Max. capacity	4 x 15 ml (conical)	4 x 15 ml
Max. RPM (speed)	3,000 min ⁻¹	3,000 min ⁻¹
Max. RCF	1,300	1,076
Cat. No.	1014	1025

HAND CENTRIFUGE

Manual centrifuge, without rotor	
Dimensions (H x W x D) in mm	285 x 140 x 175
Weight	approx. 0.9 kg
Cat. No.	1011