

HANDZENTRIFUGE

MANUAL CENTRIFUGE FOR FIELD APPLICATIONS

Manual powered for field applications with limited access to electricity. Lightweight, portable and manually powered this centrifuge is ideal for work in remote areas. Simply secure the centrifuge to a benchtop, desk, tailgate or work surface and turn by hand, no electric or battery power source required.

Performance and versatility when you need it most.

Your choice of a 4-place (15 mL) round bottom tube or conical tube rotor is available. The manual powered self-lubricating gear makes this centrifuge virtually maintenance free while generating a max RPM of 3,000 and a max RCF (G-FORCE) of 1,298.



MODEL	CONICAL ROTOR	ROUND-BOTTOM ROTOR
Max. Capacity	4 x 15 ml	4 x 15 ml
Rotor Angle	90°	50°
Max. RPM / RCF	3,000 / 1,298	3,000 / 1,077
Max. Blood Tubes*		4 (up to 10 ml)
Max. Conical Tubes*	4 (15 ml)	
Temperature Control	Ambient	Ambient
Dims. (H x W x D)	11 x 5.5 x 7 in.	11 x 5.5 x 7 in.

*Please inquire about the fit of your specific tube dimensions.

Learn more at: hettweb.com/manual

Hettich | North America

100 Cummings Center, Ste. 136L Beverly, MA 01915, U.S.A. Toll free: 1 (866) 370-4388 info@hettweb.com hettweb.com © 2019 Hettich Instruments, LP. All rights reserved. Lit. No. PS-MANUAL-EN-01-2019.

Hettich products are designed to help you achieve optimal results for your application and are built to perform to the specifications outlined in the operators manual. For application specific information and settings, please refer to your organizations standard operating procedure. As always, our Hettich representatives are here to help determine which Hettich products and/or accessories best fit your laboratory requirements.

Small Centrifuge

hettweb.com