

# MIKRO 185

MICROLITER / SPIN COLUMNS / PEDIATRIC

**Basic and economical micro centrifuge for your demanding requirements.** The MIKRO 185 delivers a variety of rotor options with or without bio-containment and is backed by a 5-year warranty. Users have the option of 12-place, 18-place, and 24-place angle rotors that fit both pediatric tubes and standard microliter tubes. An 18-place rotor for spin column kits is designed to protect the tube cap from breaking during centrifugation.

**Higher RCF values saves you time.** With its 24-place rotor, the MIKRO 185 is capable of reaching an RCF of 17,008 x g, reducing spin time without compromising sample quality.

**Flexible control panel allowing you to adjust your application parameters.** Adjust time, speed (RPM), RCF and braking with the easy to use LCD display. Pulse key available allowing you to hold and press start key for short centrifugation runs.

**Visual and audible end-of-cycle notification for a productive work-flow.** Customizable and independent control of visual and/or audio alarms can be set for end of cycle notification. This provides an across-the-room indication that your unit is ready for the next operation.



MODEL	MIKRO 185
Max. Capacity	24 x 1.5 - 2.0 ml
Max. RPM / RCF	14,000 / 18,845
Max. Spin Column Kits*	18
Max. Pediatric Tubes*	12
Temperature Control	Ambient
Dims (H x W x D)	9 x 10 x 14 in.
Warranty†	5 yr.

\*Please inquire about the fit of your specific tube dimensions. † Warranty valid in U.S. and Canada only, see warranty statement for complete details.

Learn more: [hettweb.com/products/mikro185/](http://hettweb.com/products/mikro185/)

## MIKRO 185 CENTRIFUGE PACKAGES

### MIKRO 185 MICROLITER TUBE PACKAGE - 12

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	1.5 - 2.0	11 x 38	≤ 14,000	≤ 15,558	<b>185MICRO12</b>
	<b>Includes:</b> (1) Centrifuge, (1) Fixed-angle Rotor, (12) Adapters 1.5 mL					

### MIKRO 185 MICROLITER TUBE PACKAGE - 24 - BIO-CONTAINMENT

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	24	1.5 - 2.0	11 x 38	≤ 14,000	≤ 18,845	<b>185MICRO24</b>
	<b>Includes:</b> (1) Centrifuge, (1) Fixed-angle Rotor, (24) Adapters 1.5 mL, (1) Bio-containment Lid					

### MIKRO 185 PEDIATRIC TUBE PACKAGE - 12

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	0.5	10.7 x 36	≤ 14,000	≤ 14,462	<b>185MICRO12-PED</b>
	12	0.8	8 x 45	≤ 14,000	≤ 15,558	
	<b>Includes:</b> (1) Centrifuge, (1) Fixed-angle Rotor, (12) Adapters 0.5 mL, (12) Adapters 0.8 mL					

### MIKRO 185 SPIN COLUMN PACKAGE - 18

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	18	1.5 - 2.0	11 x 38	≤ 14,000	≤ 16,654	<b>185MICROSC</b>
	<b>Includes:</b> (1) Centrifuge, (1) Fixed-angle Rotor, (18) Adapters 1.5 mL, (1) Bio-containment Lid (optional)					
						For bio-containment add <b>"-BC"</b> to package no. above

### MIKRO 185 CUSTOM PACKAGE

	<b>Supports over 29 application configurations, please inquire for more information</b>					
	<b>Tube types include:</b> Microliter (0.2 - 0.8 mL adapter combinations are available, inquire for details), Micro Spin Columns, Pediatric					

Learn more: [hettweb.com/products/mikro185/](http://hettweb.com/products/mikro185/)

**Hettich | North America**  
 100 Cummings Center, Ste. 136L  
 Beverly, MA 01915, U.S.A.  
 Toll free: 1 (866) 370-4388  
[info@hettweb.com](mailto:info@hettweb.com)  
[hettweb.com](http://hettweb.com)

© 2024 Hettich Instruments, LP. All rights reserved.  
 Lit. No. PS-MIKRO-185-EN-01-2024.

Hettich products are designed to help you achieve optimal results for your application and are built to perform to the specifications outlined in the operators manual. For application specific information and settings, please refer to your organizations standard operating procedure. As always, our Hettich representatives are here to help determine which Hettich products and/or accessories best fit your laboratory requirements.