

MIKRO 185

MICROLITER / SPIN COLUMNS / PEDIATRIC

Basic and economical micro centrifuge for your demanding requirements. The MIKRO 185 delivers a variety of rotor options with or without bio-containment and is backed by a 5-year warranty. Users have the option of 12-place, 18-place, and 24-place angle rotors that fit both pediatric tubes and standard microliter tubes. An 18-place rotor for spin column kits is designed to protect the tube cap from breaking during centrifugation.

Higher RCF values saves you time. With its 24-place rotor, the MIKRO 185 is capable of reaching an RCF of 17,008 x g, reducing spin time without compromising sample quality.

Flexible control panel allowing you to adjust your application parameters. Adjust time, speed (RPM), RCF and braking with the easy to use LCD display. Pulse key available allowing you to hold and press start key for short centrifugation runs.

Visual and audible end-of-cycle notification for a productive work-flow. Customizable and independent control of visual and/or audio alarms can be set for end of cycle notification. This provides an across-the-room indication that your unit is ready for the next operation.



MODEL	MIKRO 185
Max. Capacity	24 x 1.5 - 2.0 ml
Max. RPM / RCF	14,000 / 18,845
Max. Spin Column Kits*	18
Max. Pediatric Tubes*	12
Temperature Control	Ambient
Dims (H x W x D)	9 x 10 x 14 in.
Warranty†	5 yr.

*Please inquire about the fit of your specific tube dimensions. † Warranty valid in U.S. and Canada only, see warranty statement for complete details.

Learn more: hettweb.com/mikro185

MIKRO 185 CENTRIFUGE PACKAGES

MIKRO 185 MICROLITER TUBE PACKAGE - 12

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	1.5 - 2.0	11 x 38	≤ 13,300	≤ 14,041	185MICRO12
	Includes: (1) Centrifuge, (1) Fixed-angle Rotor, (12) Adapters 1.5 mL					

MIKRO 185 MICROLITER TUBE PACKAGE - 24 - BIO-CONTAINMENT

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	24	1.5 - 2.0	11 x 38	≤ 13,300	≤ 17,008	185MICRO24
	Includes: (1) Centrifuge, (1) Fixed-angle Rotor, (24) Adapters 1.5 mL, (1) Bio-containment Lid					

MIKRO 185 PEDIATRIC TUBE PACKAGE - 12

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	0.5	10.7 x 36	≤ 13,300	≤ 16,019	185MICRO12-PED
	12	0.8	8 x 45	≤ 13,300	≤ 17,008	
	Includes: (1) Centrifuge, (1) Fixed-angle Rotor, (12) Adapters 0.5 mL, (12) Adapters 0.8 mL					

MIKRO 185 SPIN COLUMN PACKAGE - 18

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	18	1.5 - 2.0	11 x 38	≤ 13,300	≤ 15,030	185MICROSC
	Includes: (1) Centrifuge, (1) Fixed-angle Rotor, (18) Adapters 1.5 mL, (1) Bio-containment Lid (optional)					

For bio-containment add **"-BC"** to package no. above

MIKRO 185 CUSTOM PACKAGE

	Supports over 29 application configurations, please inquire for more information					
	Tube types include: Microliter (0.2 - 0.8 mL adapter combinations are available, inquire for details), Micro Spin Columns, Pediatric					

Learn more: hettweb.com/mikro185

Hettich | North America
 100 Cummings Center, Ste. 136L
 Beverly, MA 01915, U.S.A.
 Toll free: 1 (866) 370-4388
info@hettweb.com
hettweb.com

© 2019 Hettich Instruments, LP. All rights reserved.
 Lit. No. PS-MIKRO-185-EN-01-2019.

Hettich products are designed to help you achieve optimal results for your application and are built to perform to the specifications outlined in the operators manual. For application specific information and settings, please refer to your organizations standard operating procedure. As always, our Hettich representatives are here to help determine which Hettich products and/or accessories best fit your laboratory requirements.