

MICRO CENTRIFUGES

MIKRO 185 / MIKRO 200 / MIKRO 220



Micro centrifuge models offer variety and high performance in a quality compact package. Our compact micro centrifuges are available in three

varieties and are coupled with multiple rotor choices to accommodate a diverse work-flow. Accessories are available for all micro tubes, spin column kits, conical tubes up to 50 ml, cryo tubes, pcr strips and all standard blood tubes. Autoclavable, phenol resistant, bio-containment lids are also available.

Speed and capacity across the entire MIKRO series. Whether you need a basic 24-place micro centrifuge (MIKRO 185) or a more advanced unit with greater capacity, rotor variety and optional refrigeration (MIKRO 200 / 220), our units are built for high-speed centrifugation for efficient and optimal performance. **Robust and steady temperature control systems in a compact footprint.** Our MIKRO 200 R (-10 to +40°C) and MIKRO 220 R (-20 to +40°C) centrifuges come with a quiet and reliable refrigeration system positioned at the rear of the unit. This maintains the compact width of the centrifuge and occupies minimal bench space.

Built in Germany for lasting performance. All Hettich centrifuges are made from high quality, durable materials for a long in-field service life. Each unit is backed by a limited 5-year warranty and offer robust temperature control systems, all metal containment, maintenance free motors, metal lid locks and easy to use control panels.

Reliable, quiet and safe. You can rely on our centrifuges to deliver optimal performance day-in and day-out in a safe, distraction free environment.

Powerful results at the micro level

hettweb.com



HETTICH MICRO CENTRIFUGES



MODEL	MIKR0 185	MIKRO 200 R	MIKRO 220 R
Ideal for	Run up to 24 microliter tubes, 12 pediatric tubes, or 18 spin column kits in a compact frame. High-performance with a max RCF of up to 18,845.	Up to 30 micro tubes with high RCF for best sample results. Specialty rotor designed for spin-column kits. 4-place PCR strip rotor available. Optional refrigerated model with fast pre-cool function for quick operation. Spin lipemic samples for a fraction of the cost.	Up to 48 micro tubes with bio-con- tainment. Wide rotor variety for tubes up to 50 ml including blood tubes. Optional refrigerated model available. Spin lipemic samples for a fraction of the cost.
Micro Capacity (1.5 – 2.0 ml)	24	30	48 (60)
PCR Strips		4 x 8	6 x 8
Spin Column Kits	18	24	
Cryo Tubes (1.8 ml)	12	20	20
Pediatric Tubes (0.5 ml)	12	15	15
Conical Tubes (15 / 50 ml)			6/3
Blood Tubes (Up to 10 ml)			12
Max. RPM / RCF	14,000 / 18,845	15,000 / 21,382	18,000 / 31,514
Temperature Control	Ambient	-10 °C to +40 °C	-20 °C to +40 °C
Dims (H x W x D)	9 x 10 x 14 in.	10 x 11 x 14 10 x 11 x 22 in.	12 x 13 x 16.5 12 x 13 x 25.5 in.
Warranty [†]	5 yr.	5 yr.	5 yr.

*Please inquire about the fit of your specific tube dimensions. † Warranty valid in U.S. and Canada only, see warranty statement for complete details

Please inquire about our **micro centrifuge packages**, each package comes complete with everything you need to perform your centrifugation.

Learn more: hettweb.com/microcentrifuges

Hettich | North America

100 Cummings Center, Ste. 136L Beverly, MA 01915, U.S.A. Toll free: 1 (866) 370-4388 info@hettweb.com hettweb.com © 2019 Hettich Instruments, LP. All rights reserved. Lit. No. PS-MIKRO-CENTRIFUGES-EN-01-2019.

Hettich products are designed to help you achieve optimal results for your application and are built to perform to the specifications outlined in the operators manual. For application specific information and settings, please refer to your organizations standard operating procedure. As always, our Hettich representatives are here to help determine which Hettich products and/or accessories best fit your laboratory requirements.