

ROTOFIX 32 A

BLOOD / CELL CULTURE / URINE / CYTOLOGY

Economical, rugged and simple design for uncomplicated performance. The ROTOFIX 32 A offers labs an indispensable, non-refrigerated, general-purpose benchtop centrifuge at a great value. The multi-functional and user-friendly centrifuge over-performs on everyday tasks and delivers consistent results. Ideal for all standard tubes up to 100 ml and specialty tubes such as cytology clips/funnels and chrome bath tubes.

Eliminate the need to purchase an expensive dedicated cytocentrifuge. The ROTOFIX 32 A performs multiple applications including cytology, in the same centrifuge, saving you money and bench space. Hettich's 12-place Clip Rotor is capable of spinning 0.2 to 6 ml single, double or large funnel formats with optional bio-containment, eliminating the need to purchase an expensive dedicated cytocentrifuge. Learn more about the Hettich Cyto-System for complete cytology solutions.



MODEL	ROTOFIX 32 A
Max. Capacity	4 x 100 ml, 6 x 94 ml
Max. RPM / RCF	6,000 / 4,226
Max. Blood Tubes*	40 (13 x 75 ml) [†] , 16 (16 x 100 ml)
Max. Conical Tubes*	16/4 (15/50 ml)
Temperature Control	Ambient
Dims (H x W x D)	10 x 14 x 17 in.
Warranty [†]	5 yr.

*Please inquire about the fit of your specific tube dimensions. † Warranty valid in U.S. and Canada only, see warranty statement for complete details.

Learn more: hettweb.com/rotofix32a

ROTOFIX 32 A CENTRIFUGE PACKAGES

ROTOFIX 32 A BLOOD TUBE PACKAGE

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	20	1.6 - 7	13 x 75/100	≤ 4,000	≤ 2,558	32ABLOOD
	16	4 - 10	16 x 75/100	≤ 4,000	≤ 2,522	
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 13 x 75/100 mm, (4) Adapters 16 x 75/100 mm, (4) Bio-containment Lids (optional)						For bio-containment add "-BC" to package no. above

ROTOFIX 32 A URINALYSIS PACKAGE

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	12 KOVA	17 x 102	≤ 4,000	≤ 2,451	32AURINALYSIS
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters KOVA					

ROTOFIX 32 A CELL CULTURE PACKAGE - 4

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	16	15 Conical	17 x 120	≤ 4,000	≤ 2,612	32ACELLCULTURE4
	4	50 Conical	29 x 115	≤ 4,000	≤ 2,576	
Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 15 mL, (4) Adapters 50 mL						

ROTOFIX 32 A CELL CULTURE PACKAGE - 6

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	6	15 Conical	17 x 120	≤ 4,000	≤ 2,701	32ACELLCULTURE6
	6	50 Conical	29 x 115	≤ 4,000	≤ 2,701	
Includes: (1) Centrifuge, (1) Horizontal Rotor, (6) Buckets, (6) Adapters 15 mL						

ROTOFIX 32 A CYTOLOGY PACKAGE - 12

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	0.4 - 6	--	≤ 2,000	≤ 438	32ACYTOCLIP12
	Includes: (1) Centrifuge, (1) Basic Cytology Rotor, (2) Cyto Clip 6-packs, (1) Cell Funnel, disposable 50-pack					

ROTOFIX 32 A CYTOLOGY PACKAGE - 12-S

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.
	12	0.4 - 6	--	≤ 2,000	≤ 438	32ACYTOCLIP12-S
	Includes: (1) Centrifuge, (1) Shielded Cytology Rotor, (2) Cyto Clip 6-packs, (1) Cell Funnel, disposable 50-pack					

ROTOFIX 32 A CUSTOM PACKAGE

	Supports over 210 application configurations, please inquire for more information					
	Tube types include: Microliter up to 5 mL, Blood, Pediatric, Urine, Conical and Skirted up to 50 mL, Standard up to 100 mL, Cytology					

Learn more: hettweb.com/rotofix32a

Hettich | North America
 100 Cummings Center, Ste. 136L
 Beverly, MA 01915, U.S.A.
 Toll free: 1 (866) 370-4388
info@hettweb.com
hettweb.com

© 2019 Hettich Instruments, LP. All rights reserved.
 Lit. No. PS-ROTOFIX-32-A-EN-01-2019.

Hettich products are designed to help you achieve optimal results for your application and are built to perform to the specifications outlined in the operators manual. For application specific information and settings, please refer to your organizations standard operating procedure. As always, our Hettich representatives are here to help determine which Hettich products and/or accessories best fit your laboratory requirements.