

ROTOFIX 46 | H

ASTM / SCHLENK / CHROME BATH / BABCOCK

Specialty industrial centrifuge for unique applications. The ROTOFIX 46, available with optional heating from +10 to +90 °C, offers solutions for organic and inorganic chemistry, preparatory and analytical chemistry, electroplating, and petro-chemistry. A universal swing-out rotor also accommodates standard tubes & bottles ranging from 3 mL to 290 ml.



MODEL	ROTOFIX 46	ROTOFIX 46 H (Heated)
Max. Capacity	4 x 290 ml	4 x 100 ml
Max. RPM / RCF	4,000 / 3,095	2,000 / 984
Max. Blood Tubes*	4 (up to 100 ml)	4 (up to 100 ml)
Max. Conical Tubes*	4 (165 mm)	4 (165 mm)
Temperature Control	Ambient	+10 °C to +90 °C
Dims (H x W x D)	14 x 21 x 25 in.	14 x 21 x 25 in.
Warranty†	5 yr.	5 yr.

*Please inquire about the fit of your specific tube dimensions. † Warranty valid in U.S. and Canada only, see warranty statement for complete details.

Learn more: hettweb.com/rotofix46

ROTOFIX 46 | H CENTRIFUGE PACKAGES

ROTOFIX 46 **ASTM** PACKAGE - 50

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (HEATED)
	4	50	45 x 130	≤ 2,000	≤ 984	46ASTM-50	46HASTM-50
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 50 mL						

ROTOFIX 46 **ASTM** PACKAGE - 100A

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (HEATED)
	4	100	58 x 161	≤ 2,000	≤ 984	46ASTM-100A	46HASTM-100A
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 100 mL						

ROTOFIX 46 **ASTM** PACKAGE - 100B

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (HEATED)
	4	100	44 x 168	≤ 2,000	≤ 961	46ASTM-100B	46HASTM-100B
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 100 mL						

ROTOFIX 46 **ASTM** PACKAGE - 100C

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (HEATED)
	4	100	37 x 200	≤ 2,000	≤ 961	46ASTM-100C	46HASTM-100C
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 100 mL						

ROTOFIX 46 **ASTM** PACKAGE - 100D

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (HEATED)
	4	100	45 x 165	≤ 2,000	≤ 961	46ASTM-100D	46HASTM-100D
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters 100 mL						

ROTOFIX 46 **CHROME BATH** PACKAGE

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (HEATED)
	6	30	24 x 151	≤ 2,000	≤ 917	46CHROME	46HCHROME
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (6) Buckets						

ROTOFIX 46 **SCHLENK** PACKAGE

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (HEATED)
	6	25	24 x 146.5	≤ 2,000	≤ 783	46SCHLENK	46HSCHLENK
	6	50	38 x 148.5	≤ 2,000	≤ 805		
Includes: (1) Centrifuge, (1) Horizontal Rotor, (6) Adapters							

ROTOFIX 46 **BABCOCK** PACKAGE

	Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF	PKG. NO.	PKG. NO. (HEATED)
	4	--	165	≤ 934	≤ 215	46BABCOCK	46HSBABCOCK
	Includes: (1) Centrifuge, (1) Horizontal Rotor, (4) Buckets, (4) Adapters, (4) Inserts						

ROTOFIX 46 **CUSTOM** PACKAGE

	Additional configurations are available for this unit, please inquire						
	Tube types include: ASTM, Schlenk, Babcock, Standard up to 250 mL, Conical and Skirted up to 50 mL, Bottles up to 290 mL, Specialty						

Learn more: hettweb.com/rotofix46

Hettich | North America
 100 Cummings Center, Ste. 136L
 Beverly, MA 01915, U.S.A.
 Toll free: 1 (866) 370-4388
info@hettweb.com
hettweb.com

© 2019 Hettich Instruments, LP. All rights reserved.
 Lit. No. PS-ROTOFIX-46-EN-01-2019.

Hettich products are designed to help you achieve optimal results for your application and are built to perform to the specifications outlined in the operators manual. For application specific information and settings, please refer to your organizations standard operating procedure. As always, our Hettich representatives are here to help determine which Hettich products and/or accessories best fit your laboratory requirements.