

CENTRIFUGE CABINETS

ONE AND TWO DRAWER UTILITY CABINETS

For convenient storage and positioning. Hettich base cabinets for Hettich centrifuges save space on the laboratory bench by providing enhanced flexibility.

The cabinets lower the centrifuge to an ergonomic height for easy loading and unloading. Each cabinet provides ample storage for user manuals, rotors, buckets and accessories. Important documentation and accessories remain with the centrifuge at all times. All cabinets are of sturdy design and feature casters with two wheel locks, drawer locks and 2 keys.

- Save valuable bench space
- Store accessories in one location
- Ergonomic height for loading/unloading



MODEL	LOW and NARROW	TALL and NARROW	LOW and WIDE	TALL and WIDE
Number of Drawers	1	2	1	2
Dims (H x W x D)	10 x 22 x 25.5 in.	22 x 22 x 25.5 in.	10 x 30 x 25.5 in.	22 x 30 x 25.5
Centrifuge Compatibility	ROTOFIX 32 A, UNIVERSAL 320 / 320 R, ROTINA 380 / 380 R, ROTANTA 420, ROTANTA 460, ROTOFIX 46 / 46 H		ROTANTA 420 R, ROTANTA 460 R	
Cat. No.	4612-A	4613-A	4614-K	4615-A

Learn more: hettweb.com/products/centrifuge-cabinet/

Hettich | North America
 100 Cummings Center, Ste. 136L
 Beverly, MA 01915, U.S.A.
 Toll free: 1 (866) 370-4388
info@hettweb.com
hettweb.com

© 2019 Hettich Instruments, LP. All rights reserved.
 Lit. No. PS-CENTRIFUGE-CABINETS-EN-01-2019.

Hettich products are designed to help you achieve optimal results for your application and are built to perform to the specifications outlined in the operators manual. For application specific information and settings, please refer to your organizations standard operating procedure. As always, our Hettich representatives are here to help determine which Hettich products and/or accessories best fit your laboratory requirements.