

FLOOR CENTRIFUGES

ROTANTA 460 RF / ROTIXA 500 RS / ROTO SILENTA 630 RS



Simplicity, variety and functional design. With low deck heights, front facing control panels and easy-lift lids, Hettich floor model centrifuges are designed for simple and efficient user interaction. An open chamber allows for unimpeded access for easy loading, unloading and cleaning. Versatile rotors accommodate multiple applications eliminating the need to change out rotors or move to another centrifuge. Thoughtfully designed accessories such as the HettLiner and balancing weights expedite your workflow and improve the quality of your results.

Optimal conditions for better results. Robust temperature control systems offer accurate and reliable time-saving results. Without a windshield, air moves freely throughout the chamber for greater temperature consistency. Dual sensors offer precise

“at-sample” temperature readings for greater accuracy (ROTO SILENTA 630 RS). Gain control and save time with quick pre-cool and stand-by cooling functions. Choose from multiple acceleration and deceleration modes for smooth braking and best sample results.

Built in Germany for lasting performance. All Hettich centrifuges are made from high quality, durable materials for a long in-field service life. Each unit is backed by a limited 5-year warranty and offer robust temperature control systems, all metal containment, maintenance free motors, metal lid locks and easy to use control panels.

Reliable, quiet and safe. You can rely on our centrifuges to deliver optimal performance day-in and day-out in a safe, distraction free environment.

Learn more: hettweb.com/floor-centrifuges/

Uncomplicated performance

hettweb.com

HETTICH FLOOR CENTRIFUGES



MODEL	ROTANTA 460 RF	ROTIXA 500 RS	ROTO SILENTA 630 RS
Ideal for	Floor model, refrigerated centrifuge offering high-capacity while saving bench space. Also accommodates small pediatric blood bags.	Entry into whole blood processing. Spin up to (4) 1L blood bags or 200 blood tubes. Robust and steady temperature control for best sample integrity.	Engineered for whole blood processing and large sample volumes, both single and double buckets for 6 or 12 positions within the same rotor. Non-windshielded rotor.
Blood Bags	4 x 450 ml (pediatric)	4 x 1,000 ml	12 x 750 ml 6 x 1,000 ml
Blood Tubes (Up to 7 ml / 10 ml)	196 / 148	200 / 168	180 / 180
Conical Tubes (15 / 50 ml)	96 / 40	112 / 40	138 / 48
Bottles	4 x 750 ml	4 x 1,000 ml	6 x 2,000 ml
Max. RPM / RCF	15,000 / 24,400	11,500 / 18,038	6,000 / 6,520
Temperature Control	-20 °C to +40 °C	-20 °C to +40 °C	-20 °C to +40 °C
Dims (H x W x D)	18 x 30 x 28 in.	38 x 26 x 32 in.	38 x 32 x 40 in.
Warranty†	5 yr.	5 yr.	5 yr.

*Please inquire about the fit of your specific tube dimensions. † Warranty valid in U.S. and Canada only, see warranty statement for complete details

Please inquire about our **floor centrifuge packages**, each package comes complete with everything you need to perform your centrifugation.

Learn more: hettweb.com/floor-centrifuges/



Hettich | North America
 100 Cummings Center, Ste. 136L
 Beverly, MA 01915, U.S.A.
 Toll free: 1 (866) 370-4388
info@hettweb.com
hettweb.com

© 2019 Hettich Instruments, LP. All rights reserved.
 Lit. No. PS-FLOOR-CENTRIFUGES-EN-01-2019.

Hettich products are designed to help you achieve optimal results for your application and are built to perform to the specifications outlined in the operators manual. For application specific information and settings, please refer to your organizations standard operating procedure. As always, our Hettich representatives are here to help determine which Hettich products and/or accessories best fit your laboratory requirements.