

# INDUSTRIAL CENTRIFUGES

MICRO / BENCHTOP / FLOOR STANDING



**Complete range of industrial centrifuges and accessories for maximum versatility on and off the bench.** Our compact micro centrifuges for low throughput labs are available in a variety of capacities and deliver high RCF values for optimal sample results. For larger volumes with higher throughput, our benchtop centrifuges provide a variety of accessories accommodating multiple applications to be performed in the same centrifuge. Floor models and rolling base cabinets help you reclaim valuable bench space.

**Built in Germany for lasting performance.** All Hettich centrifuges are made from high quality, durable materials for a long in-field service life. Each unit is backed by a limited 5-year warranty and offer robust temperature control systems, all metal containment, maintenance free motors and easy to use control panels.

**Reliable, quiet and safe.** You can rely on our centrifuges to deliver optimal performance day-in and day-out in a safe and distraction free environment. All units are tested and certified by leading international regulatory agencies for quality and safety.

**Pre and post sales support for the life of your product.** Not only are Hettich products built for the long-term, our dedicated support team is with you as well. We are here to keep you up and running from the centrifuges first day on the bench throughout its long life span. We even maintain parts for our centrifuges for up to 10-years after product obsolescence.

*See reverse for a sampling of what we have to offer.*

**Learn more:** [hettweb.com/product-category/centrifuges/](http://hettweb.com/product-category/centrifuges/)

## HETTICH CENTRIFUGES FOR THE INDUSTRIAL LABORATORY

We offer a complete range of centrifuges for food and beverage, petrochemical, semiconductor, automotive, agriculture, or water and wastewater industries in either the private or government sector. Ask about our centrifuge packages that make it easy to order and come complete with everything you need to perform your application.

<b>MIKRO 185</b> 		<b>MIKRO 200   R</b> 		<b>MIKRO 220   R</b> 		<b>UNIVERSAL 320   R</b> 	
<b>Max. Capacity</b>	24 x 1.5 / 2.0 mL	<b>Max. Capacity</b>	30 x 1.5 / 2.0 mL	<b>Max. Capacity</b>	60 x 1.5 / 2.0 mL	<b>Max. Capacity</b>	4 x 200 mL
<b>Max. RPM/RCF</b>	14,000/18,845	<b>Max. RPM/RCF</b>	15,000/21,382	<b>Max. RPM/RCF</b>	18,000/31,514	<b>Max. RPM/RCF</b>	15,000/21,382
<b>Temp. Control</b>	Ambient	<b>Temp. Control</b>	-10 to +40°C (200 R)	<b>Temp. Control</b>	-20 to +40°C (220 R)	<b>Temp. Control</b>	-20 to +40°C (320 R)
<b>ROTINA 380   R</b> 		<b>ROTINA 420   R</b> 		<b>ROTOFIX 46</b> 		<b>ROTOFIX 46 H</b> 	
<b>Max. Capacity</b>	4 x 290 mL	<b>Max. Capacity</b>	4 x 600 mL	<b>Max. Capacity</b>	4 x 290 mL	<b>Max. Capacity</b>	4 x 100 mL
<b>Max. RPM/RCF</b>	15,000/24,400	<b>Max. RPM/RCF</b>	15,000/24,400	<b>Max. RPM/RCF</b>	4,000/3,095	<b>Max. RPM/RCF</b>	2,000/984
<b>Temp. Control</b>	-20 to +40°C (380 R)	<b>Temp. Control</b>	-20 to +40°C (420 R)	<b>Temp. Control</b>	Ambient	<b>Temp. Control</b>	-10 to +90°C
<b>ROTANTA 460   R</b> 		<b>ROTANTA 460 RF</b> 		<b>ROTIXA 500 RS</b> 		<b>ROTO SILENTA 630 RS</b> 	
<b>Max. Capacity</b>	4 x 750 mL	<b>Max. Capacity</b>	4 x 750 mL	<b>Max. Capacity</b>	4 x 1,000 mL	<b>Max. Capacity</b>	6 x 2,000 mL
<b>Max. RPM/RCF</b>	15,000/24,400	<b>Max. RPM/RCF</b>	15,000/24,400	<b>Max. RPM/RCF</b>	11,500/18,038	<b>Max. RPM/RCF</b>	6,000/6,520
<b>Temp. Control</b>	-20 to +40°C (460 R)	<b>Temp. Control</b>	-20 to +40°C	<b>Temp. Control</b>	-20 to +40°C	<b>Temp. Control</b>	-20 to +40°C

**Note:** Max. specifications based on research product configuration not the general capability of the instrument. Please inquire for more detail on your specific tube size and type.

Learn more: [hettweb.com/product-category/centrifuges/](http://hettweb.com/product-category/centrifuges/)



**Hettich | North America**  
 100 Cummings Center, Ste. 136L  
 Beverly, MA 01915, U.S.A.  
 Toll free: 1 (866) 370-4388  
[info@hettweb.com](mailto:info@hettweb.com)  
[hettweb.com](http://hettweb.com)

© 2019 Hettich Instruments, LP. All rights reserved.  
 Lit. No. PS-INDUSTRIAL-CENTRIFUGES-EN-01-2019.

Hettich products are designed to help you achieve optimal results for your application and are built to perform to the specifications outlined in the operators manual. For application specific information and settings, please refer to your organizations standard operating procedure. As always, our Hettich representatives are here to help determine which Hettich products and/or accessories best fit your laboratory requirements.