

EBA MODEL 'S' FOR RAPID SEPARATION

OF STAT CHEMISTRY, COAG. & PEDIATRIC SAMPLES

Introducing Hettich **EBA 200 S**
and **EBA 280 S** small centrifuges

Rapid separation for a quick turnaround

These small, high-performance, centrifuges provide rapid separation for a quick turnaround of your STAT chemistry, coag. and pediatric samples. Models are compact and achieve consistent platelet counts in just minutes. Max. speeds of **6,000 RPM** in a swing-out rotor and **8,000 RPM** in fixed-angle rotor.

Made in Germany for lasting performance

Each model is made in Germany to provide years of quiet, safe and worry-free operation. Products are manufactured under strict quality assurance guidelines and are certified by all major regulatory agencies. In addition, Hettich centrifuges are backed by a limited five year warranty.

Support for the life of your product

Hettich provides free product and application support for the life of your product and maintains service parts for up to ten years after product obsolescence. We have manufactured lab equipment for more than a century and plan to continue for a century more.

End-of-cycle notification for an efficient workflow

The Hettich EBA 200 S and EBA 280 S come standard with both visual and acoustical end-of-cycle notifications for an efficient workflow. Choose from an acoustical signal that will sound in 30-sec. intervals, a visual signal via a blinking control panel or the auto-opening of the centrifuge lid.

Maintain quality control and reduce user error

User-defined programs to optimize turnaround time.



EBA 200 S AND EBA 280 S FOR RAPID SEPARATION OF STAT CHEMISTRY, COAG. & PEDIATRIC SAMPLES

Rapid separation for a quick turnaround

- STAT specimen processing
- Gel, coag. and pediatric tubes
- Clean separation in just minutes
- End-of-cycle notification
- Small and compact with choice of angle or horizontal flat separation
- Small footprint is ideal for both single and multi-unit setups
- User-defined programs to optimize turnaround time
- Free support for life of product
- Maintenance parts for up to ten years after product obsolescence
- Low maintenance
- Made in Germany
- 5-year limited warranty



RECOMMENDED CENTRIFUGE PACKAGES

EBA 200 S BLOOD TUBE PACKAGE - 8

Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
8	1.6 - 5	13 x 75	≤ 8,000	≤ 4,794
8	4 - 7	13 x 100	≤ 8,000	≤ 6,153
8	8.5 - 10	16 x 100	≤ 8,000	≤ 6,153
4	8	16 x 125	≤ 8,000	≤ 6,153
Pkg. No.		200SBL00D8		
Includes: (1) Centrifuge, (1) Built-in Angle Rotor, (8) Adapters 13 x 75 mm, (8) Adapters 13 x 100 mm				

EBA 200 S PEDIATRIC TUBE PACKAGE - 8

Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
8	0.5	10.7 x 36	≤ 8,000	≤ 3,935
Pkg. No.		200SBL00D8-PED		
Includes: (1) Centrifuge, (1) Built-in Angle Rotor, (8) Adapters 0.5 mL				

EBA 280 S BLOOD TUBE PACKAGE - 6

Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
6	1.6 - 5	13 x 75	≤ 6,000	≤ 4,266
6	4 - 7	13 x 100	≤ 6,000	≤ 5,071
6	4 - 7	16 x 75	≤ 6,000	≤ 4,669
6	8.5 - 10	16 x 100	≤ 6,000	≤ 5,071
Pkg. No.		280SBL00D6		
Includes: (1) Centrifuge, (1) Horizontal Rotor, (6) Buckets, (6) Adapters 13 x 75 mm, (6) Adapters 16 x 75 mm				

EBA 280 S PEDIATRIC TUBE PACKAGE - 6

Tube Qty.	Volume (mL)	Size (mm)	RPM	RCF
6	0.5	10.7 x 36	≤ 6,000	≤ 3,059
Pkg. No.		280SBL00D6-PED		
Includes: (1) Centrifuge, (1) Horizontal Rotor, (6) Buckets, (6) Adapters				

MODEL SPECIFICATIONS

MODEL	EBA 200 S	EBA 280 S
Rotor Type	Fixed angle	90° Swing-out
Max. Capacity	8 Blood tubes	6 Blood tubes
Max. RPM/RCF	8,000 min ⁻¹ / 6,153	6,000 min ⁻¹ / 5,071
Temperature Control	Ambient	Ambient
Dims (H x W x D)	9 x 10 x 14 in.	9.5 x 12.8 x 15.75 in.
Weight	24 lbs	24.25 lbs

Learn more at: hettweb.com/eba200

Learn more at: hettweb.com/eba280

Contact your local Hettich representative:

Hettich | North America
 100 Cummings Center, Ste. 136L
 Beverly, MA 01915, U.S.A.
 Toll free: 1 (866) 370-4388
info@hettweb.com