

AUTOMATIC TUBE STRIPPER

— In partnership with LMB

Features:

- New and redesigned tube holder for beater tube guidance
- Designed for right and left hand use
- Low weight for easy handling
- Swift and complete stripping of blood bag tubes
- Ideal for all commercial blood bag systems on the market
- Automatically centered tubing fed through pair of rollers by a guidance system
- Preventing coagulation of blood in tubing prior to preparing segments
- Power supply with stabilized output, low ripple & low interference, Short circuit & overload protection
- High efficiency and low energy consumption



SPECIFICATIONS

Overall dimensions (handle closed)

Overall dimensions (handle open)

Max stripping tube diameter

Stripping wheel diameter

Stripping wheels angle opening

Tube stripper body diameter

Max stripping speed

Weight without power supply

Operating temperature

Storage and transport temperature

Power supply voltage

Input frequency

DC motor supply voltage

Warranty†

TS08

3 x 2.25 x 7.75 in.

4.5 x 2.25 x 7.75 in.

0.25 in.

1.3 in.

25 deg

1.5 in.

15.75 in./sec

~1.75 lbs.

+10 to +40 °C

-20 to +50 °C

100-240 V AC

50 Hz / 60 Hz

12V

1 year

† Warranty valid in U.S. and Canada only, see warranty statement for complete details.

Learn more: hettweb.com

Contact your local Hettich representative:

Nick Sillings | Northeast
nick.sillings@hettweb.com

Matt Musheff | Great Lakes
matthew.musheff@hettweb.com

Kirk Smith | Southeast
kirk.smith@hettweb.com

Christy Gonzalez | West
christy.gonzalez@hettweb.com

Lynne Wilhite | Midwest
lynne.wilhite@hettweb.com

Christina D'Anna | Northwest
christina.danna@hettweb.com

Victoria Corral | Southwest
victoria.corral@hettweb.com

Ann Chainey | Canada
ann.chainey@hettweb.com

Alexandra McPartland | Sales
alexandra.mcpartland@hettweb.com

Ian Hanegraaff | Automation
ian.hanegraaff@hettweb.com

Dustin Tharp | Service Contracts
dustin.tharp@hettweb.com