

MANUAL TUBE STRIPPER

— Self-centering roller system

Designed to strip, crimp and cut tubing in the preparation of blood segments from the blood collection and processing procedures. Along with the spring loaded handle, the self-centering roller systems keeps the tube in place and prevents slipping or falling off the line.

- High quality aluminum alloy
- Anti-rust, light weight and sturdy
- Easy cleaning by wiping or soaking in IPA
- Ergonomic design for good grip and stripping
- Ease fatigue over prolonged use
- Self-centering roller system keeps tubing in position
- Smooth traction for effective stripping of blood components
- Internal spring system keeps handles in open position
- Minimize entanglement with blood tubing



SPECIFICATIONS	TS-002	TS-003
Base Material	Aluminum Alloy & Stainless Steel	Aluminum Alloy & Stainless Steel
Rollers	Special Composite Steel	Special Composite Steel
Cutters	Not Applicable	Stainless Steel
Features	Self-Centering Rollers	Self-Centering Rollers
Cutter	No	Yes
Net Weight (g) +/-5%	68	70
Warranty†	1 year	1 year

† Warranty valid in U.S. and Canada only, see warranty statement for complete details.

Learn more: hettweb.com

Contact your local Hettich representative:

Nick Sillings | Northeast
nick.sillings@hettweb.com

Matt Musheff | Great Lakes
matthew.musheff@hettweb.com

Kirk Smith | Southeast
kirk.smith@hettweb.com

Christy Gonzalez | West
christy.gonzalez@hettweb.com

Lynne Wilhite | Midwest
lynne.wilhite@hettweb.com

Christina D'Anna | Northwest
christina.danna@hettweb.com

Victoria Corral | Southwest
victoria.corral@hettweb.com

Ann Chainey | Canada
ann.chainey@hettweb.com

Alexandra McPartland | Sales
alexandra.mcpartland@hettweb.com

Ian Hanegraaff | Automation
ian.hanegraaff@hettweb.com

Dustin Tharp | Service Contracts
dustin.tharp@hettweb.com