

ROTOFIX 46 H

BENCHTOP CENTRIFUGE

— ASTM / Babcock

The **ROTOFIX 46 H** is a robust, versatile benchtop centrifuge ideal for research, industrial, chemical, and environmental testing laboratories. Applications include petrochemical and fuel, industrial fluids, and wastewater. It is available with a swing-out rotor for a maximum capacity of 4 x 100 mL and buckets for ASTM (pear-shaped) and Babcock tubes. With an adjustable maximum speed of 2,000 RPM, a maximum RCF of 984, and a low noise level of ≤ 58 dB(A), it features an impulse button for short centrifugation, easy keypad operation, and an heating range from +10 to +90 °C. Made in Germany for lasting quality, shipped and supported in North America, and backed by a five-year warranty, the ROTOFIX 46H delivers reliable performance.

Advantages:

- Heating range: +10 to +90 °C
- Made in Germany for lasting quality and performance
- Shipped and supported in North America
- Five (5) year warranty



SPECIFICATIONS	ROTOFIX 46 H (Heated)	ASTM Tubes
Max. Capacity	4 x 100 ml	
Max. RPM / RCF	2,000 / 984	
Max. ASTM Tubes*	4 (up to 100 ml)	
Temperature Control	+10 °C to +90 °C	
Dims (H x W x D)	14 x 21 x 25 in.	
Weight	132 lbs.	
Power	200 - 240 V 1 ~ / 50 - 60 Hz	
Warranty†	5 years	

*Please inquire about the fit of your specific tube dimensions. † Warranty valid in U.S. and Canada only, see warranty statement for complete details.

Learn more: hettweb.com/rotofix46h




ROTOFIX 46 H CENTRIFUGE PACKAGES


ROTOFIX 46H ASTM PACKAGE - 50

	Tube Qty	Volume (mL)	Size (mm)	RPM	RCF	ROTOFIX 46 H
	4	50	45 x 130	2,000	984	Temp. Control +10 °C to +90 °C
PKG. NO. (Packages include the centrifuge, rotor, buckets, adapters, and lids when applicable)						46HASTM-50


ROTOFIX 46H ASTM PACKAGE - 100A

	Tube Qty	Volume (mL)	Size (mm)	RPM	RCF	ROTOFIX 46 H
	4	100	58 x 161	2,000	984	Temp. Control +10 °C to +90 °C
PKG. NO. (Packages include the centrifuge, rotor, buckets, adapters, and lids when applicable)						46HASTM-100A


ROTOFIX 46H ASTM PACKAGE - 100B

	Tube Qty	Volume (mL)	Size (mm)	RPM	RCF	ROTOFIX 46 H
	4	100	44 x 168	2,000	961	Temp. Control +10 °C to +90 °C
PKG. NO. (Packages include the centrifuge, rotor, buckets, adapters, and lids when applicable)						46HASTM-100B


ROTOFIX 46H ASTM PACKAGE - 100C

	Tube Qty	Volume (mL)	Size (mm)	RPM	RCF	ROTOFIX 46 H
	4	100	37 x 200	2,000	961	Temp. Control +10 °C to +90 °C
PKG. NO. (Packages include the centrifuge, rotor, buckets, adapters, and lids when applicable)						46HASTM-100C

ROTOFIX 46H ASTM PACKAGE - 100D

	Tube Qty	Volume (mL)	Size (mm)	RPM	RCF	ROTOFIX 46 H
	4	100	45 x 165	2,000	961	Temp. Control +10 °C to +90 °C
PKG. NO. (Packages include the centrifuge, rotor, buckets, adapters, and lids when applicable)						46HASTM-100D

ROTOFIX 46H BABCOCK PACKAGE

	Tube Qty	Volume (mL)	Size (mm)	RPM	RCF	ROTOFIX 46 H
	4	--	165	934	215	Temp. Control +10 °C to +90 °C
PKG. NO. (Packages include the centrifuge, rotor, buckets, adapters, and lids when applicable)						46HSBABCOCK

Learn more: hettweb.com/rotofix46h

Contact your local Hettich representative: